


Front Riggings Attachment Kit (FRKTM6)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1107807-NLA	NLA - Tube, Support, Front Rigging, Black	1	2
2		1026900	Washer (1/4 x 23/32 x 1/16")	1	10
3		7930345	Locknut (1/4-28)	25	4
4		7930333	Locknut (1/4-20)	25	2
5		1027097	Plug Button, Black (3/4")	1	2
6		1078020-NLA	NLA - Hanger, Front Rigging, Pin Style, Black - Right	1	1
6		1078021-NLA	NLA - Hanger, Front Rigging, Pin Style, Black - Left	1	1
7		1047449-NLA	NLA - Spacer (1/4 x 3/8 x 3/16")	1	4
8		1030949	Screw, Hex Head (1/4-20 x 1-1/2")	1	2
9		1026099	Spacer, Coved (1/4 x 11/16 x 1/2")	1	4
10		1059374	Screw, Hex Head (1/4-28 x 3")	1	4
11		1107808-NLA	NLA - Kit, Front Rigging Attachment with Instruction Sheet	1	1
	A	1106633-NLA	NLA - Instruction Sheet, Front Rigging Attachment Kit	2	1

NOTE: A - Not Shown