

Battery Box Hangers

NOTE: LABAC Chairs use the same configurations as the Non-LABAC Chairs.


Chair Description	Side	14" or 15" Wide	16" thru 18" Wide	19" thru 24" Wide
4M or 5M with Standard Crossbraces Sectional Back or Manual Recliner	Right	1029690-NLA	1027935-NLA	1033959-NLA
	Left	1029691-NLA	1027936-NLA	1033960-NLA
4M or 5M with Dual Crossbraces or Depth over 20" Sectional, Non-Sectional or Manual Recliner	Right	1035848-NLA	1033961-NLA	1033959-NLA
	Left	1035849-NLA	1033962-NLA	1033960-NLA
4M with L-2P and Standard Crossbraces Sectional Back or Manual Recliner	Right	1035576-NLA	1033961-NLA	N/A
	Left	1035577-NLA	1033962-NLA	N/A
4M with L-2P and Dual Crossbraces Sectional Back	Right	1035576-NLA	1033961-NLA	1033959-NLA
	Left	1035577-NLA	1033962-NLA	1033960-NLA
4M with L-2P and Standard Crossbraces Manual Recliner	Right	1035997-NLA	1033961-NLA	1027935-NLA
	Left	1035996-NLA	1033962-NLA	1027936-NLA
4M with L-2P and Dual Crossbraces Manual Recliner	Right	Order as a Special	1033961-NLA	N/A
	Left	Order as a Special	1033962-NLA	N/A
4M with Standard Crossbraces Power Back Recliner	Right	N/A	1027935-NLA	N/A
	Left	N/A	1027936-NLA	N/A

Battery Box Hangers

NOTE: LABAC Chairs use the same configurations as the Non-LABAC Chairs.

Mounting Bracket to Raise Battery Box

Sectional Backs and Manual Recline Chairs Over 14" Wide



Battery Box Hangers

NOTE: LABAC Chairs use the same configurations as the Non-LABAC Chairs.

Description	Part Number
Barrier Plate (Not Shown) 2 each needed when raising box. Used to protect batteries from rivets.	1027604-NLA
Washer 3/16"	16348X003
Mounting Bracket	1034093-NLA
Screw	60284X006-NLA
Pop Rivet	1007094-NLA