

Frame Width Growth Kits

NOTE: Frame width growth kits listed are for power Tiger seats only. To widen an Orbit seat, see seat frame under Orbit seat option.

NOTE: Front Rigging Footplates may be needed when growing a chair.

Original Width	13" New Width	14" New Width	15" New Width	16" New Width
12"	N/A See Note(s) A	1040142-NLA	1040143-NLA See Note(s) B	1040143-NLA See Note(s) B
13"	N/A	1040142-NLA	1040143-NLA See Note(s) B	1040143-NLA See Note(s) B
14"	N/A	N/A	1040143-NLA See Note(s) B	1040143-NLA See Note(s) B
15"	N/A	N/A	N/A	N/A See Note(s) A

NOTE: A - Only wider upholstery is required.
B - Specify frame color when ordering.

Individual Hardware For Growth Kits

Part Number	Description	1040142-NLA	1040143-NLA
1037291-NLA	Spacer Tube, 4" (for 16" Wide Frame)	N/A	2
1025360-NLA	Spacer Tube 2" (for 14" Wide Handle)	2	N/A
1029930-NLA	Spacer Tube 4" (for 16" Wide Handle)	N/A	2
1029922-NLA	Release Rod	N/A	1
1037285-NLA	Slide Tube, Seat Frame	N/A	2
1038250-NLA	Slide Tube, Cane	N/A	2
1016254	Handgrip	2	2
8881037655 See Note(s) A	Seat Rail Tube (15"-16" Wide)	N/A	2
1036434-NLA	Washer (9/32 x 7/8 x 3/16")	N/A	2
1036040-NLA	Spacer (9/32 x 7/8 x 1/2")	8	N/A
1020363-NLA	Spacer (1/4 x 1/2 x 7/32")	2	N/A
1036840	Screw (1/4-20 x 2-1/4")	6	N/A
1026035-NLA	Locknut (10-24)	N/A	2
1030957-NLA	Shoulder Bolt (1/4 x 1-3/4")	N/A	2
1026033	Washer (1/4 x 1/2 x 1/16")	N/A	2
1026034-NLA	Washer (7/32 x 7/16 x 1/8")	N/A	2
1002108-NLA	Screw (1/4-20 x 2")	4	N/A

Frame Width Growth Kits

NOTE: Frame width growth kits listed are for power Tiger seats only. To widen an Orbit seat, see seat frame under Orbit seat option.

NOTE: Front Rigging Footplates may be needed when growing a chair.

Individual Hardware For Growth Kits

Part Number	Description	1040142-NLA	1040143-NLA
40157X020	Spacer, Coved (1")	N/A	2
7930301-NLA	Screw (Pkg of 25) (1/4-20 x 1-3/4")	N/A	6
1008495	Screw (1/4-20 x 1-1/2")	N/A	4

NOTE: A - Specify frame color when ordering.