



Rear Upright Tube Mounting Hardware



Rear Upright Tube Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1059271	Tip, Plastic, Black	1	1
2		1058545	Lever, Actuator	1	1
3		40157X020	Spacer, Coved (1/4 x 1 x 9/16")	1	2
4		1011511	Spring, Compression (5/8 x 1 x 3/64")	1	1
5		1058232	Plug, Latch	1	1
6	A	1055350	Screw, Hex Head (1/4-28 x 2-1/4")	1	1
6	B	1055349	Screw, Hex Head (1/4-28 x 2")	1	1
7		1059109	Plate, Backrest Mounting	1	2
8		1059127	Pin, Locking	1	1
9		1025765	Locknut (1/4-28)	1	5
10		1058832-NLA	NLA - Tube, Backrest, Black	1	1
11		1045495-NLA	NLA - Screw, Set Socket (1/4-20 x 1/4")	1	2
12		See Chart	#12 Backrest Crossmember	1	1
NOTE: A - With Transport Option B - Without Transport Option					

#12 Backrest Crossmember

12" Wide	13" Wide	14" Wide	15" Wide	16" Wide
1042508-NLA	1042509-NLA	1042510	1042511	1042512
17" Wide	18" Wide	19" Wide	20" Wide	
1050083	1050080	1050077	1050074	