

**MyOn®
United States
Junior (FXMYONJR, GRMYONJR)
Adult (MYONADLT, MYONTS, HDPG)
Ultralightweight Folding Wheelchair**



Yes, you can.

Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 11-013. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-333-6900

For Faxed In Orders Call:

1-800-678-4682 or

1-440-366-0709

For Technical Service Call:

1-800-832-4707

Literature Can Be Found at:

www.invacare.com or

1-800-828-6272

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers

© 2022 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

TABLE of CONTENTS

Description	Page Number
Accessories - Acrylic Trays	129
Accessories - Anti-Tippers (A1580, A1580P)	132
Accessories - Anti-Tippers, Long (For Model MYONADLT 19 and 20" FSTF)	135
Accessories - Crutch/Cane Carrier (A1569)	137
Accessories - Headrest and Backbrace	138
Accessories - Junior Spoke Protectors	140
Accessories - One Arm Drive (DUHARU, DUHARP)	141
Accessories - Positioning Straps	143
Accessories - Push Handles	144
Accessories - Stabilizer Kit	147
Accessories - Standard Spoke Protectors	148
Accessories - Stroller Handle (STRHNDL)	149
Arms - Cantilever - Clamp On - Desk Length (AT55D), Full Length (AT55F)	26
Arms - Conventional Adjustable Height - Desk Length (CAD), Full Length (CAF)	30
Arms - Conventional Fixed Height - Desk Length (CNAD), Full Length (CNAF)	34
Arms - Side Guards - Rigid	54
Arms - Side Guards, Clothing Guards, and Mud Guards	51
Arms - Swing Away - Front Mount - Transport Brackets	38
Arms - Swing Away - Rear Mount - Transport Brackets	40
Arms - T Style Adjustable - Desk Length (TAD), Full Length (TAF)	45
Arms - T Style Adjustable - Heavy Duty - Desk Length (TAD), Full Length (TAF)	48

TABLE of CONTENTS

Description	Page Number
Arms - T Style Adjustable - Junior and Active (TAHD, TAHJ)	42
Back Canes - Adjustable Angle Back (AABH - 8 Degree Bend)	69
Back Canes - Adjustable Angle Back (AABH - Straight)	66
Back Canes - Adjustable Height w/ Push Handles (ABHPH)	62
Back Canes - Adjustable Height w/o Push Handles (ABH)	64
Back Canes - Fixed Height w/ Push Handles (FBH - FBHF)	58
Back Canes - Fixed Height w/ Push Handles, HD (FBH - FBHF)	60
Casters - Caster, Headtube and Fork w/ Mounting Hardware	103
Crossbrace - Conventional Adjustable Depth (ADJDPTH)	13
Crossbrace - Conventional Adjustable Depth, Heavy Duty (HDPG ADJDPTH)	18
Crossbrace - Conventional Fixed Depth (FXDDPTH)	11
Crossbrace - Conventional Fixed Depth, Heavy Duty (HDPG FXDDPTH)	16
Crossbrace - Conventional Fixed Seat Depth	24
Crossbrace - Conventional Growable Seat Depth	21
Footplate Assembly - 80 Degree Supports	97
Footplate Assembly - Adult Heavy Duty (FLIPFTS, FLDPADJ, and 1PCFTPT)	95
Footplate Assembly - Adult Standard (FLIPFTS, FLDPADJ, and 1PCFTPT)	92
Footplate Assembly - Junior Elevating Leg Rests (ELVOFF, ELVOFFWO)	101
Footplate Assembly and Lift Off Pivot Hanger, 90 Degree (90STRLIFT)	99
Frame Finishes and Touch-Up Paint	6
Frame with Hardware	7

TABLE of CONTENTS

Description	Page Number
Front Riggings - Accessories	91
Front Riggings - Adult Elevating Legrest Calf Pads (ELRSS)	85
Front Riggings - Adult Elevating Legrest Hanger Assembly (ELRSS)	84
Front Riggings - Junior 80 Degree Swing Hangers	87
Front Riggings - Junior Elevating Legrest Hangers and Calf Pad (ELVOFF, ELVOFFWO)	89
Front Riggings - Swing Out Hangers - 70, 80, and 90 Degree	82
Transport (TRBKTS) Brackets	56
Upholstery - Back Brace and Hardware	71
Upholstery - Back Covers, Tensioners	74
Upholstery - Seat - Junior (Fixed and Adjustable)	80
Upholstery - Seat - Standard (FXDDPTH, ADJDPTH)	78
Wheel Locks - Hill Holder Assembly (A1598)	127
Wheel Locks - Pull to Lock (WLPULL), Pull to Lock Extended (WLPULLX)	123
Wheel Locks - Push to Lock (WLPUSH), Push to Lock Extended (WLPUSHX)	125
Wheels - Rear - Axle Plates - 20 Inch Frame	113
Wheels - Rear - Axle Plates w/ Hardware (A1573 - All Wheels)	106
Wheels - Rear - Axle Plates w/ Hardware (A1577 - Composite Wheels)	108
Wheels - Rear - Axle Plates w/ Hardware (A1577 - Spoked Wheels)	110
Wheels - Rear - Axle Plates w/ Hardware (Permanent Axle - All Wheels)	112
Wheels - Rear - Composite, Spoked, Sun Fusion, and Spinergy	116
Wheels - Rear - Handrims and Mounting Hardware	120


TABLE of CONTENTS

Description	Page Number
Wheels - Rear - Natural Fit Handrims and Hardware	122


Frame Finishes and Touch-Up Paint

Color Numbers	Color	Touch-Up Paint (.6 oz.)
24P	Wet Black	1031450-NLA
30P	Sunny Yellow	1074167-NLA
60P	Silver Metallic	1101318-NLA
61P	Electric Red	1046321-NLA
62P	Electric Blue	1046320-NLA
71P	Silver Vein	N/A
104P	Electric Teal	1101316-NLA
105P	Midnight Blue	1101326-NLA
115P	Black Prism	1133371-NLA
119P	Bubblegum Pink	1133375-NLA
120P	Grape Madness	1133376-NLA
121P	Lollypop Blue	1133377-NLA
122P	Cosmic Blue	1133378-NLA
123P	Celestial Purple	1133379-NLA
125P	Emerald Green	1133381-NLA
127P	Grasshopper Green	1133383-NLA

Frame with Hardware



Frames, Back Canes
and Forks Shown for
Reference Only



Frame with Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1172219	Kit, H-Block Hardware	1	1
3			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
4			H-Block, Front - Left	1	1
4			H-Block, Front - Right	1	1
5	B	1172246	NLA - Locknut (M6 x 1)	1	2
2			Kit, H-Block Hardware, Junior	1	1
			Locator, Seat	1	2
			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
	B	1172248	NLA - Locknut (M6 x 1)	1	2
2			Kit, H-Block Hardware, Junior (12 - 15" W)	1	1
3			Screw, Socket Head, Cap (M6 x 35mm)	1	2
5			NLA - Locknut (M6 x 1)	1	2
6	C	1172217	H-Block, Front - Right	1	1
6			H-Block, Front - Left	1	1
7			Kit, TRBKTS Mounting Hardware	1	1
8			H-Block, Front - Right	1	1
8			H-Block, Front - Left	1	1
9			Bracket, Tie Down, Front	1	2
10			Plate, Tie Down, Rear	1	2
11			nla - Screw, Hex Head w/ Patch (M6 x 1 x 50mm)	1	4
12			Screw, Hex Head Cap (M6 x 1 x 60)	1	2
13			Screw, Hex Head Cap (M6 x 1 x 50)	1	2
14			NLA - Spacer	1	2
15			NLA - Screw, Socket Head Cap (M6 x 1 x 40mm)	1	2
16			Washer (1/4 x 1/2 x 1/32")	1	12
		See Chart	Frame Assembly w/ Hardware, Standard	1	1
		See Chart	Frame Assembly w/ Hardware, Heavy Duty	1	1
		See Chart	Frame Assembly w/ Hardware - TRRO/TRBKTS	1	1

NOTE: A - Kit includes items 3, 4, and 5.
 B - Kit includes items 3, 5, and 6.
 C - Kit includes items 8 - 16.

Frame Assembly w/ Hardware, Standard

Frame Style	Seat Depth	Part Number	Hardware Part No.
MYONADLT MYONTO MYONTS	15"	8881172238	1172219
MYONADLT MYONTO MYONTS	16 - 18"	8881172244	1172219
MYONADLT MYONTO MYONTS	19 - 20"	8881172239	1172219

Frame with Hardware

Specify Frame Color when ordering, see "Frame Finishes" Category for Color choices.

Frame Assembly w/ Hardware, Heavy Duty

Frame Style	Seat Depth	Part Number	Hardware Part No.
MYONAC - HDPG	15 - 16"	8881172241	1172219
MYONAC - HDPG	17 - 18"	8881172242	1172219
MYONAC - HDPG	19 - 20"	8881172243	1172219

Specify Frame Color when ordering, see "Frame Finishes" Category for Color choices.


Frame with Hardware

Frame Assembly w/ Hardware - TRRO/TRBKTS

Frame Style	Seat Depth	Part Number	Hardware Part No.
MYONADLT MYONTO	15"	8881172238	1172216
MYONADLT MYONTO	16 - 18"	8881172244	1172216
MYONADLT MYONTO	19 - 20"	8881172239	1172216
MYONADLT MYONTS	15"	8881172238	1172217
MYONADLT MYONTS	16 - 18"	8881172244	1172217
MYONADLT MYONTS	19 - 20"	8881172239	1172217

Specify Frame Color when ordering, see "Frame Finishes" Category for Color choices.

Crossbrace - Conventional Fixed Depth (FXDDPTH)




Crossbrace - Conventional Fixed Depth (FXDDPTH)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Crossbrace Assembly w/ Hardware - Myon	1	1
3	B	1173765	Kit, Saddle Clamp w/ Hardware (Pair)	1	1
4			Clamp, Saddle, Crossbar	1	2
5			NLA - Screw, Socket Button Head (M6 x 40mm)	1	4
6			NLA - Locknut (M6 x 1)	1	4
7	C	1173771	NLA - Kit, Crossbrace Bolt w/ Hardware (Pair)	1	1
8			Washer, Plastic	1	2
9			NLA - Screw, Hex Head (M8 x 1.25 x 80mm)	1	2
10			NLA - Washer (M8 x 16mm x 1.5mm)	1	4
11			NLA - Locknut (M8 x 1.25)	1	2
12	D	1173779-NLA	NLA - Kit, Pivot Link w/ Hardware (15 - 17" W) (Pair)	1	1
13			NLA - Link, Pivot, Crossbrace (15 - 17" W)	1	2
14			NLA - Locknut (M6 x 1)	1	2
15			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12	D	1173767-NLA	NLA - Kit, Pivot Link w/ Hardware (18 - 20" W) (Pair)	1	1
13			Link, Pivot, Crossbrace (18 - 20" W)	1	2
14			NLA - Locknut (M6 x 1)	1	2
15			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
NOTE: A - Includes items 3 - 19. B - Includes items 4 - 6. C - Includes items 8 - 11. D - Includes items 13 - 15.					

#1 Crossbrace Assembly w/ Hardware - Fixed Depth

Seat Depth	15" Wide	16" Wide	17" Wide	18" Wide	19" Wide	20" Wide
15" Deep	1172249	1172250	1172251	1172252	1172253	1172254
16" Deep	1530032	1530033	1530034	1530035	1530036	1530037
17" Deep	1168733	1168734	1168735	1168736	1168737	1168738
18" Deep	1530038	1530052	1530053	1530054	1530055	1530056
19" Deep	1172255	1172256	1172257	1172258	1172259	1172260
20" Deep	1172261	1172262	1172263	1172264	1172265	1172266

Crossbrace - Conventional Adjustable Depth (ADJDPTH)



Crossbrace - Conventional Adjustable Depth (ADJDPTH)


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Crossbrace Assembly w/ Hardware - Myon	1	1
2	B	1173766	Kit, Saddle Clamp w/ Hardware (Pair)	1	1
3			Clamp, Saddle, Crossbar	1	2
4			Screw (M6 x 45mm)	1	4
5			NLA - Locknut (M6 x 1)	1	4
6	C	1173775	Kit, Crossbrace Bolt w/ Hardware (Pair)	1	1
7			Insert, Crossbrace	1	4
8			Screw, Socket Head, Cap (M6 x 55mm, 12.9)	1	2
9			Washer, Crossbrace	1	2
10			NLA - Locknut (M6 x 1)	1	4
11	D	1173768	Kit, Pivot Link w/ Hardware (18 - 20" W) (Pair)	1	1
12			Link, Pivot, Crossbrace (18-20" W)	1	2
13			Screw, Socket Head, Cap (M6 x 35mm)	1	2
14			NLA - Locknut (M6 x 1)	1	2
			Carton	1	1
			NLA - Bag, Poly Zip-Lip (8-3/4" L x 6" W)	1	1
11	D	1173780	Kit, Pivot Link w/ Hardware (15 - 17" W) (Pair)	1	1
12			Link, Pivot, Crossbrace (15-17" W)	1	2
13			Screw, Socket Head, Cap (M6 x 35mm)	1	2
14			NLA - Locknut (M6 x 1)	1	2
15	E	1172274	Kit, Seat Rail Extensions (16 - 18" D)	1	1
16			Extension, Seat Rail (2")	1	2
17			Extension, Seat Rail (1")	1	2
18			NLA - Plug Button	1	2
19			Cap, Extension, Seat Rail (1")	1	2
20			Rod, Upholstery (2")	1	2
21			Rod, Upholstery (1")	1	2
22			Guide, Seat	1	2
23			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
24			NLA - Screw, Phillips Pan Head (M3.5 x .6 x 9.5mm)	1	6
25			NLA - Locknut (M6 x 1)	1	2
			Instruction Sheet, Seat Rail Extensions	1	1
NOTE: A - Includes items 2 - 14. B - Includes items 3 - 5. C - Includes items 7 - 10. D - Includes items 12 - 14. E - Includes items 16 - 25.					

Crossbrace - Conventional Adjustable Depth (ADJDPTH)

#1 Crossbrace Assembly w/ Hardware - Adjustable Depth

Seat Width	15" W	16" W	17" W	18" W	19" W	20" W
16 - 18" D (ADJDEPTH) (SDR1618)	1530032	1530033	1530034	1530035	1530036	1530037
18 - 20" D (ADJDPTH) (SDR1820)	1530038	1530032	1530053	1530054	1530055	1530056

Crossbrace - Conventional Fixed Depth, Heavy Duty (HDPG FXDDPTH)







Crossbrace - Conventional Fixed Depth, Heavy Duty (HDPG FXDDPTH)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Crossbrace Assembly w/ Hardware - Myon Active HD	1	1
3	B	1173766	Kit, Saddle Clamp w/ Hardware (Pair)	1	2
4			Clamp, Saddle, Crossbar	1	1
5			Screw (M6 x 45mm)	1	2
6			NLA - Locknut (M6 x 1)	1	2
7	C	1173768	Kit, Pivot Link w/ Hardware (18 - 20" W) (Pair)	1	2
8			Link, Pivot, Crossbrace (18-20" W)	1	2
9			Screw, Socket Head, Cap (M6 x 35mm)	1	2
10			NLA - Locknut (M6 x 1)	1	2
			Carton	1	1
			NLA - Bag, Poly Zip-Lip (8-3/4" L x 6" W)	1	1
7	C	1173769	Kit, Pivot Link w/ Hardware (22" W)	1	2
8			Link, Pivot, Crossbrace (22" W)	1	2
9			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
10			NLA - Locknut (M6 x 1)	1	2
7	C	1173770	Kit, Pivot Link w/ Hardware (24" W)	1	2
8			Link, Pivot, Crossbrace (24" W)	1	2
9			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
10			NLA - Locknut (M6 x 1)	1	2
11	D	1173773-NLA	NLA - Kit, Seat Rail Extensions, Myon (2")	1	1
12			NLA - Plug Button	1	4
13			Extension, Seat Rail (2")	1	2
14			NLA - Screw, Phillips Pan Head (M3.5 x .6 x 9.5mm)	1	4
NOTE: A - Includes items 3 - 14. B - Includes items 4 - 6. C - Includes items 8 - 10. D - Includes items 13 - 15.					

#1 Crossbrace Assembly w/ Hardware - Fixed Depth HD

Seat Depth	20" W	22" W	24" W
16"	1488618	1488620	N/A
17"	1169097	1169098	N/A
18"	1488619	1488621	1488622
19"	N/A	1169105	1169106
20"	1169112	1169113	1169114

Crossbrace - Conventional Adjustable Depth, Heavy Duty (HDPG ADJDPTH)



Crossbrace - Conventional Adjustable Depth, Heavy Duty (HDPG ADJDPTH)


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Crossbrace Assembly - Myon Junior	1	1
2	B	1173766	Kit, Saddle Clamp w/ Hardware (Pair)	1	1
3			Clamp, Saddle, Crossbar	1	2
4			Screw (M6 x 45mm)	1	4
5			NLA - Locknut (M6 x 1)	1	4
6	C	1173774-NLA	NLA - Kit, Plug Buttons w/ Hardware	1	1
7			NLA - Plug Button	1	4
8			NLA - Screw, Phillips Pan Head (M3.5 x .6 x 9.5mm)	1	4
9	D	1173768	Kit, Pivot Link w/ Hardware (18 - 20" W) (Pair)	1	1
10			Link, Pivot, Crossbrace (18-20" W)	1	2
11			Screw, Socket Head, Cap (M6 x 35mm)	1	2
12			NLA - Locknut (M6 x 1)	1	2
			Carton	1	1
			NLA - Bag, Poly Zip-Lip (8-3/4" L x 6" W)	1	1
9	D	1173769	Kit, Pivot Link w/ Hardware (22" W)	1	1
10			Link, Pivot, Crossbrace (22" W)	1	2
11			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12			NLA - Locknut (M6 x 1)	1	2
9	D	1173770	Kit, Pivot Link w/ Hardware (24" W)	1	1
10			Link, Pivot, Crossbrace (24" W)	1	2
11			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12			NLA - Locknut (M6 x 1)	1	2
13	E	1172274	Kit, Seat Rail Extensions (16 - 18" D)	1	1
14			Extension, Seat Rail (2")	1	2
15			Extension, Seat Rail (1")	1	2
16			NLA - Plug Button	1	2
17			Cap, Extension, Seat Rail (1")	1	2
18			Rod, Upholstery (2")	1	2
19			Rod, Upholstery (1")	1	2
20			Guide, Seat	1	2
21			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
22			NLA - Screw, Phillips Pan Head (M3.5 x .6 x 9.5mm)	1	6
23			NLA - Locknut (M6 x 1)	1	2
			Instruction Sheet, Seat Rail Extensions	1	1
NOTE: A - Includes items 2 - 12. B - Includes items 3 - 5. C - Includes items 7 and 8. D - Includes items 10 - 12. E - Includes items 14 - 23.					

Crossbrace - Conventional Adjustable Depth, Heavy Duty (HDPG ADJDPTH)

#1 Crossbrace Assembly - Adjustable Depth HD

Seat Depth	20" Wide (D20)	22" Wide (D22)	24" Wide (D24)
16 - 18" Deep (SDR1618)	1488619	1488621	1488622
18 - 20" Deep (SDR1820)	N/A	1169113	1169114

Crossbrace - Conventional Growable Seat Depth



Crossbrace - Conventional Growable Seat Depth


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Crossbrace Assembly - Growable Conventional Adjustable Seat Depth	1	1
2	B	1173765	Kit, Saddle Clamp w/ Hardware (Pair)	1	1
3			Clamp, Saddle	1	2
4			NLA - Screw, Socket Button Head (M6 x 40mm)	1	4
5	C	1173777	NLA - Locknut (M6 x 1)	1	4
6			Kit, Crossbrace w/ Hardware (Pair)	1	1
7			Washer, Plastic	1	2
8			Bolt (M6 x 80mm)	1	2
9			Washer (6.6 x 22 x 2)	1	4
10	D	1173776	NLA - Locknut (M6 x 1)	1	2
11			Kit, Pivot Link w/ Hardware (12 - 15" W)	1	1
12			Link, Pivot, Crossbrace (12-15" W)	1	2
13			Screw (M6 x 45mm)	1	2
14	E	1172274	NLA - Locknut (M6 x 1)	1	4
15			Kit, Seat Rail Extensions (16 - 18" D)	1	1
16			Extension, Seat Rail (2")	1	2
17			Extension, Seat Rail (1")	1	2
18			NLA - Plug Button	1	2
19			Cap, Extension, Seat Rail (1")	1	2
20			Rod, Upholstery (2")	1	2
21			Rod, Upholstery (1")	1	2
22			Guide, Seat	1	2
23			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
24			NLA - Screw, Phillips Pan Head (M3.5 x .6 x 9.5mm)	1	6
25			NLA - Locknut (M6 x 1)	1	2
			Instruction Sheet, Seat Rail Extensions	1	1
NOTE: A - Includes items 2 - 14. B - Includes items 3 - 5. C - Includes items 7 - 10. D - Includes items 12 - 14. E - Includes items 16 - 25.					

#1 Growable Conventional Adjustable Seat Depth

Seat Depth	8 - 11" W	12 - 15" W
8 - 10" Deep	1511958	N/A
10 - 12" Deep	1511959	N/A
11 - 13" Deep	N/A	1512008
12 - 14" Deep	1513260	N/A
13 - 15" Deep	N/A	1512009

Crossbrace - Conventional Growable Seat Depth

Crossbrace - Conventional Fixed Seat Depth




Crossbrace - Conventional Fixed Seat Depth

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Crossbrace Assembly - Junior Fixed	1	1
2	B	1173765	Kit, Saddle Clamp w/ Hardware (Pair)	1	1
3			Clamp, Saddle	1	2
4			NLA - Screw, Socket Button Head (M6 x 40mm)	1	4
5			NLA - Locknut (M6 x 1)	1	4
6	C	1173776	Kit, Pivot Link w/ Hardware (12 - 15" W)	1	1
7			Link, Pivot, Crossbrace (12-15" W)	1	2
8			Screw (M6 x 45mm)	1	2
9			NLA - Locknut (M6 x 1)	1	4
10	D	1173777	Kit, Crossbrace w/ Hardware (Pair)	1	1
11			Washer, Plastic	1	2
12			Bolt (M6 x 80mm)	1	2
13			Washer (6.6 x 22 x 2)	1	4
14			NLA - Locknut (M6 x 1)	1	2
15	E	1173774-NLA	NLA - Kit, Plug Buttons w/ Hardware	1	1
16			NLA - Plug Button	1	4
17			NLA - Screw, Phillips Pan Head (M3.5 x .6 x 9.5mm)	1	4
NOTE: A - Includes items 2 - 17. B - Includes items 3 - 5. C - Includes items 7 - 9. D - Includes items 11 - 14. E - Includes items 16 - 17.					

#1 Crossbrace Assembly - Junior Fixed

Seat Depth	12" Wide Chair	13" Wide Chair	14" Wide Chair
12" Deep	1172267	N/A	N/A
13" Deep	1172268	N/A	N/A
14" Deep	1172269	1172270	1172272
15" Deep	N/A	1172271	1172273

Arms - Cantilever - Clamp On Desk Length (AT55D), Full Length (AT55F)



Arms - Cantilever - Clamp On Desk Length (AT55D), Full Length (AT55F)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	B	8881011810	Arm Pad, Desk Length w/ Law Label	1	1
1	B	8881055204	Pad, Arm, Desk Length	1	1
1	C	8881127509	Kit, Arm Pad w/Screws - Full Length	1	1
			Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	1	1
			NLA - Pad, Arm, Full Length w/ Law Label	1	1
2	A,D	1068656	Arm Assembly, Cantilever - Right	1	1
3			Bracket, Armrest Mounting, Black	1	1
4			Arm Support, Clamp On	1	1
5			Clamp, Half, Threaded	1	1
6			Clamp, Half, Non-Threaded	1	1
7			Tube, Arm, Extension, Black	1	1
8			Tube, Slide, Cantilever Arm, Black	1	1
9			Tube, Arm, Black	1	1
10			Spring, Compression (5/8 x 1 x 3/64")	1	1
11			Plug, Latch	1	1
12			Screw, Hex Head (1/4-28 x 2")	1	1
13			Tip, Plastic, Black	1	1
14			Lever, Actuator	1	1
15			NLA - Pin, Locking	1	1
16			Spacer, Coved (1/4 x 1 x 9/16")	1	1
17			Washer, Nylon (3/8 x 7/8 x 1/16")	1	1
18			NLA - Screw, Hex Head (1/4-28 x 3/4")	1	1
19			Screw, Socket Head (1/4-20 x 1")	1	2
20			NLA - Lockwasher, High Collar (1/4")	1	2
21			Locknut (1/4-28)	1	3
22			Plug Button, Black (1")	1	1
23			Screw, Phillips Pan Head (#10-32 x 1-3/8")	1	2
24			Washer, Nylon (1/4 x 1/2 x 1/32")	1	1
2	A,D	1068657	Arm Assembly, Cantilever - Left	1	1
3			Bracket, Armrest Mounting, Black	1	1
4			Arm Support, Clamp On	1	1
5			Clamp, Half, Threaded	1	1
6			Clamp, Half, Non-Threaded	1	1
7			Tube, Arm, Extension, Black	1	1
8			Tube, Slide, Cantilever Arm, Black	1	1
9			Tube, Arm, Black	1	1
10			Spring, Compression (5/8 x 1 x 3/64")	1	1
11			Plug, Latch	1	1
12			Screw, Hex Head (1/4-28 x 2")	1	1
13			Tip, Plastic, Black	1	1
14			Lever, Actuator	1	1
15			NLA - Pin, Locking	1	1
16			Spacer, Coved (1/4 x 1 x 9/16")	1	1

Arms - Cantilever - Clamp On Desk Length (AT55D), Full Length (AT55F)


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
17	A,E	1173055	Washer, Nylon (3/8 x 7/8 x 1/16")	1	1
18			NLA - Screw, Hex Head (1/4-28 x 3/4")	1	1
19			Screw, Socket Head (1/4-20 x 1")	1	2
20			NLA - Lockwasher, High Collar (1/4")	1	2
21			Locknut (1/4-28)	1	3
22			Plug Button, Black (1")	1	1
23			Screw, Phillips Pan Head (#10-32 x 1-3/8")	1	2
24			Washer, Nylon (1/4 x 1/2 x 1/32")	1	1
2			Arm Assembly, Cantilever, Clamp-On - Right	1	1
3			Bracket, Armrest Mounting, Black	1	1
4			NLA - Strap, Wheel Lock	1	1
5			Clamp, Half, Threaded	1	1
6			Clamp, Half, Non-Threaded	1	1
7			Tube, Arm, Extension, Black	1	1
8			Tube, Slide, Cantilever Arm, Black	1	1
9			Tube, Arm, Black	1	1
10			Spring, Compression (5/8 x 1 x 3/64")	1	1
11			Plug, Latch	1	1
12			Screw, Hex Head (1/4-28 x 2")	1	1
13			Tip, Plastic, Black	1	1
14			Lever, Actuator	1	1
15			NLA - Pin, Locking	1	1
16			Spacer, Coved (1/4 x 1 x 9/16")	1	1
17			Washer, Nylon (3/8 x 7/8 x 1/16")	1	1
18			NLA - Screw, Hex Head (1/4-28 x 3/4")	1	1
19			Screw, Socket Head (1/4-20 x 1")	1	2
20			NLA - Lockwasher, High Collar (1/4")	1	2
21			Locknut (1/4-28)	1	3
22			Plug Button, Black (1")	1	1
23			Screw, Phillips Pan Head (#10-32 x 1-3/8")	1	1
24			Washer, Nylon (1/4 x 1/2 x 1/32")	1	1
2	A,E	1173056	Arm Assembly, Cantilever, Clamp-On - Left	1	1
3			Bracket, Armrest Mounting, Black	1	1
4			NLA - Strap, Wheel Lock	1	1
5			Clamp, Half, Threaded	1	1
6			Clamp, Half, Non-Threaded	1	1
7			Tube, Arm, Extension, Black	1	1
8			Tube, Slide, Cantilever Arm, Black	1	1
9			Tube, Arm, Black	1	1
10			Spring, Compression (5/8 x 1 x 3/64")	1	1
11			Plug, Latch	1	1
12			Screw, Hex Head (1/4-28 x 2")	1	1
13			Tip, Plastic, Black	1	1

Arms - Cantilever - Clamp On Desk Length (AT55D), Full Length (AT55F)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
14			Lever, Actuator	1	1
15			NLA - Pin, Locking	1	1
16			Spacer, Coved (1/4 x 1 x 9/16")	1	1
17			Washer, Nylon (3/8 x 7/8 x 1/16")	1	1
18			NLA - Screw, Hex Head (1/4-28 x 3/4")	1	1
19			Screw, Socket Head (1/4-20 x 1")	1	2
20			NLA - Lockwasher, High Collar (1/4")	1	2
21			Locknut (1/4-28)	1	3
22			Plug Button, Black (1")	1	1
23			Screw, Phillips Pan Head (#10-32 x 1-3/8")	1	2
24			Washer, Nylon (1/4 x 1/2 x 1/32")	1	1
NOTE: A - Includes items 2 - 24. B - For use with AT55D. C - For use with AT55F. D - For use with ABH1418, ABH1519, ABH1620, and 8BEND Back Cane setups. E - For use with FBH14, FBH16, and FBH20 Back Cane setups.					

Arms - Conventional Adjustable Height

Desk Length (CAD), Full Length (CAF)



Arms - Conventional Adjustable Height Desk Length (CAD), Full Length (CAF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881169069	Armrest Assembly, Swingaway, Adjustable Height, Desk Length, Color - Right	1	1
3			Pad, Arm, Desk Length	1	1
5			Arm Pad Mounting Assembly - Right	1	1
6			Guard, Clothing	1	1
7			Tube, Armrest, Color - Right	1	1
8			Support, Armrest - Right	1	1
9			Lever, Release	1	1
10			Button, Spring	1	1
11			Plug, Armrest	1	1
12			Support, Armrest, Front	1	1
13			Screw (M6 x 45mm)	1	2
14			NLA - Locknut (M6 x 1)	1	2
15			Washer (6.6 x 22 x 2)	1	1
16			Bushing, Black	1	1
17			Screw, Countersunk (M4.8)	1	1
18			Screw, Socket Head (M12 x 12.5)	1	1
19			Bushing, Armrest	1	1
20			Nut, Square, Armrest	1	1
1	A	8881169070	Armrest Assembly, Swingaway, Adjustable Height, Desk Length, Color - Left	1	1
3			Pad, Arm, Desk Length	1	1
5			Arm Pad Mounting Assembly - Left	1	1
6			Guard, Clothing	1	1
7			Tube, Armrest, Color - Left	1	1
8			Support, Armrest - Left	1	1
9			Lever, Release	1	1
10			Button, Spring	1	1
11			Plug, Armrest	1	1
12			Support, Armrest, Front	1	1
13			Screw (M6 x 45mm)	1	2
14			NLA - Locknut (M6 x 1)	1	2
15			Washer (6.6 x 22 x 2)	1	1
16			Bushing, Black	1	1
17			Screw, Countersunk (M4.8)	1	2
18			Screw, Socket Head (M12 x 12.5)	1	1
19			Bushing, Armrest	1	1
20			Nut, Square, Armrest	1	1
2	B	8881169071	Armrest Assembly, Swingaway, Adjustable Height, Full Length, Color - Right	1	1
4			Arm Pad, Full Length	1	1
5			Arm Pad Mounting Assembly - Right	1	1
6			Guard, Clothing	1	1
7			Tube, Armrest, Color - Right	1	1

Arms - Conventional Adjustable Height Desk Length (CAD), Full Length (CAF)






Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
8			Support, Armrest - Right	1	1
9			Lever, Release	1	1
10			Button, Spring	1	1
11			Plug, Armrest	1	1
12			Support, Armrest, Front	1	1
13			Screw (M6 x 45mm)	1	2
14			NLA - Locknut (M6 x 1)	1	2
15			Washer (6.6 x 22 x 2)	1	1
16			Bushing, Black	1	1
17			Screw, Countersunk (M4.8)	1	3
18			Screw, Socket Head (M12 x 12.5)	1	1
19			Bushing, Armrest	1	1
20			Nut, Square, Armrest	1	1
2	B	8881169072	Armrest Assembly, Swingaway, Adjustable Height, Full Length, Color - Left	1	1
4			Arm Pad, Full Length	1	1
5			Arm Pad Mounting Assembly - Left	1	1
6			Guard, Clothing	1	1
7			Tube, Armrest, Color - Left	1	1
8			Support, Armrest - Left	1	1
9			Lever, Release	1	1
10			Button, Spring	1	1
11			Plug, Armrest	1	1
12			Support, Armrest, Front	1	1
13			Screw (M6 x 45mm)	1	2
14			NLA - Locknut (M6 x 1)	1	2
15			Washer (6.6 x 22 x 2)	1	1
16			Bushing, Black	1	1
17			Screw, Countersunk (M4.8)	1	3
18			Screw, Socket Head (M12 x 12.5)	1	1
19			Bushing, Armrest	1	1
20			Nut, Square, Armrest	1	1
21	C	1151395	Arm Pad w/Hardware, Desk Length	1	1
3			Pad, Arm, Desk Length	1	1
17			Screw, Phillips Flat Head (#4-40 x 3/8")	1	2
22	D	1151422	Arm Pad w/Hardware, Full Length	1	1
4			Pad, Arm, Full Length	1	1
17			Screw, Phillips Flat Head (#4-40 x 3/8")	1	3
23	E	1430423	Kit, Swingaway Armrest Support - Right	1	1
8			Support, Armrest - Right	1	1
13			Screw (M6 x 45mm)	1	1
14			NLA - Locknut (M6 x 1)	1	1
15			Washer (6.6 x 22 x 2)	1	1

Arms - Conventional Adjustable Height Desk Length (CAD), Full Length (CAF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
16	E	1431284	Bushing, Black	1	1
17			Screw, Countersunk (M4.8)	1	1
23			Kit, Swingaway Armrest Support - Left	1	1
8			Support, Armrest - Left	1	1
13			Screw (M6 x 45mm)	1	1
14			NLA - Locknut (M6 x 1)	1	1
15			Washer (6.6 x 22 x 2)	1	1
16			Bushing, Black	1	1
17			Screw, Countersunk (M4.8)	1	1
NOTE: A - Includes items 3, 5, and 6 - 20. B - Includes items 4 - 20. C - Kit includes items 3 and 17. D - Kit includes items 4 and 17. E - Kit includes items 8 and 13 - 17.					

Arms - Conventional Fixed Height

Desk Length (CNAD), Full Length (CNAF)



Arms - Conventional Fixed Height Desk Length (CNAD), Full Length (CNAF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881169065	Armrest Assembly, Swingaway, Fixed Height, Desk Length, Color - Right	1	1
3			Arm Pad, Desk Length	1	1
6			Guard, Clothing, Desk Length - Right	1	1
7			Lever, Release	1	1
8			Button, Spring	1	1
9			Plug, Armrest	1	1
10			Support, Armrest, Front	1	1
11			Screw (M6 x 45mm)	1	2
12			NLA - Locknut (M6 x 1)	1	2
13			Support, Armrest - Right	1	1
14			Washer (6.6 x 22 x 2)	1	1
15			Bushing, Black	1	1
16			Screw, Countersunk (M4.8)	1	1
18			Screw, Socket, Button Head (M5 x 35mm)	1	3
19			Grip, Foam Armrest	1	1
20			Tube, Armrest, Color	1	1
1	A	8881169066	Armrest Assembly, Swingaway, Fixed Height, Desk Length, Color - Left	1	1
3			Arm Pad, Desk Length	1	1
6			Guard, Clothing, Desk Length - Left	1	1
7			Lever, Release	1	1
8			Button, Spring	1	1
9			Plug, Armrest	1	1
10			Support, Armrest, Front	1	1
11			Screw (M6 x 45mm)	1	2
12			NLA - Locknut (M6 x 1)	1	2
13			Support, Armrest - Left	1	1
14			Washer (6.6 x 22 x 2)	1	1
15			Bushing, Black	1	1
16			Screw, Countersunk (M4.8)	1	1
18			Screw, Socket, Button Head (M5 x 35mm)	1	3
19			Grip, Foam Armrest	1	1
20			Tube, Armrest, Color	1	1
2	B	8881169067	Armrest Assembly, Swingaway, Fixed Height, Full Length, Color - Right	1	1
4			Arm Pad, Full Length	1	1
5			Tube, Support, Armpad, Full Length	1	1
6			Guard, Clothing, Desk Length - Right	1	1
7			Lever, Release	1	1
8			Button, Spring	1	1
9			Plug, Armrest	1	1
10			Support, Armrest, Front	1	1
11			Screw (M6 x 45mm)	1	2


Arms - Conventional Fixed Height Desk Length (CNAD), Full Length (CNAF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
12			NLA - Locknut (M6 x 1)	1	2
13			Support, Armrest - Right	1	1
14			Washer (6.6 x 22 x 2)	1	1
15			Bushing, Black	1	1
16			Screw, Countersunk (M4.8)	1	1
17			Button, Plug, Reflector (22mm)	1	1
18			Screw, Socket, Button Head (M5 x 35mm)	1	4
19			Grip, Foam Armrest	1	1
20			Tube, Armrest, Color	1	1
2	B	8881169068	Armrest Assembly, Swingaway, Fixed Height, Full Length, Color - Left	1	1
4			Arm Pad, Full Length	1	1
5			Tube, Support, Armpad, Full Length	1	1
6			Guard, Clothing, Desk Length - Left	1	1
7			Lever, Release	1	1
8			Button, Spring	1	1
9			Plug, Armrest	1	1
10			Support, Armrest, Front	1	1
11			Screw (M6 x 45mm)	1	2
12			NLA - Locknut (M6 x 1)	1	2
13			Support, Armrest - Left	1	1
14			Washer (6.6 x 22 x 2)	1	1
15			Bushing, Black	1	1
16			Screw, Countersunk (M4.8)	1	1
17			Button, Plug, Reflector (22mm)	1	1
18			Screw, Socket, Button Head (M5 x 35mm)	1	4
19			Grip, Foam Armrest	1	1
20			Tube, Armrest, Color	1	1
21	C	1182182	Kit, Armpad, Full Length (Pair) (CNAF, TAF)	1	1
3			Arm Pad, Full Length	1	2
18			Screw, Socket, Button Head (M5 x 35mm)	1	8
21	C	1162413	Kit, Armpad, Desk Length (Pair) (CNAD, TAD)	1	1
3			Arm Pad, Desk Length	1	2
18			Screw, Socket, Button Head (M5 x 35mm)	1	8
22	D	1432595	Kit, Desk to Full Armpad w/ Tube and Hardware	1	1
4			Arm Pad, Full Length	1	1
5			Tube, Support, Armpad, Full Length	1	1
17			Button, Plug, Reflector (22mm)	1	1
18			Screw, Socket, Button Head (M5 x 35mm)	1	4
23	E	1430423	Kit, Swingaway Armrest Support - Right	1	1
11			Support, Armrest - Right	1	1
12			Screw (M6 x 45mm)	1	1
13			NLA - Locknut (M6 x 1)	1	1

Arms - Conventional Fixed Height Desk Length (CNAD), Full Length (CNAF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
14	E	1431284	Washer (6.6 x 22 x 2)	1	1
15			Bushing, Black	1	1
16			Screw, Countersunk (M4.8)	1	1
23			Kit, Swingaway Armrest Support - Left	1	1
11			Support, Armrest - Left	1	1
12			Screw (M6 x 45mm)	1	1
13			NLA - Locknut (M6 x 1)	1	1
14			Washer (6.6 x 22 x 2)	1	1
15			Bushing, Black	1	1
16			Screw, Countersunk (M4.8)	1	1
NOTE: A - Includes items 3, 6-16 and 18-20. B - Includes items 4-20. C - Kit includes items 3 and 18. D - Kit includes items 4, 5, 17, and 18. E - Kit includes items 11 - 16.					


Arms - Swing Away - Front Mount - Transport Brackets



Arms - Swing Away - Front Mount - Transport Brackets

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1190148	Kit, Swingaway Arm Mounting w/ Hardware	1	1
2			Socket, Arm	1	2
3			Arm, Swing Away, Half, Padded - Right	1	1
3			Arm, Swing Away, Half, Padded - Left	1	1
4			Sleeve, Bearing	1	2
5			Screw (M6 x 1 x 60mm)	1	4
6			Bolt, Hex Head (M6 x 45mm)	1	2
7			NLA - Washer, Steel (1/4 x 1/2 x 3/64")	1	12
8		NLA - Locknut (M6 x 1)	1	6	
NOTE: A - Kit includes items 2 thru 8.					

Arms - Swing Away - Rear Mount - Transport Brackets




Arms - Swing Away - Rear Mount - Transport Brackets

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1173040	Kit, Swingaway Arm Mounting w/ Hardware	1	1
2			Socket, Arm	1	2
3			Arm, Swing Away, Half, Padded - Right	1	1
3			Arm, Swing Away, Half, Padded - Left	1	1
4			Sleeve, Bearing	1	2
5			Bolt, Hex Head (M6 x 45mm)	1	2
6			Screw (M6 x 1 x 60mm)	1	2
7			NLA - Washer, Steel (1/4 x 1/2 x 3/64")	1	8
8		NLA - Locknut (M6 x 1)	1	4	
NOTE: A - Kit includes items 2 thru 8.					

Arms - T Style Adjustable - Junior, Active (TAHJ, TAHD)


Myon, Myon Active (TAHD - Adjustable Height, Desk Length)

Myon Junior (TAHJ - Adjustable Height, Junior)



1
TAHD - Adjustable Height, Desk Length

TAHJ - Adjustable Height, Junior 13



SCM_061A

Arms - T Style Adjustable - Junior, Active (TAHJ, TAHD)

Myon, Myon Active (TAHD - Adjustable Height, Desk Length)

Myon Junior (TAHJ - Adjustable Height, Junior)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C,A	1489055	T-Arm Assembly, Desk Length, Junior - Right	1	1
2			Arm Pad, Desk Length	1	1
3			Support, T Arm, Junior	1	1
4			NLA - Support, Armrest	1	1
5			Guard, Clothing - Right/Left	1	1
6			Plug Button, Black (3/4")	1	2
7			Screw, Socket Button Head (M6 x 1 x 12mm)	1	2
8			Nut (M6)	1	2
9			Screw, Socket, Button Head (M5 x 35mm)	1	2
10			Screw, Socket Head (M6 x 50)	1	2
11			NLA - Locknut (M6 x 1)	1	2
12			Screw (M6 x 1 x 10mm)	1	1
1	C,A	1489056	T-Arm Assembly, Desk Length, Junior - Left	1	1
2			Arm Pad, Desk Length	1	1
3			Support, T Arm, Junior	1	1
4			NLA - Support, Armrest	1	1
5			Guard, Clothing - Right/Left	1	1
6			Plug Button, Black (3/4")	1	2
7			Screw, Socket Head (M6 x 50)	1	2
8			Nut (M6)	1	3
9			Screw, Socket, Button Head (M5 x 35mm)	1	2
10			Screw, Socket Head (M6 x 50)	1	2
11			NLA - Locknut (M6 x 1)	1	2
12			Screw (M6 x 1 x 10mm)	1	3
13	D,B	1192393	Kit, T-Arm Junior, Bear - Right	1	1
3			Arm Pad, Desk Length	1	1
5			NLA - Support, Armrest	1	1
6			Screw, Socket Button Head (M6 x 1 x 12mm)	1	2
8			Plug Button, Black (3/4")	1	2
9			Screw, Socket, Button Head (M5 x 35mm)	1	2
11			NLA - Locknut (M6 x 1)	1	2
12			Screw, Socket Head (M6 x 50)	1	2
14			Guard, Clothing, Junior	1	1
15			Support, Armrest, Junior - Right	1	1
16			Nut (M6)	1	2
17			Locking Pin Assembly	1	1
			Instruction Sheet, Service, Myon	1	1
13		1192394	Kit, T-Arm Junior, Bear - Left	1	1
3			Arm Pad, Desk Length	1	1
5			NLA - Support, Armrest	1	1


Arms - T Style Adjustable - Junior, Active (TAHJ, TAHD)

Myon, Myon Active (TAHD - Adjustable Height, Desk Length)

Myon Junior (TAHJ - Adjustable Height, Junior)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
6			Screw, Socket Button Head (M6 x 1 x 12mm)	1	2
8			Plug Button, Black (3/4")	1	2
9			Screw, Socket, Button Head (M5 x 35mm)	1	2
11			NLA - Locknut (M6 x 1)	1	2
12			Screw, Socket Head (M6 x 50)	1	2
14			Guard, Clothing, Junior	1	1
15			Support, Armrest, Junior, Left	1	1
16			Nut (M6)	1	2
17			Locking Pin Assembly	1	1
			Instruction Sheet, Service, Myon	1	1
NOTE: A - Includes items 2 - 12. B - Includes items 3, 5, 6, 8, 9, 11, 12, 14, 15, 16, and item 17. C - Used with Myon and Myon Active Desk Length Chairs (TAHD). D - Used with Myon Junior Desk Length Chairs (TAHJ).					

Arms - T Style Adjustable Desk Length (TAD), Full Length (TAF)




Arms - T Style Adjustable Desk Length (TAD), Full Length (TAF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1428596-0601	Armrest, T-Arm, Adjustable Height, Desk - Right	1	1
3			Arm Pad, Desk Length	1	1
5			Tube, Armrest - Right	1	1
6			Guard, Clothing - Right/Left	1	1
7			NLA - Support, Armrest	1	1
8			Screw, Socket, Button Head (M5 x 35mm)	1	2
9			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
10			NLA - Locknut (M6 x 1)	1	5
11			Plug Button, Black (3/4")	1	2
12			Screw, Socket Head (M6 x 50)	1	2
13			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
1	A	1428597-0601	Armrest, T-Arm, Adjustable Height, Desk - Left	1	1
3			Arm Pad, Desk Length	1	1
5			Tube, Armrest, Black - Left	1	1
6			Guard, Clothing - Right/Left	1	1
7			NLA - Support, Armrest	1	1
8			Screw, Socket, Button Head (M5 x 35mm)	1	2
10			Plug Button, Black (3/4")	1	2
11			Screw, Socket Head (M6 x 50)	1	2
12			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
13			NLA - Locknut (M6 x 1)	1	5
2	B	1428598-0601	Armrest, T-Arm, Adjustable Height, Full - Right	1	1
4			Arm Pad, Full Length	1	1
5			Tube, Armrest - Right	1	1
6			Guard, Clothing - Right/Left	1	1
7			NLA - Support, Armrest	1	1
8			Screw, Socket, Button Head (M5 x 35mm)	1	2
9			Button, Plug (22mm)	1	1
10			Plug Button, Black (3/4")	1	2
11			Screw, Socket Head (M6 x 50)	1	2
12			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
13			NLA - Locknut (M6 x 1)	1	5
14			Screw, Socket Head, Torx (K50 x 30mm)	1	1
15			Washer (.250 x .562 x .065/036)	1	2
16	B	1428599-0601	Tube, Support, Armpad, Full Length	1	1
2			Armrest, T-Arm, Adjustable Height, Full - Left	1	1
4			Arm Pad, Full Length	1	1
5			Tube, Armrest, Black - Left	1	1
6			Guard, Clothing - Right/Left	1	1
7			NLA - Support, Armrest	1	1
8			Screw, Socket, Button Head (M5 x 35mm)	1	2
9			Button, Plug (22mm)	1	1
10			Plug Button, Black (3/4")	1	2

Arms - T Style Adjustable Desk Length (TAD), Full Length (TAF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
11			Screw, Socket Head (M6 x 50)	1	2
12			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
13			NLA - Locknut (M6 x 1)	1	5
14			Screw, Socket Head, Torx (K50 x 30mm)	1	1
15			Washer (.250 x .562 x .065/036)	1	2
16			Tube, Support, Armpad, Full Length	1	1
17	C	1182182	Kit, Armpad, Full Length (Pair) (CNAF, TAF)	1	1
3			Arm Pad, Full Length	1	2
8			Screw, Socket, Button Head (M5 x 35mm)	1	8
	G	1162413	Kit, Armpad, Desk Length (Pair) (CNAD, TAD)	1	1
			Arm Pad, Desk Length	1	2
			Screw, Socket, Button Head (M5 x 35mm)	1	6
18	D	1432595	Kit, Desk to Full Armpad w/ Tube and Hardware	1	1
4			Arm Pad, Full Length	1	1
9			Button, Plug, Reflector (22mm)	1	1
14			Screw, Socket, Button Head (M5 x 35mm)	1	4
16			Tube, Support, Armpad, Full Length	1	1
20	E	1430421	Kit, Attaching Hardware, T-Arm	1	1
7			NLA - Support, Armrest	1	1
11			Screw, Socket Head (M6 x 50)	1	2
13			NLA - Locknut (M6 x 1)	1	2
19	F	1151387	Kit, Clothing Guard w/ Hardware	1	1
6			Guard, Clothing - Right/Left	1	1
12			Screw, Socket Button Head (M6 x 1 x 12mm)	1	3
13			NLA - Locknut (M6 x 1)	1	3
9		1151389	Package, Plug Button, Black (3/4")	10	1
10		1151390-NLA	NLA - Package, Plug Button, Black (22mm)	10	1
NOTE: A - Includes items 3, and 5 - 13. B - Includes items 4, 9, 14, 16 C - Kit includes items 3 and 8. D - Kit includes items 4, 9, 14, and 16. E - Kit includes items 7, 11, and 13. F - Kit includes items 6, 12, and 13. G - Kit includes Desk Length Arm Pad (2) and Screws (3). Not shown in illustration.					

Arms - T Style Adjustable - Heavy Duty Desk Length (TAD), Full Length (TAF)



Arms - T Style Adjustable - Heavy Duty

Desk Length (TAD), Full Length (TAF)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1478770	T-Arm Assembly w/ Hardware, Desk Length, Heavy Duty - Right	1	1
2			Arm Pad, Desk Length	1	1
3			Tube, Armrest - Right	1	1
4			NLA - Support, Armrest	1	1
5			Guard, Side, Clothing	1	1
6			Support, T Arm, Heavy Duty	1	1
7			Spacer, Nylon	1	1
8			Plug Button, Black (3/4")	1	2
9			Screw, Socket, Button Head (M5 x 35mm)	1	2
10			Screw, Phillips Flat Head (#4-40 x 3/8")	1	1
11			Screw (M40 x 14mm)	1	2
12			Nut, Double, Plastic	1	1
13			Screw, Socket Button Head (M6 x 53mm)	1	1
14			Screw, Socket Head (M6 x 50)	1	2
15			NLA - Locknut (M6 x 1)	1	4
1	A	1478771	T-Arm Assembly w/ Hardware, Desk Length, Heavy Duty - Left	1	1
2			Arm Pad, Desk Length	1	1
3			Tube, Armrest, Black - Left	1	1
4			NLA - Support, Armrest	1	1
5			Guard, Side, Clothing	1	1
6			Support, T Arm, Heavy Duty	1	1
7			Spacer, Nylon	1	1
8			Plug Button, Black (3/4")	1	2
9			Screw, Socket, Button Head (M5 x 35mm)	1	2
10			Screw, Phillips Flat Head (#4-40 x 3/8")	1	1
11			Screw (M40 x 14mm)	1	2
12			Nut, Double, Plastic	1	1
13			Screw, Socket Button Head (M6 x 53mm)	1	1
14			Screw, Socket Head (M6 x 50)	1	2
15			NLA - Locknut (M6 x 1)	1	4
1	B	1487439	T-Arm Assembly, Central w/ Hardware, Full Length, Heavy Duty - Left	1	1
2			Pad, Arm, Full Length	1	1
3			Arm Pad Mounting Assembly - Left	1	1
4			Tube, Mounting, Adjustable	1	1
5			NLA - Support, Armrest	1	1
6			Guard, Clothing	1	1
7			Support, T Arm, Heavy Duty	1	1
8			Bracket, Clothing Guard	1	2
9			Spacer, Nylon	1	1
10			Nut, Square, Armrest	1	1
11			Screw, Countersunk (M4.8)	1	2

Arms - T Style Adjustable - Heavy Duty

Desk Length (TAD), Full Length (TAF)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
12			Bushing, Armrest	1	1
13			Screw, Socket Head (M12 x 12.5)	1	1
14			Screw (M4.2 - 16 8.8)	1	4
15			NLA - Locknut (M6 x 1)	1	3
16			Screw, Socket Button Head (M6 x 53mm)	1	1
17			Screw, Socket Head (M6 x 50)	1	2
1	B	1487440	T-Arm Assembly, Central w/ Hardware, Full Length, Heavy Duty - Right	1	1
2			Pad, Arm, Full Length	1	1
3			Arm Pad Mounting Assembly - Right	1	1
4			Tube, Mounting, Adjustable	1	1
5			NLA - Support, Armrest	1	1
6			Guard, Clothing	1	1
7			Support, T Arm, Heavy Duty	1	1
8			Bracket, Clothing Guard	1	2
9			Spacer, Nylon	1	1
10			Nut, Square, Armrest	1	1
11			Screw, Countersunk (M4.8)	1	2
12			Bushing, Armrest	1	1
13			Screw, Socket Head (M12 x 12.5)	1	1
14			Screw (M4.2 - 16 8.8)	1	4
15			NLA - Locknut (M6 x 1)	1	3
16			Screw, Socket Button Head (M6 x 53mm)	1	1
17			Screw, Socket Head (M6 x 50)	1	2
NOTE: A - Includes items 2 - 15. B - Includes items 2 - 17.					

Side Guards, Clothing Guards, and Mud Guards






Side Guards, Clothing Guards, and Mud Guards

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1512647	Clothing Guard Assembly, Junior - Right	1	1
2			Guard, Clothing, Junior	1	1
3			NLA - Support, Armrest	1	1
4			Support, Side Guard	1	1
5			Screw, Socket Head (M6 x 50)	1	2
6			Screw, Socket Button Head (M6 x 1 x 12mm)	1	2
7			Screw (M6 x 1 x 10mm)	1	1
8			Cap, Guard Support	1	1
9			Rivet	1	1
10			Nut (M6)	1	3
11			NLA - Locknut (M6 x 1)	1	2
1	A	1513214	Clothing Guard Assembly, Junior - Left	1	1
2			Guard, Clothing, Junior	1	1
3			NLA - Support, Armrest	1	1
4			Support, Side Guard	1	1
5			Screw, Socket Head (M6 x 50)	1	2
6			Screw, Socket Button Head (M6 x 1 x 12mm)	1	2
7			Screw (M6 x 1 x 10mm)	1	1
8			Cap, Guard Support	1	1
9			Rivet	1	1
10			Nut (M6)	1	3
11			NLA - Locknut (M6 x 1)	1	2
12	B,C	1492738	Side Guard Assembly w/ Hardware, Black - Right	1	1
3			NLA - Support, Armrest	1	1
5			Screw, Socket Head (M6 x 50)	1	2
7			Screw (M6 x 1 x 10mm)	1	1
10			Nut (M6)	1	1
11			NLA - Locknut (M6 x 1)	1	2
13			Guard, Side - Right	1	1
14			Rivet (4.8mm x 12mm)	1	4
15			Plate, Reinforcing	1	1
16			Support, Side Guard	1	1
12	B,C	1492739	Side Guard Assembly w/ Hardware, Black - Left	1	1
3			NLA - Support, Armrest	1	1
5			Screw, Socket Head (M6 x 50)	1	2
7			Screw (M6 x 1 x 10mm)	1	1
10			Nut (M6)	1	1
11			NLA - Locknut (M6 x 1)	1	2
13			Guard, Side - Left	1	1
14			Rivet (4.8mm x 12mm)	1	4
15			Plate, Reinforcing	1	1
16			Support, Side Guard	1	1
12	B,D	1486041	Side Guard Assembly w/ Hardware, Forward - Right	1	1
3			NLA - Support, Armrest	1	1

Side Guards, Clothing Guards, and Mud Guards

Item Number	Note	Part Number	Description	Quantity Per			
				Package	Assembly		
5	B,D	1486042	Screw, Socket Head (M6 x 50)	1	2		
7			Screw (M6 x 1 x 10mm)	1	1		
10			Nut (M6)	1	1		
11			NLA - Locknut (M6 x 1)	1	2		
13			Guard, Side - Right	1	1		
14			Rivet (4.8mm x 12mm)	1	4		
15			Plate, Reinforcing	1	1		
16			Support, Side Guard	1	1		
12			Side Guard Assembly w/ Hardware, Forward - Left	1	1		
3			NLA - Support, Armrest	1	1		
5			Screw, Socket Head (M6 x 50)	1	2		
7			Screw (M6 x 1 x 10mm)	1	1		
10			Nut (M6)	1	1		
11			NLA - Locknut (M6 x 1)	1	2		
13			Guard, Side - Left	1	1		
14			Rivet (4.8mm x 12mm)	1	4		
15			Plate, Reinforcing	1	1		
16			Support, Side Guard	1	1		
NOTE: A - Includes items 2 - 11. Used with model RESK B - Includes items 3, 5, 7, 10, 11, and 13 - 16. C - Used for models (A1571 w/ 16-17" Seat Depth), and model REMG D - Used for models (A1571 w/ 18-20" Seat Depth)							


Rigid Side Guards



Rigid Side Guards

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1036342	Kit, Side Guard Hardware (Pair)	1	1
2		1032249	Clamp, Half, Threaded	1	1
3		1018692-NLA	NLA - Clamp, Half, Non-Threaded	1	1
4		1013971	Washer, Nylon (1/4 x 7/16 x 5/64")	1	1
5		1003545	Screw, Socket Head (1/4-20 x 1-1/4")	1	1
6		1036344	Guard, Side, Large w/o Quad Hole - Right (7-3/4" H x 12-3/8" D)	1	1
6		1036345	Guard, Side, Large w/o Quad Hole - Left (7-3/4" H x 12-3/8" D)	1	1
6		1036346	Guard, Side, Small w/o Quad Hole - Right (6-1/4" H x 11-1/4" D)	1	1
6		1036347	Guard, Side, Small w/o Quad Hole - Left (6-1/4" H x 11-1/4" D)	1	1
7		1036405	Guard, Side, Large w/ Quad Hole - Right (7-3/4" H x 12-3/8" D)	1	1
7		1036406	Guard, Side, Large w/ Quad Hole - Left (7-3/4" H x 12-3/8" D)	1	1
7		1036407	Guard, Side, Small w/ Quad Hole - Right (6-1/4" H x 11-1/4" D)	1	1
7		1036408	Guard, Side, Small w/ Quad Hole - Left (6-1/4" H x 11-1/4" D)	1	1
NOTE: A - Includes Items 2-5					

Transport (TRBKTS) Brackets




Transport (TRBKTS) Brackets

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1193532	Kit, Hardware, Mounting, Transfer Brackets, Swing Away Arms	1	1
2			Plate, Tie Down, Rear	1	2
3			Bracket, Tie Down, Front	1	2
4			H-Block, Front - Right	1	1
4			H-Block, Front - Left	1	1
5			Screw, Hex Head Cap (M6 x 1 x 60)	1	4
6			nla - Screw, Hex Head w/ Patch (M6 x 1 x 50mm)	1	4
7			NLA - Screw, Socket Head Cap (M6 x 1 x 40mm)	1	2
8			NLA - Spacer	1	2
9			Spacer (6.1 x 10 x 2.4)	1	4
10			Washer (1/4 x 1/2 x 1/32")	1	12
11	NLA - Locknut (M6 x 1)	1	6		
NOTE: A - Kit includes items 2 thru 11.					

Fixed Height w/ Push Handles (FBH - FBHF)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.



Fixed Height w/ Push Handles (FBH - FBHF)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1172973	Kit, Cane Hardware w/o Lumbar Support (Fixed or Foldable) (Pair)	1	1
3			Bracket, Cane, Rear	1	2
4			Grip, Cane	1	2
5			Button, Plug (22mm)	1	2
6			Screw, Phillips Pan Head (M4.8 X 16mm)	1	2
7			Screw, Cap, Socket Head (M6 x 35mm)	1	2
8			NLA - Locknut (M6 x 1)	1	2
2	B	1172976	Kit, Cane Bracket Attaching Hardware (Pair)	1	1
8			NLA - Locknut (M6 x 1)	1	4
9			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
10			Screw, Socket Button Head (M6 x 53mm)	1	2
11		1423408-0010	Back Cane, Fixed, Zero Degree (16" H)	1	2
11		1425358-0010	Back Cane, Fixed, Zero Degree (20" H)	1	2
11		1485917-0010	Back Cane, Fixed, 8 Degree (14" H)	1	2
11		1423411-0010	Back Cane, Fixed, 7 Degree (16" H)	1	2
11		1423412-0010	Back Cane, Fixed, 7 Degree (20" H)	1	2
11		1434492	Back Cane, Fold Down (16" H)	1	2
11		1545612	Back Cane, Fold Down (17" H)	1	2
11		1548860	Back Cane, Fold Down (20" H)	1	2
NOTE: A - Kit Includes items 3 - 8. B - Kit Includes items 8, 9, and 10.					

Back Canes - Fixed Height w/ Push Handles (FBH-FBHF)

Back Cane Style	Item #1 Hardware Part No.	Item #2 Back Cane Attaching Hardware Part No.	Item #11 Back Cane Part No.
FBH14 (8 Degree)	1172973	1172976	1485917-0010
FBH16 (Straight)	1172973	1172976	1423408-0010
FBH16 (8 Degree)	1172973	1172976	1423411-0010
FBH20 (Straight)	1172973	1172976	1425358-0010
FBH20 (8 Degree)	1172973	1172976	1423412-0010
FBH16F (Fold Down)	1172973	1172976	1434492
FBH17F (Fold Down)	1172973	1172976	1545612
FBH20F (Fold Down)	1172973	1172976	1548860

Fixed Height w/ Push Handles, HD (FBH - FBHF)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.



Fixed Height w/ Push Handles, HD (FBH - FBHF)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1172974	Kit, Cane Hardware w/ Lumbar Support (Fixed or Foldable) (Pair)	1	1
2			Bracket, Cane, Rear	1	2
3			Screw, Cap, Socket Head (M6 x 35mm)	1	2
4			NLA - Locknut (M6 x 1)	1	2
5			Grip, Cane	1	2
6			Lever, Indexable	1	2
7			Insert, Rigidizer	1	2
8			Screw, Phillips Pan Head (M4.8 X 16mm)	1	2
9	B	1172976	Kit, Cane Bracket Attaching Hardware (Pair)	1	1
4			NLA - Locknut (M6 x 1)	1	4
10			Screw, Socket Button Head (M6 x 53mm)	1	2
11		1423408-0010	NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12			Back Cane, Fixed, Zero Degree (16" H)	1	2
12			Back Cane, Fixed, Zero Degree (20" H)	1	2
12			Back Cane, Fold Down (16" H)	1	2
12			Back Cane, Fold Down (17" H)	1	2
12			Back Cane, Fold Down (20" H)	1	2
12			Back Cane, Fixed, 8 Degree (14" H)	1	2
12			Back Cane, Fixed, 7 Degree (16" H)	1	2
12		1423412-0010	Back Cane, Fixed, 7 Degree (20" H)	1	2


NOTE: A - Kit includes items 2 - 8.
B - Kit includes items 4, 10, and 11.

Back Canes - Fixed Height w/ Push Handles, HD (FBH-FBHF)

Back Cane Style	Item #1 Hardware Part No.	Item #9 Back Cane Attaching Hardware Part No.	Item #12 Back Cane Part No.
FBH14 (8 Degree)	1172974	1172976	1485917-0010
FBH16 (Straight)	1172974	1172976	1423408-0010
FBH16 (8 Degree)	1172974	1172976	1423411-0010
FBH20 (Straight)	1172974	1172976	1425358-0010
FBH20 (8 Degree)	1172974	1172976	1423412-0010
FBH16F (Straight)	1172974	1172976	1434492
FBH17F (Folding)	1172974	1172976	1545612
FBH20F (Folding)	1172974	1172976	1548860

Adjustable Height w/ Push Handles (ABHPH)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.



Adjustable Height w/ Push Handles (ABHPH)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1172972	Kit, Cane Hardware w/o Lumbar Support (Adjustable Height) (Pair)	1	1
2			Bracket, Cane, Rear	1	2
3			Bushing, Flanged	1	2
4			Grip, Hand, Foam (4")	1	2
5			Insert, Backbrace	1	2
6			Screw, Cap, Socket Head (M6 x 35mm)	1	2
7			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
8			Screw, Phillips Pan Head (M4.8 X 16mm)	1	2
9			Panel Fastener, Snap In	1	2
10			NLA - Locknut (M6 x 1)	1	4
11	B	1172976	Kit, Cane Bracket Attaching Hardware (Pair)	1	1
10			NLA - Locknut (M6 x 1)	1	4
12			Screw, Socket Button Head (M6 x 53mm)	1	2
13			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
14		1169322	Insert, Back Cane, Adjustable Height	1	2
15		1172979	Back Cane, Adjustable Height, Straight (12 - 16" H)	1	2
15		1169323	Back Cane, Adjustable Height, Straight (15 - 19" H)	1	2
15		1172992	Back Cane, Adjustable Height, 8 Degree (12 - 16" H)	1	2
15		1169324	Back Cane, Adjustable Height, 8 Degree (15 - 19" H)	1	2


NOTE: A - Kit includes items 2 - 10.
B - Kit includes items 10, 12, and 13.

Back Canes - Adjustable Height w/ Push Handles (ABHPH)

Back Cane Style	Item #1 Hardware Part No.	Item #11 Cane Attaching Hardware Part No.	Item #15 Back Cane Part No.
ABH1216PH (8 Degree)	1172972	1172976	1172992
ABH1216PH (Straight)	1172972	1172976	1172979
ABH1519PH (8 Degree)	1172972	1172976	1169324
ABH1519PH (Straight)	1172972	1172976	1169323
ABH1620PH (8 Degree)	1172972	1172976	N/A
ABH1620PH (Straight)	1172972	1172976	1172981

Adjustable Height w/o Push Handles (ABH)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.



Adjustable Height w/o Push Handles (ABH)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1172977	Kit, Cane Hardware w/o Push Handle (Adjustable Height)	1	1
2			Bracket, Cane, Rear	1	2
3			Bushing, Flanged	1	2
4			Button, Plug (22mm)	1	2
5			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
6			Screw, Cap, Socket Head (M6 x 35mm)	1	2
7			NLA - Locknut (M6 x 1)	1	2
8	B	1172976	Kit, Cane Bracket Attaching Hardware (Pair)	1	1
7			NLA - Locknut (M6 x 1)	1	4
9			Screw, Socket Button Head (M6 x 53mm)	1	2
10			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
11		1169322	Insert, Back Cane, Adjustable Height	1	2
12		1172983	Back Cane, Adjustable Height, Straight (12 - 16")	1	2
12		1172985	Back Cane, Adjustable Height, Straight (15 - 19")	1	2
12		1172996	Back Cane, Adjustable Height, 8 Degree (12 - 16")	1	2
12		1172998	Back Cane, Adjustable Height, 8 Degree (15 - 19")	1	2
NOTE: A - Kit includes items 2 - 7. B - Kit includes items 7, 9, and 10.					

Back Canes - Adjustable Height w/o Push Handles (ABH)

Back Cane Style	Item #1 Hardware Part No.	Item #8 Cane Attaching Hardware Part No.	Item #12 Back Cane Part No.
ABH1216 (8 Degree)	1172977	1172976	1172996
ABH1216 (Straight)	1172977	1172976	1172983
ABH1519 (8 Degree)	1172977	1172976	1172998
ABH1519 (Straight)	1172977	1172976	1172985
ABH1620 (8 Degree)	1172977	1172976	N/A
ABH1620 (Straight)	1172977	1172976	N/A

Adjustable Angle Back (AABH - Straight)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.



SCM_443A

Adjustable Angle Back (AABH - Straight)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1172988	Kit, Adjustable Angle Back Canes w/ Hardware (12 - 16") (Pair)	1	1
2			Back Cane, Adjustable Height, Straight (12 - 16" H)	1	2
3			Insert, Back Cane, Adjustable Height (14 - 18" H)	1	2
4			Plate, Adjustable Angle, Inner - Right	1	1
4			Plate, Adjustable Angle, Inner - Left	1	1
5			Plate, Adjustable Angle, Outer - Right	1	1
5			Plate, Adjustable Angle, Outer - Left	1	1
6			Spacer, Coved (1/4 x 1 x 13/32")	1	8
7			Support, Seat Tube, Rear	1	2
8			Backrest, Excenter	1	2
9			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
10			NLA - Screw, Socket Head Cap (M6 x 1 x 40mm)	1	4
11			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12			Screw, Socket Button Head (M6 x 53mm)	1	2
13			Washer (1/4 x 9/16 x 1/16")	1	2
14			NLA - Locknut (M6 x 1)	1	8
1	A	1172989	Kit, Adjustable Angle Back Canes w/ Hardware (15 - 19") (Pair)	1	1
2			Back Cane, Adjustable Height, Straight (15 - 19" H)	1	2
3			Insert, Back Cane, Adjustable Height (14 - 18" H)	1	2
4			Plate, Adjustable Angle, Inner - Right	1	1
4			Plate, Adjustable Angle, Inner - Left	1	1
5			Plate, Adjustable Angle, Outer - Right	1	1
5			Plate, Adjustable Angle, Outer - Left	1	1
6			Spacer, Coved (1/4 x 1 x 13/32")	1	8
7			Support, Seat Tube, Rear	1	2
8			Backrest, Excenter	1	2
9			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
10			NLA - Screw, Socket Head Cap (M6 x 1 x 40mm)	1	4
11			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12			Screw, Socket Button Head (M6 x 53mm)	1	2
13			Washer (1/4 x 9/16 x 1/16")	1	2
14			NLA - Locknut (M6 x 1)	1	8
2			Back Cane, Adjustable Height, Straight (16 - 20" H)	1	2
3			Insert, Back Cane, Adjustable Height (14 - 18" H)	1	2
4			Plate, Adjustable Angle, Inner - Right	1	1
4			Plate, Adjustable Angle, Inner - Left	1	1
5			Plate, Adjustable Angle, Outer - Right	1	1
5			Plate, Adjustable Angle, Outer - Left	1	1
6			Spacer, Coved (1/4 x 1 x 13/32")	1	8
7			Support, Seat Tube, Rear	1	2


Adjustable Angle Back (AABH - Straight)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
8			Backrest, Excenter	1	2
9			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
10			NLA - Screw, Socket Head Cap (M6 x 1 x 40mm)	1	4
11			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12			Screw, Socket Button Head (M6 x 53mm)	1	2
13			Washer (1/4 x 9/16 x 1/16")	1	2
14			NLA - Locknut (M6 x 1)	1	8
15	B	1172976	Kit, Cane Bracket Attaching Hardware (Pair)	1	1
11			Screw, Socket Button Head (M6 x 53mm)	1	2
12			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
14			NLA - Locknut (M6 x 1)	1	4
NOTE: A - Kit includes items 2 - 14. B - Kit includes items 11, 12, and 14.					

Adjustable Angle Back (AABH - 8 Degree Bend)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.




Adjustable Angle Back (AABH - 8 Degree Bend)

For Back Brace and Indexable Lever, see category titled "Upholstery - Back Brace and Hardware" page.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1173001	Kit, Adjustable Angle Back Canes w/ Hardware, 8 Degree (12 - 16") (Pair)	1	1
2			Back Cane, Adjustable Height, 8 Degree (12 - 16" H)	1	2
3			Insert, Back Cane, Adjustable Height (14 - 18" H)	1	2
4			Plate, Adjustable Angle, Inner - Right	1	1
4			Plate, Adjustable Angle, Inner - Left	1	1
5			Plate, Adjustable Angle, Outer - Right	1	1
5			Plate, Adjustable Angle, Outer - Left	1	1
6			Spacer, Coved (1/4 x 1 x 13/32")	1	8
7			Support, Seat Tube, Rear	1	2
8			Backrest, Excenter	1	2
9			Screw, Socket Head Cap (M6 x 1 x 30mm)	1	2
10			NLA - Screw, Socket Head Cap (M6 x 1 x 40mm)	1	4
11			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12			Screw, Socket Button Head (M6 x 53mm)	1	2
13			Washer (1/4 x 9/16 x 1/16")	1	2
14	A	1169695	NLA - Locknut (M6 x 1)	1	8
1			Kit, Adjustable Angle Back Canes w/ Hardware, 8 Degree (15 - 19") (Pair)	1	1
2			Back Cane, Adjustable Height, 8 Degree (15 - 19" H)	1	2
3			Insert, Back Cane, Adjustable Height (14 - 18" H)	1	2
4			Plate, Adjustable Angle, Inner - Right	1	1
4			Plate, Adjustable Angle, Inner - Left	1	1
5			Plate, Adjustable Angle, Outer - Right	1	1
5			Plate, Adjustable Angle, Outer - Left	1	1
6			Spacer, Coved (1/4 x 1 x 13/32")	1	8
7			Support, Seat Tube, Rear	1	2
8			Backrest, Excenter	1	2
9			Screw, Socket Head Cap (M6 x 1 x 30)	1	2
10			NLA - Screw, Socket Head Cap (M6 x 1 x 40mm)	1	4
11			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
12			Screw, Socket Button Head (M6 x 53mm)	1	2
13			Washer (1/4 x 9/16 x 1/16")	1	2
14	B	1172976	NLA - Locknut (M6 x 1)	1	8
15			Kit, Cane Bracket Attaching Hardware (Pair)	1	1
11			Screw, Socket Button Head (M6 x 53mm)	1	2
12			NLA - Screw, Socket Button Head (M6 x 40mm)	1	2
14			NLA - Locknut (M6 x 1)	1	4

NOTE: A - Kit includes items 2 - 14.
B - Kit includes items 11, 12, and 14.

Upholstery - Back Brace and Hardware



Upholstery - Back Brace and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1170661	Kit, Hardware, Back Canes, Adjustable Angle w/ Backbrace	1	1
2			Grip, Hand, Foam (4")	1	1
3			Insert, Backbrace	1	1
4			Lever, Indexable	1	1
5			Screw, Phillips Pan Head (M4.8 X 16mm)	1	1
6			Panel Fastener, Snap In	1	1
7		See Chart	#7. Back Brace	1	1

#7 Back Brace


Back Widths	Back Cane Style (Back Brace and Indexable Lever comes standard on model picks listed below no matter what Back Upholstery is selected).	ADJTENS OPTION Selected (Back Brace and Indexable Lever will be added on Back Cane Model Picks listed below.)	Part Number
15" W	HDPG AABH1216 AABH1519 AABH1620	FBH/FBHF, HD FBH/FBHF, ABHPH, ABH, AABH-Straight, AABH-8 Degree Bend	1586520-0601
16" W	HDPG AABH1216 AABH1519 AABH1620	FBH/FBHF, HD FBH/FBHF, ABHPH, ABH, AABH-Straight, AABH-8 Degree Bend	1586521-0601
17" W	HDPG AABH1216 AABH1519 AABH1620	FBH/FBHF, HD FBH/FBHF, ABHPH, ABH, AABH-Straight, AABH-8 Degree Bend	1586522-0601
18" W	HDPG AABH1216 AABH1519 AABH1620	FBH/FBHF, HD FBH/FBHF, ABHPH, ABH, AABH-Straight, AABH-8 Degree Bend	1586523-0601

Upholstery - Back Brace and Hardware


#7 Back Brace

Back Widths	Back Cane Style (Back Brace and Indexable Lever comes standard on model picks listed below no matter what Back Upholstery is selected).	ADJTENS OPTION Selected (Back Brace and Indexable Lever will be added on Back Cane Model Picks listed below.)	Part Number
19" W	HDPG AABH1216 AABH1519 AABH1620	FBH/FBHF, HD FBH/FBHF, ABHPH, ABH, AABH-Straight, AABH-8 Degree Bend	1586524-0601
20" W	HDPG AABH1216 AABH1519 AABH1620	FBH/FBHF, HD FBH/FBHF, ABHPH, ABH, AABH-Straight, AABH-8 Degree Bend	1586525-0601
22" W	HDPG AABH1216 AABH1519 AABH1620	FBH/FBHF, HD FBH/FBHF, ABHPH, ABH, AABH-Straight, AABH-8 Degree Bend	1586526-0601
24" W	HDPG AABH1216 AABH1519 AABH1620	FBH/FBHF, HD FBH/FBHF, ABHPH, ABH, AABH-Straight, AABH-8 Degree Bend	1586527-0601


Upholstery - Back Covers, Tensioners




Adult Back Cover
(Fixed and Adjustable)



**Junior Adjustable Tension
Back Cover (ADJWHB)**



Slip-On Tensioner
(SLIPBACK)



Adjustable Tensioner
(ADJTENS)

SCM_043A

Upholstery - Back Covers, Tensioners

#1. Adult Back Cover (Fixed and Adjustable)

	Back Cane Style	Back Cane Style	Back Cane Style
Back Widths	FBH16 FBH17F ABH1519PH ABH1519 AABH1519	FBH20 FBH20F ABH1620PH ABH1620 AABH1620	FBH14 ABH1216PH ABH1216 AABH1216
15" W	1168653	1169724	1168645
16" W	1168654	1169725	1168646
17" W	1168655	1169726	1168647
18" W	1168656	1169727	1168648
19" W	1168657	1169728	1168649
20" W	1168658	1169729	1168650
22" W	1168659	1169730	1168651
24" W	1168660	1169731	1168652

Upholstery - Back Covers, Tensioners

#2. Slip On Back (SLIPBACK) Upholstery

	Back Cane Style	Back Cane Style	Back Cane Style
Seat Widths	FBH14 ABH1216PH ABH1216 AAHB1216	FBH16 FBH20 FBH16F FBH17F FBH20F	ABH1519PH ABH1519 AABH1519 ABH1620PH ABH1620 AABH1620
15" W	1532247	1532100	1532100
16" W	1532248	1532101	1532101
17" W	1532249	1532102	1532102
18" W	1532250	1532103	1532103
19" W	1532251	1532104	1532104
20" W	1532252	1532105	1532105
22" W	1532253	1532245	1532245
24" W	1532254	1532246	1532246

#3. Junior Adjustable Upholstery Cover (ADJWHB)

	Cover Type
Back Widths	ADJWHB
8 - 11" W	1170352
12 - 15" W	1170353
Adjustable Back Cuff (2 ea.)	1506257
	See Note(s) A

NOTE: A - Not shown on illustration.


When ordering Adjustable Upholstery Cover (ADJWHB) you must also order a Adjustable Tensioner (ADJTENS).

Upholstery - Back Covers, Tensioners

#4. Adjustable Tensioner (ADJTENS)

Back Widths	Part Number
15 - 18" W	1498000
19 - 24" W	1498022

Upholstery - Seat - Myon Standard (FXDDPTH, ADJDPTH)



Seat Upholstery - Fixed Depth (FXDDPTH), Sail Cloth (SAILCL)

Chair Depth	15" W	16" W	17" W	18" W	19" W	20" W
15" D	1172904	1172905	1172906	1172907	1172908	1172909
16" D	1168739	1168740	1168741	1168742	1168743	1168744
17" D	1172910	1172911	1172912	1172913	1172914	1172915
18" D	1169020	1169021	1169022	1169023	1169024	1169025
19" D	1172916	1172917	1172918	1172919	1172920	1172921
20" D	N/A	N/A	1169099	1169100	1169101	1169102
Seat Upholstery Assembly includes Item 1 (Seat Upholstery) and Item 2 (Rods)						

Upholstery - Seat - Myon Standard (FXDDPTH, ADJDPTH)

Seat Upholstery - Adjustable Depth (ADJDPTH), Sail Cloth (SAILCL)

Chair Depth	15" W	16" W	17" W	18" W	19" W	20" W
16" - 18" D	1168627	1168628	1436008	1436009	1436010	1168629
18" - 20" D	1172962	1172963	1172964	1172965	1172966	1172967


Seat Upholstery - Heavy Duty - Fixed Depth (FXDDPTH), Sail Cloth (SAILCL)

Chair Depth	20" W (HDPG)	22" W (HDPG)	24" W (HDPG)
16" D	1168744	1172922	N/A
17" D	1169019	1172923	N/A
18" D	1169025	1172924	1172925
19" D	N/A	1169116	1169117
20" D	1169102	1169103	1169104

Seat Upholstery - Heavy Duty - Adjustable Depth (ADJDPTH), Sail Cloth (SAILCL)

Chair Depth	20" W (HDPG)	22" W (HDPG)	24" W (HDPG)
16" - 18" D	1168629	1168632	1168633
18" - 20" D	N/A	1168638	1168639

Upholstery - Seat - Junior (Fixed and Adjustable)



Upholstery - Seat - Myon Junior - Sail Cloth (SAILCL)

Chair Depth	12" W	13" W	14" W
12" D	1172897	N/A	N/A
13" D	1172898	N/A	N/A
14" D	1172899	1172900	1172901
15" D	N/A	1172902	1172904


Seat Upholstery Assembly includes Item 1 (Seat Upholstery) and Item 2 (Rods).

Upholstery - Seat - Junior (Fixed and Adjustable)

Upholstery - Seat - Myon Junior - Growable (ADJWDS), Sail Cloth (SAILCL)

Chair Depth	8" - 11" W	12" - 15" W
8" - 10" D	1506248	N/A
10" - 12" D	1508441	N/A
12" - 14" D	1168637	N/A
11" - 13" D	N/A	1506256
13" - 15" D	N/A	1508442
Seat Upholstery Assembly includes Item 1 (Seat Upholstery) and Item 2 (Rods).		

Front Rigging - Swing Out Hangers 70, 80, and 90 Degree




SCM_037A

Front Rigging - Swing Out Hangers 70, 80, and 90 Degree


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881173077	Hanger, Front Rigging, 70 Degree, Color	1	1
2			Tube, Rigging, 70 Degree, Color	1	1
3			Housing, Locking Pin	1	1
4			Support, Hanger, Footrest	1	1
5			Bracket, Front Rigging	1	1
6			Calf Support Strap	1	1
7			Screw, Socket Head Cone PT Set (M10 x 1.5 x 16)	1	1
8			Screw, Socket Head, Cap (M6 x 35mm)	1	1
9			Screw, Countersunk (M4.8)	1	1
10			NLA - Locknut (M6 x 1)	1	1
11	B	8881169050	Hanger, Front Rigging, 80 Degree, Color	1	1
4			Support, Hanger, Footrest	1	1
5			Bracket, Front Rigging	1	1
8			Screw, Socket Head, Cap (M6 x 35mm)	1	1
9			Screw, Countersunk (M4.8)	1	1
10			NLA - Locknut (M6 x 1)	1	1
12			Tube, Upper, Front Rigging, Color	1	1
13			Housing, Legrest, 80 Degree	1	1
14			Screw, Stop (M10 x 1.5 x 23mm)	1	1
15			Screw, Stop (M10 x 1.5 x 12mm)	1	1
16	C	8881172207	Key, Cap	1	1
17			Hanger, Front Rigging, 90 Degree, Color	1	1
5			Bracket, Front Rigging	1	1
8			Screw, Socket Head, Cap (M6 x 35mm)	1	1
10			NLA - Locknut (M6 x 1)	1	1
18			Tube, 90 Degree Footrest, Color	1	1
NOTE: A - Includes items 2 - 10. B - Includes items 4, 5, 8, 9, 10, and 12 - 16. C - Includes items 5, 8, 10, and 18.					

Adult Elevating Legrest Hanger Assembly (ELRSS)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1508945	Hanger Assembly, Elevating Legrest - Left	1	1
1		1508948	Hanger Assembly, Elevating Legrest - Right	1	1




Adult Elevating Leg Rest Calf Pads (ELRSS)



Adult Elevating Leg Rest Calf Pads (ELRSS)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1169052	Calf Pad Assembly w/ Hardware (15 - 17" W)	1	1
2			Pad, Calf (15 - 17" W)	1	1
3			Plate, Attaching, Calf Pad (15 - 17" W)	1	1
4			Plate, Articulation, Female	1	1
5			Plate, Adjustment Angle	1	1
6			Plate, Articulation	1	1
7			Knob, Locking, Star	1	1
8			Pin, Spring (M8 x 55mm)	1	1
9			Pin, Spring (M8 x 50mm)	1	1
10			Clamp, Calf Pad Plate	1	4
11			Screw, Socket, Button Head (M6 x 16mm)	1	1
12			Spacer, Sleeve	1	1
13			Washer (2.4mm x 15.6mm)	1	1
14	A	1169053	Ring, Locking	1	1
15			Screw, VBA (M3.5 x 15)	1	4
1			Calf Pad Assembly w/ Hardware (18 - 20" W)	1	1
2			Pad, Calf (18 - 20" W)	1	1
3			Plate, Attaching, Calf Pad (18 - 20" W)	1	1
4			Plate, Articulation, Female	1	1
5			Plate, Adjustment Angle	1	1
6			Plate, Articulation	1	1
7			Knob, Locking, Star	1	1
8			Pin, Spring (M8 x 50mm)	1	1
9			Pin, Spring (M8 x 55mm)	1	1
10			Clamp, Calf Pad Plate	1	4
11	Screw, Socket, Button Head (M6 x 16mm)	1	1		
12	Spacer, Sleeve	1	1		
13	Washer (2.4mm x 15.6mm)	1	1		
14	Ring, Locking	1	1		
15	Screw, VBA (M3.5 x 15)	1	4		
NOTE: A - Includes items 2 - 15.					

Front Riggings - Junior 80 Degree Swing Hangers




Front Riggings - Junior 80 Degree Swing Hangers


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881514392	Hanger Assembly w/ Hardware, 80 Degree, Long - Right	1	1
3			Hanger Bar, 80 Degree, Long, Color - Right	1	1
4			Screw (M6 x 45mm)	1	1
5			Screw (M5 x .8 x 25mm)	1	1
6			Pivot, Hanger, 80 Degree	1	1
7			Cap, Cylinder, Hanger	1	1
8			Screw (M5 x .8 x 23mm)	1	3
9			Screw (M5 x .8 x 40mm)	1	1
10			Collar, Anti-Rattle	1	1
11			Washer (M6 x 12mm x 0.8mm)	1	1
12			Spring, Double	1	1
13			Lock, Latch Release, Female	1	1
14			Lock, Latch Release, Male	1	1
15			Housing, Latch	1	1
16			Bushing, Latch Housing	1	1
17			Nut, Lock (M5 x .8)	1	5
18			Bushing, Latch	1	1
1	A	8881514393	Hanger Assembly w/ Hardware, 80 Degree, Long - Left	1	1
3			Hanger Bar, 80 Degree, Long, Color - Left	1	1
4			Screw (M6 x 45mm)	1	1
5			Screw (M5 x .8 x 25mm)	1	1
6			Pivot, Hanger, 80 Degree	1	1
7			Cap, Cylinder, Hanger	1	1
8			Screw (M5 x .8 x 23mm)	1	3
9			Screw (M5 x .8 x 40mm)	1	1
10			Collar, Anti-Rattle	1	1
11			Washer (M6 x 12mm x 0.8mm)	1	1
12			Spring, Double	1	1
13			Lock, Latch Release, Female	1	1
14			Lock, Latch Release, Male	1	1
15			Housing, Latch	1	1
16			Bushing, Latch Housing	1	1
17			Nut, Lock (M5 x .8)	1	5
18			Bushing, Latch	1	1
20	B	1169143	Kit, Hardware, Block Sector	1	1
21			Screw (M6 x 1 x 60mm)	1	1
22			Position Latch	1	1

NOTE: A - Includes items 3, and 4-19.
B - Includes items 21 and 22.

Junior Elevating Legrest Hangers and Calf Pad (ELVOFF, ELVOFFWO)



Elevating Legrest Hanger



Elevating Legrest Calfpad


Junior Elevating Legrest Hangers and Calf Pad (ELVOFF, ELVOFFWO)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1450950	Hanger Assembly, Elevating Legrest (12" W)	1	2
2			Screw (M8 x 35mm)	1	1
3			Support, Hanger, Elevated Legrest	1	1
4			Indexer, Hanger, 90 Degree - Myon Junior	1	1
5			Tube, Upper, Elevating Legrest (12" W)	1	1
6			Tube, Lower, Elevating Legrest (12" W)	1	1
7			Screw, CBHC (M8 x 20mm)	1	1
8			Screw, Ball (7.5 mm)	1	1
9			Bushing, Armrest	1	1
10			Nut, Square, Armrest	1	1
1		1428808	Hanger Assembly, Elevating Leg Rest (13 - 15" W)	1	2
2			Screw (M8 x 35mm)	1	1
3			Support, Hanger, Elevated Legrest	1	1
4			Indexer, Hanger, 90 Degree - Myon Junior	1	1
5			Tube, Upper, Elevating Legrest (13 - 14" W)	1	1
6			Tube, Lower, Elevating Legrest (13 - 14" W)	1	1
7			Screw, CBHC (M8 x 20mm)	1	1
8			Screw, Ball (7.5 mm)	1	1
9			Bushing, Armrest	1	1
10			Nut, Square, Armrest	1	1
12		1403981	Calf Pad Assembly w/ Hardware (12 - 15" W)	1	1
13			Calf Pad (15 - 17" W)	1	1
14			Plate, Attaching, Calf Pad (15 - 17" W)	1	1
15			Mounting, Calf Pad Back	1	1
16			Button, Tightening	1	1
17			Screw, VBA (M3.5 x 15)	1	4
18			Clamp, Calf Pad Plate	1	2
19			Pin, Spring (M8 x 50mm)	1	1

Front Riggings - Accessories

Front Rigging	Description	Part Number
A1337	Calf Strap	1032243
		See Note(s) A
HLOOP FLDPADJ	Heel Loop	1492315
		See Note(s) B
HLOOP FLDPADJ	Heel Loop	1478467
		See Note(s) C
A9342	"H" Strap, 12-13" W	1032214
A9342	"H" Strap, 14-18" W	1032215
FRMPRBK	Neoprene Frame Protector Front - Black	1163422
FRMPRBK	Neoprene Frame Protector Back - Black	1163421
FRMPRBK	Neoprene Frame Protector Front - Blue	1159758
FRMPRBK	Neoprene Frame Protector Back - Blue	1159756
FRMPRBK	Neoprene Frame Protector Front - Red	1159759
FRMPRBK	Neoprene Frame Protector Back - Red	1159757
MYCVR	Monkey Cover	1511844
BRCVR	Bear Cover	1511843
NOTE: A - Used on all wheelchair models. B - Used on wheelchair models MYONADLT and MYONTS. C - Used on wheelchair models FXMYONJR and GRMYONJR.		

Footplate Assembly - Adult Standard (FLIPFTS, FLDPADJ, and 1PCFTPT)



SCM_028A

Footplate Assembly - Adult Standard (FLIPFTS, FLDPADJ, and 1PCFTPT)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Footplate Assembly, Flip Up (FLIPFTS)	1	2
2		See Chart	#2 Footplate Assembly, Flip Up Depth and Adjustable Angle (FLDPADJ)	1	2
3		See Chart	#3 One Piece Footplate Assembly (1PCFTPT)	1	1

#1 Footplate Assembly, Flip Up (FLIPFTS)

Footplate Type	Seat to Floor Height	15 - 16" W Left	15 - 16" W Right	17 - 18" W Left	17 - 18" W Right	19 - 20" W Left	19 - 20" W Right
FLIPFTS	14"	1169119	1169120	1169121	1169122	1169123	1169124
FLIPFTS	15 - 16"	1530620	1530614	1530622	1530616	1530623	1530617
FLIPFTS	17 - 20"	1530633	1530630	1530634	1530631	1530635	1530632

#2 Footplate Assembly, Flip Up Depth and Adjustable Angle (FLDPADJ)

Footplate Type	Chair Width	Seat to Floor Height FSTF14"	Seat to Floor Height FSTF15-16"	Seat to Floor Height FSTF17-20"
FLDPADJ	15" W - Left	1169131	1428657	1440816
FLDPADJ	15" W - Right	1169125	1428651	1440810
FLDPADJ	16" W - Left	1169132	1428658	1440817
FLDPADJ	16" W - Right	1169126	1428652	1440811
FLDPADJ	17" W - Left	1169133	1428659	1440818
FLDPADJ	17" W - Right	1169127	1428653	1440812
FLDPADJ	18" W - Left	1169134	1428660	1440819
FLDPADJ	18" W - Right	1169128	1428654	1440813
FLDPADJ	19" W - Left	1169135	1428661	1440820
FLDPADJ	19" W - Right	1169129	1428655	1440814
FLDPADJ	20" W - Left	1169136	1428662	1440821
FLDPADJ	20" W - Right	1169130	1428656	1440815




Footplate Assembly - Adult Standard

(FLIPFTS, FLDPADJ, and 1PCFTPT)

#3 One Piece Footplate Assembly (1PCFTPT)

Footplate Type	Seat to Floor Height	15" W	16" W	17" W	18" W	19" W	20" W
1PCFTPT	15 - 20"	1496945	1496946	1496947	1496948	1496949	1496953

Footplate Assembly - Adult Heavy Duty (FLIPFTS, FLDPADJ, and 1PCFTPT)



SCM_027A

Footplate Assembly - Adult Heavy Duty (FLIPFTS, FLDPADJ, and 1PCFTPT)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Footplate Assembly, Flip Up (FLIPFTS, FLDPADJ)	1	2
2		1532963	#2 Footplate Assembly, Flip Up Depth and Adjustable Angle	1	1
2		1532964	Footplate Assembly w/ Hardware, Adjustable Depth, Short - Left (20 - 24" W)	1	1
3		See Chart	#3 One Piece Footplate Assembly (1PCFTPT)	1	1




#1 Footplate Assembly, Flip Up (FLIPFTS, FLDPADJ)

Footplate Type	Seat to Floor Height	20" W Left	20" W Right	22 - 24" W Left	22 - 24" W Right
FLIPFTS	14"	1169123	1169124	1169137	1169138
FLIPFTS	15 - 16"	1530623	1530617	1488624	1488623
FLIPFTS	17 - 20"	1530635	1530632	1488626	1488625
FLDPADJ	15 - 20"	1532963	1532964	1532963	1532964

#3 One Piece Footplate Assembly (1PCFTPT)

Footplate Type	Seat to Floor Height	20" W (HD)	22" W (HD)	24" W (HD)
1PCFTPT	15 - 20"	1496953	1510382	1510383

Footplate Assembly - 80 Degree Supports



Footplate Assembly - 80 Degree Supports

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Footplates - 80 Degree Supports (FLIPFTS)	1	1
2		See Chart	#2 Footplates - 80 Degree Supports (FLPDADJ)	1	1
3		1528511	Footboard Assembly w/ Hardware, 80 Degree (8 - 11" W)	1	1
3		1528515	Footboard Assembly w/ Hardware, 80 Degree (12 - 15" W)	1	1






#1 Footplates - 80 Degree Supports (FLIPFTS)

Footrest Support Type	Footplate Type	Chair Width 12-13" W	Chair Width 14-15" W
80STRSHRT0 80STRSHRT30 80STRLONG0 80STRLONG30	FLIPFTS Right	1532606	1532608
80STRSHRT0 80STRSHRT30 80STRLONG0 80STRLONG30	FLIPFTS Left	1532607	1532609

#2 Footplates - 80 Degree Supports (FLPDADJ)

Footrest Support Type	Chair Width 10" W	Chair Width 11-13" W	Chair Width 14" W	Chair Width 15" W
80STRSHRT0 80STRSHRT30 80STRLONG0 80STRLONG30	1512257	1494144	1512255	1512259



Footplate Assembly and 90 Degree Lift Off Pivot Hanger (90STRLIFT w/ FLDADJ Footrests or 1PCFTPT)



Footplate Assembly and 90 Degree Lift Off Pivot Hanger (90STRLIFT w/ FLDADJ Footrests or 1PCFTPT)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1512289	Footplate Assembly w/ Hardware, Adjustable (10 - 11" W)	2	1
1	A	1512291	Footplate Assembly w/ Hardware, Adjustable (12 - 13" W)	2	1
1	A	1512287	Footplate Assembly w/ Hardware, Adjustable (14" W)	2	1
1	A	1512410	Footplate Assembly w/ Hardware, Adjustable (15" W)	2	1
2		1528773	Footboard Assembly w/ Hardware, 90 Degree (8 - 11" W)	1	1
2		1528818	Footboard Assembly w/ Hardware, 90 Degree (12 - 15" W)	1	1
3	A,B	1513492	Lift Off Pivot Hanger Assembly, 90 Degree	2	1
4	C	1497559	Kit, Insert, Legrest	1	1
5			Screw, Socket Button Head (M6 x 1 x 12mm)	1	1
6			Insert, Legrest Orienter	1	1
7	A	1511844	Cover, Hanger, 90 Degree, Monkey	2	1
8	A	1511843	Cover, Hanger, 90 Degree, Panda	2	1
NOTE: A - Accommodates both the right and left side of the wheelchair. B - Includes 2 ea. Screws, Supports and Indexers C - Includes items 5 and 6					


Footplates - Elevating Leg Rests (ELVOFF, ELVOFFWO)



Footplates - Elevating Leg Rests (ELVOFF, ELVOFFWO)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1450909	Footplate Assembly w/ Hardware, Composite, ELR - Left (12-13" W)	1	1
1	A	1450910	Footplate Assembly w/ Hardware, Composite, ELR - Right (12-13" W)	1	1
1	A	1449013	Footplate Assembly w/ Hardware, Composite, ELR - Left (14-15" W)	1	1
1	A	1449014	Footplate Assembly w/ Hardware, Composite, ELR - Right (14-15" W)	1	1
10	B	1507029	Footplate Assembly, Dual w/ Hardware, Adjustable, ELR (12 - 13" W)	1	1
10	B	1513240	Footplate Assembly, Dual w/ Hardware, Adjustable, ELR (15" W)	1	1
10	B	1513241	Footplate Assembly, Dual w/ Hardware, Adjustable, ELR (14" W)	1	1
NOTE: A - Includes items 2 - 7. B - Includes items 11 - 17.					

Casters - Caster, Headtube and Fork w/ Mtg. Hdwr.



Casters - Caster, Headtube and Fork w/ Mtg. Hdwr.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Casters	1	1
2		See Chart	#2 Headtube Forks w/ Hardware	1	1
3		1169075	Kit, Caster Attaching Hardware, Narrow Fork (Pair)	1	1
6			Bolt, Cap, Hex Head (M8 x 55mm)	1	2
7			NLA - Locknut (M8 x 1.25)	1	2
8	B		Spacer (5/16 x 15/32 x 1/4")	1	4
4		1173020	Kit, Caster Attaching Hardware, Wide (Pair)	1	1
5			NLA - Screw, Hex Head (M8 x 1.25 x 80mm)	1	2
7			NLA - Locknut (M8 x 1.25)	1	2
9			Spacer, Caster	1	4
10	C,E	1169073	Kit, Headtube Attaching Hardware (Pair)	1	1
11			Cap, Head Tube, Caster Housing	1	2
12			Plate, Mounting Tube Head	1	2
13			Plate, Mounting Caster, Inner	1	2
14			Insert, Support Fork	1	2
15	D		Screw, (M6 x 1 x 37mm)	1	4
10		1172240	Kit, Headtube Attachment Hardware (TRBKTS/TRRO)	1	1
11			Cap, Head Tube, Caster Housing	1	2
12			Plate, Mounting Tube Head	1	2
13			Plate, Mounting Caster, Inner	1	2
14			Insert, Support Fork	1	2
NOTE: A - Includes items 6, 7, and 8. B - Includes items 5, 7, and 9. C - Includes items 11 - 15. D - Includes items 11 - 14. E - Please review your Owners Manual for Torque Specification and Torque Label when replacing parts from this kit.					

#1 Casters

Caster Type	Wheel/Tire Type	Tire Size	Part No.	Mounting Hardware
5ALUM	Aluminum Urethane	5"	1026952	1169075
AT901	Composite Urethane	5"	1029496	1169075
AT5SOFT	Soft Roll	5 x 1-1/2"	1157706	1173020
AT902	Composite Urethane	6"	1029500	1169075
6PNEU	Composite Pneumatic	6 x 1-1/4"	1044040	1169075
AT903	Composite Urethane	8 x 1-1/4"	1166026	1169076
AT1292	Composite Pneumatic	8 x 1-1/4"	1166028	1169076

Casters - Caster, Headtube and Fork w/ Mtg. Hdwr.

#1 Casters

Caster Type	Wheel/Tire Type	Tire Size	Part No.	Mounting Hardware
1492	Pneumatic w/ FF Insert	8 x 1-1/4"	1166027	1169076

#2 Headtube Forks w/ Hardware

Caster Style	(FRKL1) Super Low Single Hole	(FRKS1) Normal Range Single Hole	(FRKS2) Normal Range Two Hole	(FRKS3) Normal Range Three Hole
5ALUM	1542007 See Note(s) A	1542007	1538091	1538090
AT901	1542007 See Note(s) A	1542007	1538091	1538090
AT5SOFT	1542007 See Note(s) A	1542007	1538092	1449113
AT902	1542007	1488399	1538092	1538090
6PNEU	1542007	1488399	1538092	1538090
AT903	1542007	N/A	N/A	1449113
AT1292	1542007	N/A	N/A	1449113
1492	1542007	N/A	N/A	1449113


NOTE: A - Also order Part Number 1449245 (2) when order replacement parts.

Axle Plates w/ Hardware (A1573 - All Wheels)


Standard and Heavy Duty

The camber is built into the axle plate. If the flat side is out it is 3 degree.
If the round surface is outside it is 0 degree.

A1573 - Standard Wheels
Quad Release Axle and Hardware



A1573 - Heavy Duty Wheels
Quad Release Axle and Hardware



Axle Plates w/ Hardware (A1573 - All Wheels)

Standard and Heavy Duty



The camber is built into the axle plate. If the flat side is out it is 3 degree.
If the round surface is outside it is 0 degree.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1173024	Kit, Axle Plate, Quad Release Hardware (Pair)	1	1
3			Axle, Quad Release	1	2
4			Bushing, Axle, Threaded	1	2
5			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	2
6			NLA - Locknut (M6 x 1)	1	2
7			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
8			Nut, Hex, Jam (3/4-16)	1	2
9			Plate, Quick Release Axle	1	2
2			B	1173026	Kit, Axle Plate, Quad Release Hardware, Heavy Duty (Pair)
3	Axle, Quad Release	1			2
4	Bushing, Axle, Threaded	1			2
7	Washer, Serrated (M20.3 x 30mm x 1.5mm)	1			2
8	Nut, Hex, Jam (3/4-16)	1			2
10	Bracket, Axle, Quick Release - Right	1			1
10	Bracket, Axle, Quick Release - Left	1			1
11	Screw, Scoket Head Cap (M6 x 50)	1			4
12	Nut, w/ Insert, HFR (M6)	1			4
NOTE: A - Kit includes items 3 - 9. B - Kit includes items 3, 4, 7, 8, 10, 11 and 12.					

Axle Plates w/ Hardware (A1577 - Composite Wheels)

Standard and Heavy Duty

The camber is built into the axle plate. If the flat side is out it is 3 degree.
If the round surface is outside it is 0 degree.



SCM_086A

Axle Plates w/ Hardware (A1577 - Composite Wheels)

Standard and Heavy Duty



The camber is built into the axle plate. If the flat side is out it is 3 degree.
If the round surface is outside it is 0 degree.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1169077	Kit, Hardware, Quick Release Axle Plate (Pair)	1	1
2			Axle, Quick Release Suspension	1	2
3			Bushing, Axle, Threaded	1	2
4			Nut, Hex, Jam (3/4-16)	1	4
5			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	4
6			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	2
7			NLA - Locknut (M6 x 1)	1	2
8			Plate, Quick Release Axle	1	2
9	B	1523612	Axle Release Plate Assembly w/ Hardware, Quick Extended (Pair)	1	1
2			Quick Release Suspension Axle	1	2
3			Bushing, Axle, Threaded	1	2
4			Nut, Hex, Jam (3/4-16)	1	4
5			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	4
10			Screw, Socket Head (M6 x 50)	1	4
11			Nut, w/ Insert, HFR (M6)	1	4
12			Bracket, Axle, Quick Release - Right	1	1
12			Bracket, Axle, Quick Release - Left	1	1
NOTE: A - Kit includes items 2 - 8. B - Kit includes items 2, 3, 4, 5, and items 10, 11, and 12					

Axle Plates w/ Hardware (A1577 - Spoked Wheels)

Standard and Heavy Duty

The camber is built into the axle plate. If the flat side is out it is 3 degree.
If the round surface is outside it is 0 degree.



SCM_087A

Axle Plates w/ Hardware (A1577 - Spoked Wheels)

Standard and Heavy Duty


The camber is built into the axle plate. If the flat side is out it is 3 degree.

If the round surface is outside it is 0 degree.

Item Number	Note	Part Number	Description	Quantity Per			
				Package	Assembly		
1	A	1173080	Kit, Axle Plate, Quick Release Hardware, Standard - Spoked Wheels (Pair)	1	1		
2			Axle, Quick Release Suspension	1	2		
3			Bushing, Axle Thread	1	2		
4			Nut, Hex, Jam (3/4-16)	1	4		
5			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	4		
10			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	2		
11			NLA - Locknut (M6 x 1)	1	2		
12			Plate, Quick Release Axle	1	2		
9			B	1173079	Kit, Axle Plate, Quick Release Hardware - Heavy Duty - Spoked Wheels (Pair)	1	1
2					Axle, Quick Release Suspension	1	2
3					Bushing, Axle Thread	1	2
4					Nut, Hex, Jam (3/4-16)	1	4
5	Washer, Serrated (M20.3 x 30mm x 1.5mm)	1			4		
6	Bracket, Axle, Quick Release - Left	1			1		
6	Bracket, Axle, Quick Release - Right	1			1		
7	Screw, Scocket Head Cap (M6 x 50)	1			4		
8		Nut, w/ Insert, HFR (M6)	1	4			
NOTE: A - Kit includes items 2 - 5, and items 10, 11, and 12. B - Kit includes items 2 - 8.							






Axle Plates w/ Hardware (Permanent Axle - All Wheels)

The camber is built into the axle plate. If the flat side is out it is 3 degree. If the round surface is outside it is 0 degree.



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1169078	Kit, Hardware, Fixed Axle Plate (Pair)	1	1
2			Screw, Cap, Hex Head (M12 x 1.75 x 115mm)	1	2
3			Spacer, Axle	1	2
4			Washer (M22 x 12mm x 2.5mm)	1	2
5			Nut, Hex (M12 X 1.75)	1	2
6			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	2
7			NLA - Locknut (M6 x 1)	1	2
8			Plate, Axle, Fixed	1	2
NOTE: A - Kit includes items 2 - 8.					

Axle Plates - 20 Inch Frame



SCM_017B

Axle Plates - 20 Inch Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1188533	Kit, Axle Plate Assembly, Dual Amputee, Quick Release	1	1
2			Support, Rear Wheel	1	4
3			Plate, Support, Axle, Dual Amputee	1	2
4			Axle, Quick Release Suspension	1	2
5			Bushing, Axle, Threaded	1	2
7			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	4
8			Screw, Socket Button Head (M6 x 53mm)	1	4
9			Spacer (13 x 16 x 5)	1	6
10			Bolt, Hex Head Cap (M12 x 1.75 x 40mm)	1	6
11			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	4
12			Nut, Hex, Jam (3/4-16)	1	4
13			NLA - Locknut (M6 x 1)	1	8
14	B	1188534	Kit, Axle Plate Assembly, Dual Amputee, Quick Release, Spoked Wheel	1	1
2			Support, Rear Wheel	1	4
3			Plate, Support, Axle, Dual Amputee	1	2
7			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	4
8			Screw, Socket Button Head (M6 x 53mm)	1	4
9			Spacer (13 x 16 x 5)	1	6
10			Bolt, Hex Head Cap (M12 x 1.75 x 40mm)	1	6
11			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	4
12			Nut, Hex, Jam (3/4-16)	1	4
13			NLA - Locknut (M6 x 1)	1	8
15			Axle, Quick Release Suspension	1	2
22			Bushing, Axle Thread	1	2
16	C	1188535	Kit, Axle Plate Assembly, Quick Release	1	1
2			Support, Rear Wheel	1	2
4			Axle, Quick Release Suspension	1	2
7			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	2
11			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	4
12			Nut, Hex, Jam (3/4-16)	1	4
13			NLA - Locknut (M6 x 1)	1	2
23	D	1188536	Bushing, Axle, Threaded	1	2
18			Kit, Axle Plate Assembly, Quick Release, Spoked Wheels	1	1
2			Support, Rear Wheel	1	2
7			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	2
11			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	4
12			Nut, Hex, Jam (3/4-16)	1	4
13			NLA - Locknut (M6 x 1)	1	2
15	E	1188538	Axle, Quick Release Suspension	1	2
25			Bushing, Axle Thread	1	2
20			Kit, Axle Plate Assembly, Quad Release	1	1
2			Support, Rear Wheel	1	2

Axle Plates - 20 Inch Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
5			Bushing, Axle, Threaded	1	2
7			Screw, C-Sunk Head (M6 x 1 x 40mm)	1	2
11			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
12			Nut, Hex, Jam (3/4-16)	1	2
13			NLA - Locknut (M6 x 1)	1	2
21			Axle, Quad Release	1	2
NOTE: A - Kit includes items 2 - 5, and 7 - 13. B - Kit includes items 2, 3, and items 7 - 13, 15, and item 22. C - Kit includes items 2, 4, 7, 11 - 13, and item 24. D - Kit includes items 2, 7, items 11 - 13, item 15, and item 25. E - Kit includes items 2, 4, 6, 7, 11, 12, 14, and item 21.					

Rear Wheels

Composite, Spoked, Sun Fusion, and Spinergy

Rear Wheels - Composite (COM)

Wheel Style	Tire Type	Complete Wheel Assembly
18COM (18")	Pneumatic w/ Flat Free	1137256
20COM (20")	Pneumatic	1026422
20COM (20")	Pneumatic w/ Flat Free	1026423
20COM (20")	Urethane	1038364
22COM (22")	Pneumatic	1026420
22COM (22")	Pneumatic w/ Flat Free	1026421
22COM (22")	Urethane	1032116
24COM7 (24")	Pneumatic	1026362-NLA
24COM7 (24")	Pneumatic w/ Flat Free	1026363
24COM7 (24")	Urethane	1032118

Rear Wheels

Composite, Spoked, Sun Fusion, and Spinergy

Rear Wheels - Spoked (ITESPK)

Wheel Style	Tire Type	Complete Wheel Assembly (w/ Tire, Tube, and/or Rim Strip)	Tire	Tube/ Insert	Rim Strip
22ITESPK (22")	Pneumatic	8881159332 See Note(s) A	1018544	1018545	1018546
22ITESPK (22")	Pneumatic w/ Flat Free	8881159332 See Note(s) A	1018544	1023323	N/A
24ITESPK (24")	Pneumatic	8881159332	40176X060	40176X061	40176X008
24ITESPK (24")	Pneumatic w/ Flat Free	8881159332 See Note(s) A	40176X060	1010906	N/A
24ITESPK (24")	Primo	8881159332 See Note(s) A	1026138	1032240	40176X008
24ITESPK (24")	Schwalbe	8881159332 See Note(s) A	1139813	1032240	40176X008
25ITESPK (25")	Primo	8881159332 See Note(s) A	1044947	1029471	40176X008
25ITESPK (25")	Schwalbe	8881159332 See Note(s) A	1139814	1029471	40176X008
26ITESPK (26")	Pneumatic	8881159332 See Note(s) A	1026142	1026140	1026141-NLA
26ITESPK (26")	Primo	8881159332 See Note(s) A	1029470	1029471	1026141-NLA
26ITESPK (26")	Schwalbe	8881159332 See Note(s) A	1139815	1029471	1026141-NLA

NOTE: A - If you need to order these Wheels without a Handrim, it must be ordered as a special thru the configurator.

Rear Wheels

Composite, Spoked, Sun Fusion, and Spinergy

Rear Wheels - KIK Tire Spoked (ITESPK)

Wheel Style	Tire Type	Wheel Assembly	Tire	Tube/Insert
24ITESPK (24")	KIK	1172782	1062852	N/A

Rear Wheels - Sun Fusion (FN)

Wheel Style	Tire Type	Wheel Assy. (w/ Brgs.)	Tire Number	Tube Number	Rim Strip
FN1624 (24")	Pneumatic	1163457	40176X060	40176X061	40176X008
FN1624 (24")	Primo	1163457	1026138	1105941	40176X008
FN1624 (24")	KIK Black	1163457	1062852	N/A	N/A
FN1624 (24")	Schwalbe	1163457	1139813	1105941	40176X008
FN1625 (25")	Primo	1163458-NLA	1044947	1104326	40176X008
FN1625 (25")	Primo Treaded	1163458-NLA	1130527	1104326	40176X008
FN1625 (25")	KIK Black	1163458-NLA	1144435	N/A	N/A
FN1625 (25")	Schwalbe	1163458-NLA	1139814	1104326	40176X008
FN1626 (26")	Pneumatic	1163459	1026142	1026140	1026141-NLA
FN1626 (26")	Primo	1163459	1029470	1104326	1026141-NLA
FN1626 (26")	Schwalbe	1163459	1139815	1104326	1026141-NLA

Rear Wheels

Composite, Spoked, Sun Fusion, and Spinergy

Rear Wheel w/o Tire - Spinergy Spox Wheels (SPIN)

Wheel w/ Bearings

Wheel Size	Black (BLK)	Yellow (YEL)	White (WHT)	Blue (BLU)	Red (RED)
SPIN24 (24")	1144427	1144428	1159103-NLA	1159102-NLA	1159101-NLA
SPIN25 (25")	1144429	1144430-NLA	1159100-NLA	1159099-NLA	1159098-NLA
SPIN26 (26")	1144431-NLA	1144432	1159097	1159096	1159095

Rear Wheel Tires, Tubes, and Rim Strip- Spinergy Spox Wheels (SPIN)

Tire, Tube, and Rim Strip

Wheel Style	Tire Type	Tire Part No.	Tube Part No.	Rim Strip Part No.
SPIN24 (24")	Primo	1026138	1104326	40176X008
SPIN24 (24")	Schwalbe	1139813	1105941	40176X008
SPIN25 (25")	Primo	1044947	1104326	40176X008
SPIN25 (25")	Primo Treaded	1130527	1104326	40176X008
SPIN25 (25")	Schwalbe	1139814	1104326	40176X008
SPIN26 (26")	Primo	1029470	1104326	1026141-NLA
SPIN26 (26")	Schwalbe	1139815	1104326	1026141-NLA

Handrims and Mounting Hardware

Rear Wheels - Handrims

Composite Wheels (COM)

Handrim	18COM	20COM	22COM	24COM7
Aluminum Anodized (ALUMPC)	N/A	1033374	1033373	1033372
Black Anodized (BLACKANO)	N/A	N/A	N/A	N/A
Plastic Coated (A1512)	N/A	1115400	1115401	1115397
Plastic - 8" Oblique (AT1224-8)	N/A	N/A	N/A	1115396
Plastic - 8" Vertical (AT1225-8)	N/A	N/A	N/A	1115394
Plastic - 12" Vertical (AT1225-12)	N/A	N/A	N/A	1115395
Aluminum Handrims (AHDRH)	1059369	N/A	N/A	N/A

Rear Wheels - Handrims

Spoked Wheels (ITESPK)

Handrim	22ITESPK	24ITESPK	25ITESPK	26ITESPK
Complete Wheel Assembly	8881159332	8881159332	8881159332	8881159332
Aluminum Anodized (ALUMPC)	1043443	1026136	1029493	1026137
Black Anodized (BLACKANO)	N/A	1159659	N/A	1159661
Plastic Coated (A1512)	1043689-NLA	1035612	1043249	1035613

Handrims and Mounting Hardware

Rear Wheels - Handrims

Sun Fusion (FN)

Handrim Style	FN1624	FN1625	FN1626
Aluminum Anodized (ALUMPC)	1114300	1114301-NLA	1114302
Black Anodized (BLACKANO)	1160753-NLA	1160754-NLA	1160755
Plastic Coated (A1512)	1168964	1043249	1168965-NLA

Rear Wheels - Handrims

Spinergy Spox (SPIN)

Handrim Material	SPIN24 (24")	SPIN25 (25")	SPIN26 (26")
Aluminum Anodized (ALUMPC)	1114300	1114301-NLA	1114302
Black Anodized (BLACKANO)	1160753-NLA	1160754-NLA	1160755
Plastic Coated (A1512)	1168964	1043249	1168965-NLA

Rear Wheel - Handrim Hardware

Handrim Style	Kit Part Number	Quantity per Chair
AT1225-8	1026380	2
AT1224-8	1026381	2
AT1225-12	1026382	2
A1512	1026383	2
ALUMPC and CHDRM	1026384	2
P1512	1026388	2
COM	1065790	2

Handrims - Natural Fit Handrims

Wheel Style	Natural Fit LT Standard Grip (NFSDLT)	Natural Fit LT Super Grip (NFSPLT)	Natural Fit Standard Grip (NFITSD)	Natural Fit Super Grip (NFITSP)	Mounting Hardware
24COM7 (24")	1163239	1163239	1144488	1144488	1145680-NLA
24ITESPK (24")	1163241	1163241	1146264	1146264	1160547
25ITESPK (25")	1163242	1163242	1146265	1146265	1160547
26ITESPK (26")	1163243	1163243	1146266-NLA	1146266-NLA	1160547
FN1624 (24")	1163241	1163241	1146264	1146264	1160547
FN1625 (25")	1163242	1163242	1146265	1146265	1160547
FN1626 (26")	1163243	1163243	1146266-NLA	1146266-NLA	1160547
SPIN24 (24")	1163241	1163241	1146264	1146264	1160547
SPIN25 (25")	1163242	1163242	1146265	1146265	1160547
SPIN26 (26")	1163243	1163243	1146266-NLA	1146266-NLA	1160547

Diagram illustrating the installation of the front and rear seat frames onto the side frame.

Top Diagram (Front Seat Frame):

- 1: Front seat frame.
- 3: Side frame (shown for reference only).
- 4: Bracket used to attach the front seat frame to the side frame.
- 5: Screws used to secure the bracket.

Bottom Diagram (Rear Seat Frame):

- 2: Rear seat frame.
- 3: Side frame (shown for reference only).
- 4: Bracket used to attach the rear seat frame to the side frame.
- 5: Screws used to secure the bracket.


Wheel Locks - Pull to Lock

Pull to Lock (WLPULL), Pull to Lock Extended (WLPULLX)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1538079	Wheel Lock Assembly w/ Hardware, Pull to Lock - Right	1	1
1		1538078	Wheel Lock Assembly w/ Hardware, Pull to Lock - Left	1	1
2		1541957	Wheel Lock Assembly w/ Hardware, Pull to Lock w/ Extensions - Right	1	1
2		1541956	Wheel Lock Assembly w/ Hardware, Pull to Lock w/ Extensions - Left	1	1
3		1423107	Kit, Hardware, Clamp Brake (Pair)	1	1
4			Screw (M6 x 1 x 35mm)	1	4
5			Mounting Clamp, Wheel Lock	1	2
NOTE: A - Includes items 4 and 5.					

Wheel Locks - Push to Lock

Push to Lock (WLPUSH), Push to Lock Extended (WLPUSHX)




Wheel Locks - Push to Lock

Push to Lock (WLPUSH), Push to Lock Extended (WLPUSHX)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1425634	Wheel Lock Assembly, Push to Lock - Left	1	1
1		1425635	Wheel Lock Assembly, Push to Lock - Right	1	1
2		1425602	Wheel Lock Assembly, Push to Lock, Extended - Left	1	1
2		1425603	Wheel Lock Assembly, Push to Lock, Extended - Right	1	1
	A	1425642	Wheel Lock Assembly w/ Hardware, Dual Amputee - Left	1	1
	A	1425643	Wheel Lock Assembly w/ Hardware, Dual Amputee - Right	1	1
3	B	1423107	Kit, Hardware, Clamp Brake (Pair)	1	1
4			Screw (M6 x 1 x 35mm)	1	4
5			Mounting Clamp, Wheel Lock	1	2
NOTE: A - Not shown in illustration. B - Includes items 4 and 5.					




Wheel Locks - Hill Holder Assembly (A1598)



Wheel Locks - Hill Holder Assembly (A1598)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1039323-NLA	NLA - Hill Holder Assembly w/ Hardware - Right	1	1
3			NLA - Strap, Wheel Lock	1	1
4			NLA - Screw, Button Head Cap (1/4-28 x 7/8")	1	1
5			NLA - Screw, Socket Button Head (1/4-28 x 1-1/4")	1	1
6			NLA - Spacer (1/4 x 1/2 x 1/2")	1	1
7			NLA - Spacer (1/4 x 1/2 x 5/32")	1	1
8			Locknut (1/4-28)	1	5
9			NLA - Bracket Assembly, Black - Right	1	1
10			NLA - Link, Stop, Black	1	1
11			NLA - Handle, Wheel Lock, Black - Right	1	1
12			Tip, Plastic, Black	1	1
13			Screw, Socket Button Head (1/4-28 x 5/8")	1	2
14			Shoe, Wheel Lock	1	1
15			Washer, Nylon (1/4 x 7/16 x 5/64")	1	1
16			Washer, Nylon (1/4 x 1/2 x 1/32")	1	7
17			NLA - Spring, Torsion - Right	1	1
1			B	1039324-NLA	NLA - Hill Holder Assembly w/ Hardware - Left
3	NLA - Strap, Wheel Lock	1			1
4	NLA - Screw, Button Head Cap (1/4-28 x 7/8")	1			1
5	NLA - Screw, Socket Button Head (1/4-28 x 1-1/4")	1			1
6	NLA - Spacer (1/4 x 1/2 x 1/2")	1			1
7	NLA - Spacer (1/4 x 1/2 x 5/32")	1			1
8	Locknut (1/4-28)	1			5
9	NLA - Bracket Assembly, Black - Left	1			1
10	NLA - Link, Stop, Black	1			1
11	NLA - Handle, Wheel Lock, Black - Left	1			1
12	Tip, Plastic, Black	1			1
13	Screw, Socket Button Head (1/4-28 x 5/8")	1			2
14	Shoe, Wheel Lock	1			1
15	Washer, Nylon (1/4 x 7/16 x 5/64")	1			1
16	Washer, Nylon (1/4 x 1/2 x 1/32")	1			7
17	NLA - Spring, Torsion - Left	1			1
2		1039343			Kit, Clamp Wheel Lock Hardware (Pair)
18			Clamp, Wheel Lock, Threaded (7/8")	1	2
19			NLA - Clamp, Wheel Lock, Non-Threaded (7/8")	1	2
20			Screw, Allen Head (1/4-20 x 1")	1	4
NOTE: A - Includes items 3 - 17. B - Includes items 18 - 20					

Acrylic Trays



Acrylic Trays

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Acrylic Tray - Standard Trays - Before 10/17/18	1	1
1		See Chart	#1 Acrylic Tray - Standard Trays - After 10/16/18	1	1
2		See Chart	#2 Acrylic Tray - Swing Out	1	1
3		L10024	Knob, Tray	1	1

#1 Acrylic Tray - Standard Trays - Before 10/17/18

Chair Width	Tray (ACTRY) Part Number
12" Wide	1492757
13" Wide	1492758
14" Wide	1492759
15" Wide	1450921
16" Wide	1450922
17" Wide	1450923
18" Wide	1450924
19" Wide	1450925
20" Wide	1450926
Item includes Tray and Mat.	

#1 Acrylic Tray - Standard Trays - After 10/16/18





Chair Width	Tray (ACTRY) Part Number
15-16" Wide	1601454
17-18" Wide	1601458
19-20" Wide	1601459
Item includes Tray and Mat.	

Acrylic Trays

#2 Acrylic Tray - Swing Out

Models	8 - 9" Wide Chair	10 - 11" Wide Chair	12 - 13" Wide Chair	14 - 15" Wide Chair
ACTRSA ACTLSA	1513253	1513254	1513255	1513256
Item includes Tray, Knob, Support and Screws				

Anti-Tippers (A1580, A1580P)



SCM_033A

Anti-Tippers (A1580, A1580P)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1147653	Anti-Tip Assembly (Pair)	2	1
			Tube, Anti Tip Assembly	1	2
			NLA - Label, Warning, Anti-Tip	1	2
			NLA - Label, Warning, Anti-Tip	1	2
2	B	1173073	Anti-Tip Assembly, Low Seat to Floor (Pair)	2	1
			Tube, Anti-Tip Assembly, Low Seat to Floor	1	2
			NLA - Label, Warning, Anti-Tip	1	2
			NLA - Label, Warning, Anti-Tip	1	2
3	C	1173027	Anti-Tip Assembly, Super Low Seat to Floor (Pair)	2	1
			NLA - Tube, Anti-Tip Assembly, Super Low	1	2
			NLA - Label, Warning, Anti-Tip	1	2
			NLA - Label, Warning, Anti-Tip	1	2
4	D	1173028	Anti-Tip Assembly w/ Tipping Aid (Pair)	2	1
			Tube, Anti-Tip Assembly	1	2
			Tube, Anti-Tip Assembly w/ Tipping Aid	1	2
			Cap, Foot Support	1	2
			Screw, Hex Head Cap (M5 x .8 x 45mm)	1	2
			Wheel, 37.6	1	4
			Nut, Lock (M5 x .8)	1	2
			Spring, Torsion with Button Lock	1	2
			Button, Spring	1	2
			NLA - Anti-Rattle (1")	1	2
			NLA - Label, Warning, Anti-Tip	1	2
			NLA - Label, Warning, Anti-Tip	1	2
		1147676	Cap Kit, Step Tube	1	1
			Cap, Tube 25	1	2
		See Chart	Anti-Tipper Cross Reference	1	1



NOTE: A - Used with FBH, ABH, and AABH Series Back Canes, except Chairs with FRKL1 and FRKS1 Forks.
B - Used with FBH, ABH, and AABH Series Back Canes and FRKS1 Forks.
C - Used with FBH, ABH, and AABH Series Back Canes and FRKL1 Forks.
D - Not shown in illustration, includes Tube Cap (2).

Anti-Tipper Cross Reference

	FBH, ABH, and AABH Series Except for FRKL1 (RSTF14-15) Except for FRKS1 (RSTF15-16)	FBH, ABH, and AABH Series For FRKS1 (RSTF15-16)	FBH, ABH, and AABH Series For FRKL1 (RSTF14-15)	A1580P
Anti-Tipper Part No.	1147653	1173073	1173027	1173028

Anti-Tippers **(A1580, A1580P)**

Anti-Tippers, Long (For 19 and 20" FSTF) (A1580, A1580P)




Anti-Tippers, Long (For 19 and 20" FSTF) (A1580, A1580P)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	60108489	Anti-Tip Assembly, Long (Pair)	2	1
			Anti-Tip Assembly, Long	1	2
			NLA - Label, Warning, Anti-Tip	1	2
			NLA - Label, Warning, Anti-Tip	1	2
2	A	60108492	Anti-Tip Assembly w/ Tipping Aid, Long (Pair)	2	1
			Tube, Grey Anti Rattle	1	2
			Tube, Anti-Tip w/ Tipping Aid Long	1	2
			Cap, Foot Support	1	2
			Screw, Hex Head Cap (M5 x .8 x 45mm)	1	2
			Wheel, 37.6	1	4
			Nut, Lock (M5 x .8)	1	2
			Spring, Torsion with Button Lock	1	2
			Button, Spring	1	2
			NLA - Anti-Rattle (1")	1	2
			NLA - Label, Warning, Anti-Tip	1	2
			NLA - Label, Warning, Anti-Tip	1	2


NOTE: A - Used for 19" and 20" Front Seat To Floor (FSTF) heights only.

Crutch/Cane Carrier (A1569)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	0572-02056	Cane Holder Assembly w/ Hardware	1	1
2			Strap, Hook Loop, Cane Holder	1	1
3			Holder, Cane	1	1
4			Clamp, Cane Holder	1	2
5			Spacer, Cane Clamp Holder	1	2
6			Bolt, Hex Head (M6 x 45mm)	1	2
7			NLA - Locknut (M6 x 1)	1	2
NOTE: A - Includes items 2 - 7.					

Headrest and Backbrace



SCM_065A

Headrest and Backbrace

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1192337	Kit, Headrest and Backbrace, Seat (8 - 11", 1" Back Canes)	1	1
2			Headrest Assembly	1	1
3			Backbrace Assembly, TLSCP 205-280	1	1
			Instruction Sheet, Service, Myon	1	1
1		1192338	Kit, Headrest and Backbrace, Seat (8 - 11", 7/8" Back Canes)	1	1
2			Headrest Assembly	1	1
3			Backbrace Assembly, TLSCP 205-280	1	1
			Instruction Sheet, Service, Myon	1	1
1		1192339	Kit, Headrest and Backbrace, Seat (12 - 15", 1" Back Canes)	1	1
2			Headrest Assembly	1	1
3			Backbrace Assembly, TLSCP 205-280	1	1
			Instruction Sheet, Service, Myon	1	1
1		1192340	Kit, Headrest and Backbrace, Seat (12 - 15", 7/8" Back Canes)	1	1
2			Headrest Assembly	1	1
3			Backbrace Assembly, TLSCP 205-280	1	1
			Instruction Sheet, Service, Myon	1	1

Headrest and Backbrace - Where Used


Back Type	Part Number
ABH1216PH, ABH1519PH ABH1620PH, AABH1216 AABH1519, AABH1620 ABH1216, ABH1519, ABH1620	1192337
FBH16, FBH20, FBH16F, FBH20F	1192338
ABH1216PH, ABH1519PH ABH1620PH, AABH1216 AABH1519, AABH1620 ABH1216, ABH1519, ABH1620	1192339
FBH16, FBH20, FBH16F, FBH20F	1192340

Junior Spoke Protectors


Spoke Protector Style	20" Wheels	22" Wheels
A1307M	1492138	1522862
A1307SP	1513257	1522863
A1307FL	1513258	1522864

One Arm Drive

Urethane (DUHALU, DUHARU), Pneumatic (DUHARP, DUHARL)



Solid Tire



Pneumatic Tire

One Arm Drive

Urethane (DUHALU, DUHARU), Pneumatic (DUHARP, DUHARL)


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1513269	One Arm Drive Assembly, Pneumatic w/ Hardware (15 - 17" W)	1	1
2			Shaft, Telescopic, One Arm Drive (15 - 17" W)	1	1
3			Wheel, Active, Pneumatic, DHR (24")	1	1
4			Wheel, Passive, Pneumatic, DHR (24")	1	1
5			Bushing, Quick Release	1	2
6			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
7			Nut, One Arm Drive	1	2
1	A	1513270	One Arm Drive Assembly, Pneumatic w/ Hardware (18 - 20" W)	1	1
2			Shaft, Telescopic, One Arm Drive (18 - 20" W)	1	1
3			Wheel, Active, Pneumatic, DHR (24")	1	1
4			Wheel, Passive, Pneumatic, DHR (24")	1	1
5			Bushing, Quick Release	1	2
6			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
7			Nut, One Arm Drive	1	2
1	A	1513271	One Arm Drive Assembly, Solid w/ Hardware (15 - 17" W)	1	1
2			Shaft, Telescopic, One Arm Drive (15 - 17" W)	1	1
3			Wheel, Solid, Active, DHR (24")	1	1
4			Wheel, Solid, Passive (24")	1	1
5			Bushing, Quick Release	1	2
6			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
7			Nut, One Arm Drive	1	2
1	A	1513272	One Arm Drive Assembly, Solid w/ Hardware (18 - 20" W)	1	1
2			Shaft, Telescopic, One Arm Drive (18 - 20" W)	1	1
3			Wheel, Solid, Active, DHR (24")	1	1
4			Wheel, Solid, Passive (24")	1	1
5			Bushing, Quick Release	1	2
6			Washer, Serrated (M20.3 x 30mm x 1.5mm)	1	2
7			Nut, One Arm Drive	1	2

NOTE: A - Includes items 3 - 7.

Positioning Straps

Chair Widths	A1321 Hook and Loop Seat Positioning Straps	A1515 Auto Style Seat Positioning Straps	9515 Hook and Loop Chest Positioning Straps	7321 Auto Style Chest Positioning Straps	TRRO Transport Lap Belt
15 - 19" W	1147654	1147656	N/A	N/A	N/A
20" W 22 - 24" W (HD)	1147655	1147657	N/A	N/A	N/A
11 - 15" W	1147654	1147656	1004310	1016945-NLA	N/A
11 - 15" W	N/A	N/A	N/A	1016946-NLA	N/A
12 - 15" W	N/A	N/A	N/A	N/A	1162581
16 - 20" W	N/A	N/A	N/A	N/A	1162582
Includes Buckle, Strap and Fasteners.					

Push Handles



SCM_060A

Push Handles

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1192341	Kit, Handle, Stroller 7/8" Cane	1	1
			Instruction Sheet, Service, Myon	1	1
	B	1557041	Kit, Push Handles w/ Hardware (7/8" Cane)	1	1
			Push Bar Assembly	1	2
			Swivel Coupling	1	2
			Ring, Push, Reduction Bar (7/8")	1	4
			Lever, Indexable, Push Handle	1	2
			Screw (M8 x 35 A2K 8.8 DIN931)	1	2
			Washer (22 x 8.4 x 1.5 ZP NFE25513)	1	2
			Spring, Torsion with Button Lock	1	2
2			Button, Plug (22mm)	1	2
	C	1192342	Kit, Handle, Stroller 1" Cane	1	1
			Instruction Sheet, Service, Myon	1	1
	D	1192345	Kit, Push Handles w/ Hardware (1" Cane)	1	1
			Push Bar Assembly	1	2
			Swivel Coupling	1	2
			Ring, Push, Reduction Bar (7/8")	1	2
			Ring, Push, Reduction Bar (1")	1	2
			Lever, Indexable, Push Handle	1	2
			Screw (M8 x 35 A2K 8.8 DIN931)	1	2
3			Washer (22 x 8.4 x 1.5 ZP NFE25513)	1	2
			Spring, Torsion with Button Lock	1	2
			Button, Plug (22mm)	1	2
	E	1557061	Kit, Swivel Coupling w/ Lever	1	1
			Swivel Coupling	1	2
			Ring, Push, Reduction Bar (7/8")	1	4
			Lever, Indexable, Push Handle	1	2
			Screw (M8 x 35 A2K 8.8 DIN931)	1	2
			Washer (22 x 8.4 x 1.5 ZP NFE25513)	1	2
	F	1557062	Kit, Push Arm Assembly w/ Hardware	1	1
4			Push Bar Assembly	1	2
			Spring, Torsion with Button Lock	1	2
			Button, Plug (22mm)	1	2

NOTE: A - Kit includes 1557041 and Instruction Sheet.
 B - Kit includes items 2 - 4, and 6 - 10.
 C - Kit includes 1192345 and Instruction Sheet.
 D - Item 5 used with 1192345 only.
 E - Kit includes items 3, 4, 6, 7, 8.
 F - Kit includes items 2, 9, and 10.

Push Handles - Where Used


Back Type	Part Number
-----------	-------------

Push Handles

Push Handles - Where Used

Back Type	Part Number
ABH1216PH, ABH1519PH ABH1620PH, AABH1216 AABH1519, AABH1620 ABH1216, ABH1519, ABH1620	1192341
FBH16, FBH20, FBH16F, FBH20F	1192342

Accessories - Stabilizer Kit




Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	1169697	Kit, Stabilizer w/ Hardware (Pair)	1	1
2			Stabilizer, Fork, Myon Junior	1	2
3			Screw, (M6 x 1 x 37mm)	1	4
4			Screw, Socket Head, Cap (M6 x 1 x 37mm)	1	2
5			Nut (M6)	1	6
6			Washer (M6 x 12mm x 1.2mm)	1	6
NOTE: A - Includes items 2 - 6. B - Review Owner's Manual for specific installation instructions before replacing these parts on the chair.					

Standard Spoke Protectors

Wheel Style	20" Wheels	22" Wheels	24" Wheels	25" Wheels	26" Wheels
20COM	1032146-NLA	N/A	N/A	N/A	N/A
22COM 22ITESPK	N/A	1032147-NLA	N/A	N/A	N/A
24COM7	N/A	N/A	1028442-NLA	N/A	N/A
24ITESPK FN1624	N/A	N/A	1028441-NLA	N/A	N/A
25ITESPK FN1625 SPIN25	N/A	N/A	N/A	1028443-NLA	N/A
26ITESPK FN1626 SPIN26	N/A	N/A	N/A	N/A	1028443-NLA

Stroller Handle (STRHNDL)



Stroller Handle (STRHNDL)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1513242	Push Handle Assembly w/ Hardware (8 - 11" W)	1	1
2			Bar, Push, Attendant (8 - 11" W)	1	1
3			Bracket, Attaching, Pushbar	1	2
4			Lever, Indexable, Push Handle	1	2
5			Adapter, Brake Handle, Push Bar	1	2
6			Interface, Push Bar	1	2
7			Cap, Plastic (17 mm)	1	2
1	A	1513243	Push Handle Assembly w/ Hardware (12 - 15" W)	1	1
2			Bar, Push, Attendant (12 - 15" W)	1	1
3			Bracket, Attaching, Pushbar	1	2
4			Lever, Indexable, Push Handle	1	2
5			Adapter, Brake Handle, Push Bar	1	2
6			Interface, Push Bar	1	2
7			Cap, Plastic (17 mm)	1	2
NOTE: A - Includes items 2 - 7.					

Invacare Corporation

www.invacare.com



One Invacare Way
Elyria, Oh USA
44036-2125
800-333-6900

Technical Support
800-832-4707
Customer Service
800-333-6900

Literature Orders
800-333-6900

© 2022 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. (Revised 04/06/22)

Form No. 11-013