


Crossbrace Assembly 6626Y2

This model was discontinued 10/4/00

Model 6626Y2 (Discontinued)
(36" - 69" Wide)

Model 6626Y2 Crossbrace can be used to convert models 6625, 6628, or 6658 to a Home Bed Rail.



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		16356X000	Barrel Assembly, Right Head End / Left Foot End (Includes: Plunger Knob and Pin)	1	2
1		16326X000	Barrel Assembly, Left Head End / Right Foot End (Includes: Plunger Knob and Pin)	1	2
2		48806X011-NLA	NLA - Button, Snap	1	4
3		40188X005-NLA	NLA - Spring, Torsion	1	4
4	A	06626X003-NLA	NLA - Tube, Channel, Chrome - Outside	1	2
4	A	06626X004-NLA	NLA - Tube, Channel, Chrome - Inside	1	2
5		1011239-NLA	NLA - Label, Caution	1	2
6		60290X002	Locknut (5/16-18)	1	4
7		16308X000	Washer, Flat (5/16 x 11/16 x 1/16")	1	4
8	B	06626X002-NLA	NLA - Crossbrace Assembly	1	2

NOTE: A - Channel Tube is No Longer Available. Order Item number 8.
B - Includes Items 1-7.