


Seat Frames & Frame Hardware



Seat Frames & Frame Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		15680X013	Anti-Rattle (1")	1	4
2		1029501	Plug Button, Black (1")	1	2
3		1025704	Washer, Flat (5/16 x 5/8 x 1/16")	1	4
4		1028769	Locknut (5/16-18)	1	4
5		1046494	Screw, Socket Button Head (5/16-18 x 7/8")	1	4
6		1046497-NLA	NLA - Washer, Nylon (5/16 x 7/8 x 1/64")	1	4
7		1046496-NLA	NLA - Pad, Foam	1	2
8		40117X021	Bearing (1/2" I.D. x 1-1/8" O.D.)	1	4
9		1024749	Clamp	1	2
10		1026911-NLA	NLA - Screw, Socket Set Cone Point (1/4-20 x 1/2")	1	2
11		See Chart	#11 Seat Frames and Frame Hardware	1	1

#11 Seat Frames

			Seat Depth				
Seat Width	14"	15"	16"	17"	18"	19"	
12"	8881048722-NL A	--	--	--	--	--	
13"	8881048723-NL A	--	--	--	--	--	
14"	8881048724-NL A	8881048731-NL A	8881048739-NL A	8881048748-NL A	8881048755-NL A	--	
15"	8881048725-NL A	8881048732-NL A	8881048740-NL A	8881048749-NL A	8881048756-NL A	8881048762-NL A	
16"	8881048726-NL A	8881048733-NL A	8881048741-NL A	8881048750-NL A	8881048757-NL A	8881048763-NL A	
17"	8881048727-NL A	8881048734-NL A	8881048742-NL A	8881048751-NL A	8881048758-NL A	8881048764-NL A	
18"	8881048728-NL A	8881048735-NL A	8881048743-NL A	8881048752-NL A	8881048759-NL A	8881048765-NL A	
19"	--	--	8881048744-NL A	8881048753-NL A	8881048760-NL A	8881048766-NL A	
20"	--	--	8881048745-NL A	8881048754-NL A	8881048761-NL A	8881048767-NL A	
21"	--	--	--	--	8881082651-NL A	--	
22"	--	--	--	--	8881082652-NL A	--	

Specify frame color when ordering.