




# Side Frames without Hardware - Invacare Models




Invacare MG and  
Invacare MGI with  
Fixed Frame  
(Before 3/99)  
Inside Back Tube  
Mounted Axle Lug




Invacare MG and  
Invacare MGI with  
Fixed Frame  
(After 2/99)  
Outside Mounted  
Axle Lug




Invacare MG with  
Removable Arm  
(Before 3/99)  
"S" Frame Style



Invacare MG and  
Invacare MGI with  
Removable Arm  
(After 2/99)  
"Box" Frame Style



Invacare MG and  
Invacare MGI with  
Permanent Arm  
(After 2/99)  
"Box" Frame Style



Invacare MG with  
Permanent Arm  
(Before 3/99)  
"S" Frame Style

# Side Frames without Hardware - Invacare Models

## Invacare MGI (International Accounts Only)

Part Number	Description
1086145-NLA See Note(s) A	NLA - Frame, Box Style with Hardware, Permanent Arm, Chrome - Right
1086146-NLA See Note(s) A	NLA - Frame, Box Style with Hardware, Permanent Arm, Chrome - Left
1086149-NLA See Note(s) A	NLA - Frame, Box Style with Hardware, Removable Arm, Chrome - Right
1086150-NLA See Note(s) A	NLA - Frame, Box Style with Hardware, Removable Arm, Chrome - Left
1086147-NLA See Note(s) A	NLA - Frame with Hardware, Outside Mounted Lug, Fixed Frame - Right
1086148-NLA See Note(s) A	NLA - Frame with Hardware, Outside Mounted Lug, Fixed Frame - Left

NOTE: A - Frames include head tube bearings.

## Invacare MG (After 2/99)

Part Number	Description
1086145-NLA See Note(s) A	NLA - Frame, Box Style with Hardware, Permanent Arm, Chrome - Right
1086146-NLA See Note(s) A	NLA - Frame, Box Style with Hardware, Permanent Arm, Chrome - Left
1086149-NLA See Note(s) A	NLA - Frame, Box Style with Hardware, Removable Arm, Chrome - Right
1086150-NLA See Note(s) A	NLA - Frame, Box Style with Hardware, Removable Arm, Chrome - Left
1086147-NLA See Note(s) A	NLA - Frame with Hardware, Outside Mounted Lug, Fixed Frame - Right
1086148-NLA See Note(s) A	NLA - Frame with Hardware, Outside Mounted Lug, Fixed Frame - Left

NOTE: A - Frames include head tube bearings.

# Side Frames without Hardware - Invacare Models

## Invacare MG (Before 3/99)

Part Number	Description
1092382-NLA See Note(s) A	NLA - Frame, "S" Style with Hardware, Permanent Arm, Chrome - Right
1092383-NLA See Note(s) A	NLA - Frame, "S" Style with Hardware, Permanent Arm, Chrome - Left
1092384-NLA See Note(s) A	NLA - Frame, "S" Style with Hardware, Removable Arm, Chrome - Right
1092385-NLA See Note(s) A	NLA - Frame, "S" Style with Hardware, Removable Arm, Chrome - Left
1092386-NLA See Note(s) A	NLA - Frame with Hardware, Inside Back Tube Mounted Lug Fixed Frame, Chrome - Right
1092387-NLA See Note(s) A	NLA - Frame with Hardware, Inside Back Tube Mounted Lug Fixed Frame, Chrome - Left

NOTE: A - Frames include head tube bearings.