


High Mount Push to Lock Wheel Locks

High Mount Push to Locks



High Mount Push to Lock Wheel Locks

High Mount Push to Locks

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1046406-NLA	NLA - Wheel Lock, High Mount, Push to Lock - Right	1	1
1	A	1046407-NLA	NLA - Wheel Lock, High Mount, Push to Lock - Left	1	1
1		1125459	Wheel Lock Assembly, Push/Pull - Right	1	1
1		1125460	Wheel Lock Assembly, Push/Pull - Left	1	1
2		1068210	Kit, Wheel Lock Attaching Hardware (1" Tubing)	1	1
			Screw, Allen Head (1/4-20 x 1")	1	2
			Clamp, Half, Non-Threaded	1	1
			Clamp, Half, Threaded	1	1

NOTE: A - No Longer Available. Order Push/Pull Wheel Locks 1125459 RT and 1125460 LT as a pair to match.