


# High Mount Push to Lock Wheel Locks

## High Mount Push to Locks



# High Mount Push to Lock Wheel Locks

## High Mount Push to Locks

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1046406-NLA	NLA - Wheel Lock, High Mount, Push to Lock - Right	1	1
1	A	1046407-NLA	NLA - Wheel Lock, High Mount, Push to Lock - Left	1	1
1		1125459	Wheel Lock Assembly, Push/Pull - Right	1	1
1		1125460	Wheel Lock Assembly, Push/Pull - Left	1	1
2		1068210-NLA	Kit, Wheel Lock Attaching Hardware (1" Tubing)	1	1
			Screw, Allen Head (1/4-20 x 1")	1	2
			Clamp, Half, Non-Threaded	1	1
			Clamp, Half, Threaded	1	1
NOTE: A - No Longer Available. Order Push/Pull Wheel Locks 1125459 RT and 1125460 LT as a pair to match.					