


Rigid Tubular Footplate Tapered Front Only



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1026911	Screw, Socket Set Cone Point (1/4-20 x 1/2")	1	2
2		1024749	Clamp	1	2
3		1026124	Collar, Anti-Rattle	1	2
4		See Chart	#4 Rigid Tubular Footplate	1	1

Rigid Tubular Footplate Tapered Front Only

#4 Rigid Tubular Footplates

Chair Width	Part Number
14"	1043272-NLA
15"	1043273-NLA
16"	1043274-NLA
17"	1043275-NLA
18"	1043276-NLA
19"	1043277-NLA
20"	1043278-NLA
21"	1082974
22"	1082980