

Diagram illustrating the assembly of the backrest crossmember. The components are numbered as follows:

- 1: Main vertical support tube.
- 2: Upper connector or bush.
- 3: Lower connector or bush.
- 4: Small cylindrical component, likely a pin or spacer.
- 5: Long horizontal tube, likely the backrest crossmember.
- 6: Small rectangular component, likely a pin or spacer.
- 7: Small cylindrical component, likely a pin or spacer.
- 8: Small cylindrical component, likely a pin or spacer.
- 9: The entire assembly, labeled "backrest crossmember reference only".

(Not available on chairs with AT99 arms)

Push Handles

Bolt On and Quick Release

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1032209-NLA	NLA - Handle, Push, Black - Right	1	1
1		1032210-NLA	NLA - Handle, Push, Black - Left	1	1
2		1032211-NLA	NLA - Clamp, Threaded, Top	1	1
3		1032212-NLA	NLA - Clamp, Non-Threaded, Bottom	1	1
4		1027097	Plug Button, Black (3/4")	1	1
5		1028431-NLA	NLA - Grip, Foam, Black (5")	1	1
6		1032217-NLA	NLA - Screw, Set Flat Point (1/4-20 x 1")	1	1
7		1026907	Washer (1/4 x 1/2 x 1/16")	1	2
8		1025713-NLA	NLA - Screw, Hex Head (1/4-20 x 1")	1	2
9		1032213-NLA	NLA - Push Handle Assembly w/ Hardware, Bolt On (Pair)	2	1
10		1006229	Grip, Foam (8")	1	1
11		1045470	Pin, Quick Release with Lanyard	1	1
12	B	16348X000-NLA	NLA - Rivet (3/16 x 7/16")	1	1
13		1046290-NLA	NLA - Handle, Push with Quick Release Pin, Black	1	1
NOTE: A - Includes items 1-8 B - Includes Push Handle and items 10-12					