


Back Canes with Push Handles



Back Canes with Push Handles

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1006229	Grip, Foam (8")	1	1
2		1027097	Plug Button, Black (3/4")	1	1
3		1025195	NLA - Locknut (1/4-20)	1	1
4		1008432	Screw, Socket Button Head (1/4-20 x 1-1/4")	1	1
5	A,E	1054773-NLA	NLA - Cane, Back, Straight w/ Push Handle and Hardware (10-14" H)	1	1
			Grip, Foam (8")	1	1
			Plug Button, Black (3/4")	1	1
			NLA - Back Cane w/ Push Handle, Straight (10-14" H)	1	1
			Screw, Socket Button Head (1/4-20 x 1-1/4")	1	1
			NLA - Locknut (1/4-20)	1	1
5	B,E	1054774-NLA	NLA - Cane, Back, Straight w/ Push Handle and Hardware (12-16" H)	1	1
			Grip, Foam (8")	1	1
			Plug Button, Black (3/4")	1	1
			NLA - Back Cane w/ Push Handle, Straight (12-16" H)	1	1
			Screw, Socket Button Head (1/4-20 x 1-1/4")	1	1
			NLA - Locknut (1/4-20)	1	1
5	C,E	1027850-NLA	NLA - Cane, Back, Straight w/ Push Handle and Hardware (12-16" and 14-18" H)	1	1
			Grip, Foam (8")	1	1
			NLA - Cane, Back, Straight, Black - Right/Left	1	1
			Screw w/ Washer (#10 x 7/8")	1	1
			Screw, Socket Button Head (1/4-20 x 1-1/4")	1	1
			NLA - Locknut (1/4-20)	1	1
5	D,E	1044031-NLA	NLA - Back Cane w/ Hardware, Straight w/ Push Handle (14-18" and 16-20" H)	1	1
			NLA - Back Cane, Straight w/ Push Handle	1	1
			Grip, Foam (8")	1	1
			Screw, Socket Button Head (1/4-20 x 1-1/4")	1	1
			NLA - Locknut (1/4-20)	1	1
NOTE: A - 10-14" Back Height B - 12-16" Back Height C - 14-18" Back Height D - 16-20" Back Height E - Includes items 1-4					