

Pronto® 31 **Powered Mobility Vehicle**



Yes, you can.

Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 13-682. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-333-6900

For Faxed In Orders Call:

1-800-678-4682 or

1-440-366-0709

For Technical Service Call:

1-800-832-4707

Literature Can Be Found at:

www.invacare.com or

1-800-828-6272

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.


Name	Phone Numbers

© 2020 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

TABLE of CONTENTS

Description	Page Number
01. Controller and Frame Hardware	3
02. Anti-Tip Wheel Assembly	5
03. Drive Wheels, Casters and Forks	7
04. Motors	9
05. Base Shroud and Battery Box	11
06. Seat Post	13
07. Seat and Seat Base	14
08. Armrests	16
09. Positioning Strap	18
10. Joystick	20
11. Footplate and Footplate Bracket	22
12. Batteries, Cables, and Chargers	24

Controller and Frame Hardware




SCP_004A

Controller and Frame Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186175-NLA	NLA - Kit, Headtube Bearing w/ Hardware (Pair) Includes: (4) ea. Bearing (2) ea. Spacer (2) ea. Nylon Nut (2) ea. Cap End	1	1
2		1186176	Kit, Non-Ex Controller w/ Hardware Includes: (1) ea. PG nVR2 Controller (2) ea. Phillips Head Screw	1	1
	A	1186839	Programmer, Controller	1	1
NOTE: A - Not shown in illustration.					

Anti-Tip Wheel Assembly




SCP_006A

Anti-Tip Wheel Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186177-NLA	NLA - Kit, Anti-Tip Wheel Assembly w/ Hardware Includes: (2) ea. Wheel (2) ea. Washer (2) ea. Nylon Nut (2) ea. Screw	1	1



Drive Wheels, Casters and Forks



Drive Wheels, Casters and Forks

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186182	Kit, Drive Wheel FF Assembly w/ Hardware Includes: (1) ea. Drive Wheel (1) ea. M10 Nut (1) ea. Washer	1	1
2		1186178	Kit, Fork, Caster w/ Hardware Includes: (1) ea. Fork (1) ea. Fork Stud (1) ea. Nylon Nut	1	1
3		1186179	Kit, Solid Caster Wheel Assembly Includes: (1) ea. Wheel Assembly (2) ea. Washer (1) ea. Screw (1) ea. Nylon Nut	1	1



Motors



Motors

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186180	Kit, Motor/Gearbox - Left Includes: (1) ea. Left Motor and Gearbox (3) ea. Screws	1	1
1		1186181	Kit, Motor/Gearbox - Right Includes: (1) ea. Right Motor and Gearbox (3) ea. Screws	1	1
2		1186350	Kit, Drive Motor Hardware Includes: (3) ea. Screw (1) ea. M10 Nut (1) ea. Washer	1	1

Base Shroud and Battery Box




SCP_010A

Base Shroud and Battery Box

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186184	Kit, Black Shroud Includes: (1) ea. Black Shroud	1	1
1		1186185	Kit, Blue Shroud Includes: (1) ea. Blue Shroud	1	1
2		1186183	Kit, Battery Box Includes: (1) ea. Battery Box (6) ea. Screw (1) ea. Reuseable Hardware	1	1
3	A	1186186	Kit, Fuse Includes: (2) ea. Fuse	1	1
NOTE: A - Not shown in illustration.					


Seat Post



SCP_011A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186187	Kit, Seat Post w/ Hardware Includes: (1) ea. Seat Post (1) ea. Hex Head Bolt (2) ea. Concave Washer (1) ea. Nylon Nut	1	1

Seat and Seat Base




SCP_020A

Seat and Seat Base

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186188	Seat Assembly Includes: (1) ea. Rear Base Plate (1) ea. "L" Type Release Lever (1) ea. Right Lever Tip, Black (1) ea. Spring Lever (3) ea. Knob (1) ea. Square Plug (1) ea. Right Locating Lever (1) ea. Set Screw (1) ea. Spacer (1) ea. Hex Screw (1) ea. Spring Washer (1) ea. Push Tie-mount (1) ea. Headrest Assembly (1) ea. Headrest Cover, Black (1) ea. Backrest Cover, Black (1) ea. Seat Cover, Black	1	1
2	A	1186189	Kit, Replacement Knob Includes: (3) ea. Knob	1	1
		1186351	Kit, Seat Refurbish Includes: 2 (ea.) Armrest Pad 2 (ea.) Knob 4 (ea.) Countersunk Screws 1 (ea.) Headrest Cover 1 (ea.) Backrest Cover 1 (ea.) Seat Cover 1 (ea.) Instruction Sheet	1	1
NOTE: A - Not shown on illustration.					

Armrests




SCP_002A

Armrests

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186191-NLA	NLA - Kit, Armrest - Right Includes: (1) ea. Right Armrest Assembly	1	1
2		1186190-NLA	NLA - Kit, Armrest - Left Includes: (1) ea. Left Armrest Assembly	1	1
3		1186192	Kit, Armpad Replacement Includes: (1) ea. Knob (1) ea. Armrest Pad (2) ea. Socket Head Screw	1	1

Positioning Strap




SCP_014A

Positioning Strap

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186193	Kit, Positioning Strap w/ Hardware Includes: (1) ea. Seat Belt (2) ea. Washer (2) ea. Bushing (2) ea. Washer (2) ea. Screw	1	1


Joystick



Joystick

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186196	Kit, Joystick Mount w/ Hardware Includes: (1) ea. Tube (2) ea. Plug, Black (2) ea. Screw (2) ea. Screw (2) ea. Nut (1) ea. Joystick Bracket	1	1
2		1186197	Kit, Joystick Replacement Includes: (1) ea. Joystick	1	1

Footplate and Footplate Bracket




Footplate and Footplate Bracket

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1186194	Kit, Footplate w/ Hardware Includes: (1) ea. Footplate (1) ea. Brass Sleeve (1) ea. Hex Head Bolt (2) ea. Washer (1) ea. Nylon Nut	1	1
2		1186195	Kit, Footplate Bracket w/ Hardware Includes: (1) ea. Screw (1) ea. Nut (1) ea. Footplate Bracket (2) ea. Screw (2) ea. Nylon Nut	1	1

Batteries, Cables, and Chargers

Item 3 is Temporarily Not Available, Update to follow.



SCP_021B

Batteries, Cables, and Chargers

Item 3 is Temporarily Not Available, Update to follow.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1148830	Charger, Off-Board, 24 DC (5 AMP)	1	1
2		1144798	Cable, Charger	1	1
4		1186199	Kit, Battery Harness Includes: (1) ea. Battery Wire	1	1
NOTE: A - Cable is not included with this item, please order 1144798 when replacing the Charger.					

Invacare Corporation

www.invacare.com



Yes, you can.

One Invacare Way
Elyria, Oh USA
44036-2125
800-333-6900

Technical Support
800-832-4707
Customer Service
800-333-6900

Literature Orders
800-333-6900

© 2020 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.
(Revised 04/14/20)

Form No. 13-682