

Front Riggings Reference Chart

For AT5544, 60HD, 70NHD, PW93, PW93E, 60MFX

NOTE: Calf pads with attaching hardware not included on legrest.

NOTE: Pivot/ Slide tube attaching hardware not included.

NOTE: Pivot/Slide tube and footplates do not include heel loops. See footrest and legrest accessories.

1651 = Composite footplate with 9-7/8" curved pivot/slide tube.
AT5543 = Adjustable angle footplate with 9-7/8" straight pivot/slide tube.
1450 = Deluxe aluminum footplate with "L" style pivot/slide tube.
1350 = Extra large footplate with "L" style pivot/slide tube.

Support Assemblies

Pivot/Slide tube and footplates are interchangeable with the support assemblies listed.

Rigging Type	AT5544	60HD	70NHD	PW93	PW93E	60MFX
Right	8881045291 See Note(s) A	8881053158 See Note(s) A	8881053149 See Note(s) A	8881070776 See Note(s) A	8881070768-NL A See Note(s) A	8881026998 See Note(s) A
Left	8881045292 See Note(s) A	8881053157 See Note(s) A	8881053148 See Note(s) A	8881070777 See Note(s) A	8881070769-NL A See Note(s) A	8881026999 See Note(s) A

NOTE: A - Specify frame color when ordering, see page titled "Frame Finishes and Touch-Up Paints".

Pivot/Slide Tube and Footplates

Footplate Type	1651 13-14" W	1651 15-16" W	1651 17-24" W	1350 17-24" W	AT5543 12-13" W	AT5543 14-15" W	AT5543 16-24" W
Right	1057859	1057869	1057879	1026395	1060442	1060444	1062756
Left	1057860	1057870	1057880	1026396	1060443	1060445	1062757

Front Riggings Reference Chart

For AT5544, 60HD, 70NHD, PW93, PW93E, 60MFX

2" Longer Pivot/Slide Tube and Footplates

Footplate Type	1651 13-14" W	1651 15-16" W	1651 17-24" W	1350 17-24" W	AT5543 12-13" W	AT5543 14-15" W	AT5543 16-24" W
Right	1057861	1057871	1057881	1086515	N/A	1062758	1062760
Left	1057862-NLA	1057872	1057882	1086517	N/A	1062759	1062761

4" Longer Pivot/Slide Tube and Footplates

Footplate Type	1651 13-14" W	1651 15-16" W	1651 17-24" W	1350 17-24" W	AT5543 12-13" W	AT5543 14-15" W	AT5543 16-24" W
Right	1058419-NLA	1058421	1058423	1086516	N/A	N/A	N/A
Left	1058420-NLA	1058422	1058424	1086518	N/A	N/A	N/A

Support Assembly

Rigging Type	AT5544
Right	8881045291 See Note(s) A
Left	8881045292 See Note(s) A
NOTE: A - Specify frame color when ordering, see page titled Frame Finishes and Touch-Up Paints.	

Pivot/Slide Tube and Footplates

Footplate Type	1450 16-24" W
Right	1028723-NLA
Left	1028724-NLA