


Seat Upholstery



Seat Upholstery Hardware Kit

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A,B	1026933-NLA	NLA - Package, Phillips Head Screw (#10-32 x 3/4")	14	1
	A	See Chart	Seats	1	1

NOTE: A - Not Shown
B - Includes 14 ea. Screws (10-32 x 3/4")

Chair Width	14" Deep	15" Deep	16" Deep	17" Deep	18" Deep	19" Deep
14"	1032080	1032081	1032082	1032083	1032084	--
15"	1032085	1032086	1032087	1032088	1032089	--
16"	1032090	1032091	1032092	1032093	1032094	1046478-NLA
17"	1032095	1032096	1032097	1032098	1032099	1046479-NLA
18"	1032100	1032101	1032102	1032103	1032104	1046480-NLA
19"	--	--	1032105	1032106	1032107	1046481-NLA
20"	--	--	1032108	1032109	1032110	1046482-NLA
21"	--	--	--	--	1083022-NLA	--
22"	--	--	--	--	1083028-NLA	--