



Half Arms - AT99



Half Arms - AT99

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881053570	Kit, Arm Pad w/ Screws, Desk Length (Pair)	2	1
1		1032244	Grip, Foam Hand (9" L)	1	1
2		1055205	Package, Screw Phillips Pan Head (#10-32 x 1-3/8")	2	1
3		1032202	Half Arm Tube - Right	1	1
3		1032203	Half Arm Tube - Left	1	1
4		1032206-NLA	NLA - Socket, Arm, Rear - Right	1	1
4		1032207-NLA	NLA - Socket, Arm, Rear - Left	1	1
5		1018687	Sleeve, Bearing	1	1
6		1055348	Screw, Hex Head (1/4-28 x 1-3/4")	1	1
7		1065362	Locknut (1/4-28)	1	3
8		1026907	Washer (1/4 x 1/2 x 1/16")	1	2
9		1026900	Washer (1/4 x 23/32 x 1/16")	1	2
10	B	1055349	Screw, Hex Head (1/4-28 x 2")	1	1
11		40157X020	Spacer, Coved (1/4 x 1 x 9/16")	1	2
12		1032195	Kit, Half Arm Mounting Hardware (Pair)	2	1
	C	1059656	Instruction Sheet, Half Arms	1	1
NOTE: A - Specify U67 upholstery color when ordering. B - Includes items 4-11 and Instruction Sheet. C - Not Shown					