

Wheels without Handrims - 1/2" Axle (Chairs built Before 6/24/99)

Wheel Style	Wheel Size	Tire Type	Single Wheel Assembly	Tire	Tube Insert	Bearing	Bearing Spacer
U1763BK	20 x 1-1/4"	Urethane	1038364 See Note(s) B, A	1037774	N/A	40117X022	1015756
1763BK	20 x 1-3/8"	Pneumatic	1026422 See Note(s) A	1017519	1017520	40117X022	1015756
4763BK	20 x 1-3/8"	Pneumatic w/ Flat Free Insert	1026423 See Note(s) B, A	1017519	1023322	40117X022	1015756
U1764BK	22 x 1-1/4"	Urethane	1032116 See Note(s) A	1028250	N/A	40117X022	1015756
1764BK	22 x 1-3/8"	Pneumatic	1026420 See Note(s) A	1018544	1018545	40117X022	1015756
4764BK	22 x 1-3/8"	Pneumatic w/ Flat Free Insert	1026421 See Note(s) B, A	1018544	1023323	40117X022	1015756
7 Spoke U2244BK	24 x 1-1/4"	Urethane	1032118 See Note(s) A	1026904	N/A	40117X022	1015756
2244BK	24 x 1-3/8"	Pneumatic	1026362-NLA	40176X060	40176X061	40117X022	1015756
4244BK	24 x 1-3/8"	Pneumatic w/ Flat Free Insert	1026363 See Note(s) B	40176X060	1010906	40117X022	1015756
U2222 9-Spoke	24 x 1"	Urethane	1039320-NLA See Note(s) C	1039304	N/A	40117X022	1015756
U2222 9-Spoke	24 x 1"	Urethane	1080380-NLA See Note(s) D, B	1080400	N/A	40117X022	1015756

NOTE: A - When ordering 20" or 22" wheels, a wheel lock shoe spacer may also be required. See Wheel Locks.
B - Composite Wheels with Flat Free Inserts or Urethane Tires have a plug which caps an opening in the rim.
C - Before 11/20/98
D - Between 11/20/98 to 6/23/99