


# Pivot/Slide Tube w/ Footplate Assy. - Composite - 1651 (Not Including Heel Loops)



## 60 Degree (13-19") Curved End

Right/Left	Small Long 14"	Small 2" Longer 14"	Medium Long 15"-16"	Medium 2" Longer 15"-16"	Large Long 17"-22"	Large 2" Longer 17"-22"
Right	1028391	1032627-NLA	1028387	1032625-NLA	1028383	1032623-NLA
Left	1028392	1032628-NLA	1028388	1032626-NLA	1028384	1032624-NLA

# Pivot/Slide Tube w/ Footplate Assy. - Composite - 1651 (Not Including Heel Loops)

**60 Degree MFX (9-3/4-14") Curved End**

Right/Left	Small Short 14"	Medium Short 15"-16"	Large Short 17"-22"
Right	1028389-NLA	1028385-NLA	1028381-NLA
Left	1028390-NLA	1028386-NLA	1028382-NLA

## 70 Degree (14-19") Straight

Right/Left	Small Short 14"	Small 2" Longer 14"	Medium Short 15"-16"	Medium 2" Longer 15"-16"	Large Short 17"-22"	Large 2" Longer 17"-22"
Right	1042784-NLA	1042790	1042786	1042792	1042788	1042794
Left	1042785-NLA	1042791	1042787	1042793	1042789	1042795

## 70 Degree Tapered (14-1/2-19-1/2") Straight

Right/Left	Small Short 16"-17"	Small 2" Longer 16"-17"	Medium Short 18"-19"	Medium 2" Longer 18"-19"	Large Short 20"-22"	Large 2" Longer 20"-22"
Right	1042784-NLA	1042790	1042786	1042792	1042788	1042794
Left	1042785-NLA	1042791	1042787	1042793	1042789	1042795

# Pivot/Slide Tube w/ Footplate Assy. - Composite - 1651 (Not Including Heel Loops)

## Fixed Front End

Right/Left	Small Short 14"	Medium Short 15"-16"	Large Short 17"-22"
Right	1061968-NLA	1061970-NLA	1061972-NLA
Left	1061969-NLA	1061971-NLA	1061973-NLA

## Elevating Legrest (16-20") Curved End

Right/Left	Small Long 14"	Small 2" Longer 14"	Medium Long 15"-16"	Medium 2" Longer 15"-16"	Large Long 17"-22"	Large 2" Long 17"-22"
Right	1028391	1032627-NLA	1028387	1032625-NLA	1028383	1032623-NLA
Left	1028392	1032628-NLA	1028388	1032626-NLA	1028384	1032624-NLA

After 11/93 Straight Slide Tubes are used on 70 degree Tapered and 70 degree Non-Tapered Supports. All chairs made before 12/93 had Slide Tubes with the Curved End, use the 60 degree chart for reference.