


Crossbraces & Hardware

Chairs with heavy duty packages have dual crossbraces. Before ordering crossbraces for these chairs please contact specials dept.



Crossbraces & Hardware

Chairs with heavy duty packages have dual crossbraces. Before ordering crossbraces for these chairs please contact specials dept.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1056328	Plug Button, Black (1")	1	4
2		See Chart	#2 Crossbrace Pair Kits (10-12" Deep)	1	1
2		See Chart	#2 Crossbrace Pair Kits (14-18" Deep)	1	1
2		See Chart	#2 Single Crossbraces (10-14" W x 10-14" D)	1	1
2		See Chart	#2 Single Crossbraces (14-20" W x 14-18" D)	1	1
3		1026106-NLA	NLA - Washer (1/4 x 5/8 x 1/16")	1	4
4		7930334-NLA	NLA - Locknut (1/4-28)	25	2
5		1016263	Spacer, Coved (3/4 x 1 x 1/2")	1	4
6		1025704	Washer, Flat (5/16 x 5/8 x 1/16")	1	2
7		7930330-NLA	NLA - Locknut (5/16-24)	25	1
8	A	1025642-NLA	NLA - Link, Pivot, Black (10 - 12" W)	1	2
8	B	1025645	Link, Pivot, Black (13 - 15" W)	1	2
8	C	1025647	Link, Pivot, Black (16 - 17" W)	1	2
8	D	1025649	Link, Pivot, Black (18 - 20" W)	1	2
9		1055348	Screw, Hex Head (1/4-28 x 1-3/4")	1	2
10		1025705-NLA	NLA - Screw, Hex Head (5/16-24 x 3-1/4")	1	1
11		7930313	Package, Washer, Nylon (5/16 x 7/8 x 1/8")	25	1

NOTE: A - 10-12" Wide
 B - 13-15" Wide
 C - 16-17" Wide
 D - 18-20" Wide

#2 Crossbrace Pair Kits (10-12" Deep)

Chair Width	10" Deep	11" Deep	12" Deep
10"	1031540-NLA	1033810-NLA	1031542-NLA
12"	1031541-NLA	1033811-NLA	1031543-NLA

Crossbraces replaced on chairs built before 9/91 need the following retrofit kits which include: the right and left crossbrace and center bolt with hardware.
 Specify Frame Color when ordering.

#2 Crossbrace Pair Kits (14-18" Deep)

Chair Width	14" Deep	15" Deep	16" Deep	17" Deep	18" Deep
14"	1031509-NLA	1031510-NLA	1031511-NLA	1031512-NLA	1031513-NLA
15"	1031514-NLA	1031515-NLA	1031516-NLA	1031517-NLA	1031518-NLA
16"	1031519-NLA	1031520-NLA	1031526-NLA	1031521-NLA	1031522-NLA
17"	1031523-NLA	1031524-NLA	1031525-NLA	1031527-NLA	1031528-NLA
18"	1031529-NLA	1031530-NLA	1031531-NLA	1031532-NLA	1031533-NLA

Crossbraces & Hardware

Chairs with heavy duty packages have dual crossbraces. Before ordering crossbraces for these chairs please contact specials dept.

#2 Crossbrace Pair Kits (14-18" Deep)

Chair Width	14" Deep	15" Deep	16" Deep	17" Deep	18" Deep
19"	N/A	N/A	1031534-NLA	1031535-NLA	1031536-NLA
20"	N/A	N/A	1031537-NLA	1031538-NLA	1031539-NLA

Crossbraces replaced on chairs built before 9/91 need the following retrofit kits which include: the right and left crossbrace and center bolt with hardware.
Specify Frame Color when ordering.

Crossbraces & Hardware

Chairs with heavy duty packages have dual crossbraces. Before ordering crossbraces for these chairs please contact specials dept.

#2 Single Crossbraces (10-14" W x 10-14" D)

Chair Width	10" Deep	11" Deep	12" Deep	13" Deep	14" Deep
10" Right	1025657-NLA	1028241-NLA	1025659-NLA	N/A	N/A
10" Left	1025658-NLA	1028242-NLA	1025660-NLA	N/A	N/A
12" Right	1025661-NLA	1028243-NLA	1025663-NLA	1032150-NLA	1032152-NLA
12" Left	1025662-NLA	1028244-NLA	1025664-NLA	1032151-NLA	1032153-NLA
13" Right	1032154-NLA	1032156	1032158-NLA	1032160-NLA	1032162-NLA
13" Left	1032155-NLA	1032157-NLA	1032159-NLA	1032161-NLA	1032163-NLA
14" Right	1032164-NLA	1032166-NLA	1032168-NLA	1032170-NLA	1025665
14" Left	1032165-NLA	1032167-NLA	1032169-NLA	1032171-NLA	1025666

Crossbraces & Hardware

Chairs with heavy duty packages have dual crossbraces. Before ordering crossbraces for these chairs please contact specials dept.

#2 Single Crossbraces (14-20" W x 14-18" D)

Chair Width	14" Deep	15" Deep	16" Deep	17" Deep	18" Deep
14" Right	1025665	1028203-NLA	1025667-NLA	1028205-NLA	1028207-NLA
14" Left	1025666	1028204-NLA	1025668-NLA	1028206-NLA	1028208-NLA
15" Right	1025669-NLA	1028209-NLA	1025671-NLA	1028211-NLA	1028213-NLA See Note(s) C
15" Left	1025670-NLA	1028210-NLA	1025672-NLA	1028212-NLA	1028214-NLA See Note(s) C
16" Right	1025673-NLA See Note(s) A	1028215	1025675	1028217	1028219
16" Left	1025674-NLA See Note(s) A	1028216	1025676	1028218	1028220
17" Right	1025677-NLA See Note(s) B	1028221	1025679	1028223	1028225
17" Left	1025678-NLA See Note(s) B	1028222	1025680	1028224	1028226
18" Right	1025681-NLA	1028227	1025683	1028229	1028231
18" Left	1025682-NLA	1028228-NLA	1025684	1028230	1028232
19" Right	N/A	N/A	1025685	1028233	1028235
19" Left	N/A	N/A	1025686	1028234	1028236
20" Right	N/A	N/A	1025687	1028237	1028239
20" Left	N/A	N/A	1025688	1028238	1028240

NOTE: A - No longer available as an individual component. Use 1031519 as the replacement part, see "Crossbrace Pair Kits" to find the replacement part number.
 B - No longer available as an individual component. Use 1031523 as the replacement part, see "Crossbrace Pair Kits" to find the replacement part number.
 C - No longer available as an individual component. Use 1031518 as the replacement part, see "Crossbrace Pair Kits" to find the replacement part number.