


Seat Height Adjustment Slide Tubes - Pronto M71 with PTO



Seat Height Adjustment Slide Tubes - Pronto M71 with PTO

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1118688	Slide Tube Assembly, Seat Height Adjustment, Rear, Black - Pronto M71	1	1
2		1118675-NLA	NLA - Bracket, Telescoping PTO, Black - Rear	1	2
3		1088932	Screw, Socket Button Head (1/4-20 x 3/8")	1	4
4		1118369-NLA	NLA - Decal, PTO Prop Rod	1	1
5		1118674-NLA	NLA - Tube, Slide, Seat Height Adjustment, Front - Pronto M71	1	2
6		1026158	Locknut (5/16-18)	1	4
7		1021344	Screw, Socket Head (5/16-18 x 1-1/2")	1	4
8		1105026	Plug Button, Head Tube, Black	1	2
9	B	1118691-NLA	NLA - Seat Height Adjustment Slide Tubes w/ Hardware, Pronto M71	1	1

NOTE: A - Includes items 2-4
B - Includes items 1-8