


# Frame Hardware

## Xtra/Jr.

(Built After 3/03)



## Frame Hardware Xtra/Jr. (Built After 3/03)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1024786	Bushing, Pivot	1	12
2		1025703-NLA	NLA - Screw, Phillips Oval Head (#10-32 x 1")	1	2
3		1007541	Clamp, Seat	1	2
4		1059284-NLA	NLA - Screw, Socket Button Head (1/4-28 x 1-1/4")	1	8
5		1025765	Locknut (1/4-28)	1	8
6		1024749	Clamp	1	2
7		1026911-NLA	NLA - Screw, Socket Set Cone Point (1/4-20 x 1/2")	1	2
8		1013870-NLA	NLA - Anti-Rattle (1")	1	2
9	A	1107510	Hardware, Frame, Front (1 Pair)	1	1
10		1029501	Plug Button, Black (1")	1	2
11		1026888	Washer, Nylon (3/4 x 1 x 1/32")	1	4
12		1024787-NLA	NLA - Tube, Insert, Upper	1	2
13		1024788	Tube, Insert, Lower	1	2
14		1052957	Screw, Hex Head w/ Patch (1/4-20 x 2")	1	2
15		1026907	Washer (1/4 x 1/2 x 1/16")	1	6
16		1036627-NLA	NLA - Plug Button (1")	1	2
17		1055366	Locknut (1/4-28)	1	4
18	D	1054571	Kit, Sector Block - Right	1	1
18	D	1054572	Kit, Sector Block - Left	1	1
19		1009861	Spacer, Coved (1/4 x 11/16 x 1/2")	1	2
20		1026099	Spacer, Coved (1/4 x 11/16 x 1/2")	1	2
21		1113501	Screw, Eccentric (1/4-28 x 1-11/16")	1	4
22	B	1116626-NLA	NLA - Hardware, Frame (Pair)	2	1
22	C	1116627-NLA	NLA - Hardware, Frame (Pair)	2	1

NOTE: A - Fixed Front Frames. Includes Insert Tubes and items 1-8.  
B - 14-18" Deep Swingaway Frames. Includes items 1-5, 10-21.  
C - 12" Deep Frames. Includes items 1-5, 12,13,15,17, 19-21. Does not include items 10,14,16, & 18.  
D - Includes sector block and hardware, items 14 and 15.