


# Bed Ends - Manual & Semi Electric

## 5100

### Fixed Height

discontinued 6/03



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		51000M087-NLA	NLA - Panel, Decorative - Head End	1	2
1		51000M076-NLA	NLA - Panel, Decorative, Foot End	1	2
1	D	1116647	Panel, High Impact, Head End	1	2
1	D	1116648	Panel, High Impact, Foot End	1	2
2		51000M091	NLA - Trim, Decorative - Head End	1	1
2		51000M080-NLA	NLA - Trim, Decorative, Foot End	1	1
3		51000M082	Socket, Caster	1	2
4		51000M581	Guide, Leg	1	4
5		51000M908-NLA	NLA - Screw (#6 x 5/16")	1	12
6	A	51000M000-NLA	NLA - Caster Assembly, Swivel w/ Brake (3")	1	1
6	B	54000M446	Caster Assembly, Swivel w/ Brake (3")	1	1
6	A	51000M001	Caster Assembly, Swivel (3")	1	1
6	B,G	54000M445-NLA	NLA - Caster Assembly, Swivel (3")	1	1
6	A,E	5100M008-NLA	NLA - Package, Caster	1	1
7		52000M214-NLA	NLA - Fixed Height Head End Assembly without Casters	1	1
7		52000M213-NLA	NLA - Foot End Assembly, Fixed Height without Casters	1	1
	C	51000M426	S-Clip, Panel	1	1

NOTE: A - Before 10/94  
 B - After 9/94  
 C - Not Shown  
 D - Two (2) Panels are need PER Bed End. Head End takes two (2) panels and Foot End takes two (2) panels.  
 E - Includes: Two (2) Casters with Brakes and Two (2) Casters without Brakes.  
 G - No Longer Available, replace with 54000M446.