


# Careguard by Invacare

## Single Drive Bed Motors and Components



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		7001M267-NLA	NLA - Motor, Head	1	1
2		7001M263-NLA	NLA - Foot, Drive Tube	1	1
3		54000M515	Screw, Slotted Hex Head Tap (1/4-20 x 3/4")	1	3
4		5300M054	Mount, Motor	1	1
5		54000M721-NLA	NLA - Cover, Front	1	1
6		54000M641	Screw, Hex Head (5/16-18 x 1")	1	4
7		54000M686-NLA	NLA - Bushing, Sleeve	1	2
8		54000M493	Clamp, Cable, Steel (9/16")	1	1
9		53000M434	Locknut (5/16-18)	1	2
10		54000M492-NLA	NLA - Guard, Motor (5-3/4" Tall)	1	1