



Drive Shafts - 80" and 84"

Both Drive Shafts have the same part number. Drive Shafts manufactured before 4/15/03 can ONLY be used on older style bed models; Drive Shafts manufactured after 4/14/03 can be used on EITHER the old style or the new, IVC models.



Drive Shafts - 80" and 84"

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	7000M713	Drive Shaft Assembly (84")	1	1
	A	7000M585	Drive Shaft Assembly (80")	1	1
		40188X009	Button, Single Snap	1	1
		53000M734	Cap, End, Protective	1	2
<p>NOTE: A - Both Drive Shafts have the same part number (before and after dates). The shafts made before 4/15/03 can only be used on older style bed models; shafts made after 4/14/03 can be used on EITHER the old or the NEW, IVC models.</p> <p>Note: These items are not shown.</p>					