


Back Cane and Spreader Bar Assembly

Recline Only (2GR) and Tilt/ Recline (2GTR)

NOTE: When changing back cane heights, new Back Assembly Mounting Hardware may be required. See Back Assembly Mounting Hardware page



Back Cane and Spreader Bar Assembly

Recline Only (2GR) and Tilt/ Recline (2GTR)

NOTE: When changing back cane heights, new Back Assembly Mounting Hardware may be required. See Back Assembly Mounting Hardware page

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Spreader Bar	1	1
2		See Chart	#2 Back Canes	1	2
3		1078214	Screw, Hex Head w/ Patch (1/4-20 x 1/2")	1	4
4		1026907	Washer (1/4 x 1/2 x 1/16")	1	4
5		1082699	T-Nut (1/4-20)	1	4
6		1086587-NLA	NLA - Link, Curved, Black	1	2
7		1031673-NLA	NLA - Washer, Nylon (35/64 x 1 x 1/16")	1	4
8		1086323-NLA	NLA - Screw, Shoulder Hex Socket Head (1/2 x 5/8", 3/8-16)	1	2
9		1021851-NLA	NLA - Bolt, Socket Head Shoulder (1/2 x 1-1/2")	1	2
10		1062604	Screw, Shoulder Hex Socket (3/8 x 1-1/4", 5/16-18)	1	2
11		1025704	Washer, Flat (5/16 x 5/8 x 1/16")	1	6
12		1070398-NLA	NLA - Screw, Hex Head (5/16-18 x 1-1/2")	1	2
13		1082725-NLA	NLA - Block, Slide	1	2
14		1082724-NLA	NLA - Rod, Slide	1	2
15		1039860	Screw, Hex Head w/ Patch (5/16-18 x 3/4")	1	2
16		1026158	Locknut (5/16-18)	1	4
17		40363X000	Hex Nut (3/8"-16)	1	4
18		1065214	Washer, Flat (3/8 x 11/16 x 1/16")	1	2
19		1082829-NLA	NLA - Plug Button, Back Cane, Black	1	2
20	A	1087001-NLA	NLA - Kit, Back Cane and Spreader Bar Mounting Hardware	2	1

NOTE: A - Includes items 3-19

Item #1 Spreader Bar

16" Wide	17" Wide	18" Wide	19" Wide	20" Wide	21" Wide	22" Wide
1086101-NLA	See Note See Note(s) A	1086103-NLA	See Note See Note(s) A	See Note See Note(s) A	1086106-NLA	See Note See Note(s) A

NOTE: A - No longer available as a service part.

Item #2 Back Canes

Back Height

Side	20-21" High	22-23" High	24-25" High	26" High
Right	1086948-NLA	1086950-NLA	1086952-NLA	1086954-NLA
Left	1086949-NLA	1086951-NLA	1086953-NLA	1086955