


Stroller Handles



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1133698	Grip, Foam Hand (9" L)	1	1
2		1056328	Plug Button, Black (1")	1	1
3		1043722-NLA	NLA - Handle, Stroller, Black - Right	1	1
3		1043723-NLA	NLA - Handle, Stroller, Black - Left	1	1
4		16348X000-NLA	NLA - Rivet (3/16 x 7/16")	1	1
5		1045470	Pin, Quick Release with Lanyard	1	1
6		15680X013	Anti-Rattle (1")	1	1
7		1044714-NLA	NLA - Stroller Handle Assembly	2	1