



# Shrouds and Attaching Hardware

After 7/13/03



**Side Shrouds**



**Battery Shrouds**

# Shrouds and Attaching Hardware

## After 7/13/03

### Magic Black

Right	Left	Pair
1112817-NLA See Note(s) A	1112818-NLA See Note(s) A	1112825-NLA
NOTE: A - No longer available. Use 1112825 as its replacement.		

### Dark Blue Pearl

Right	Left	Pair
1112811-NLA See Note(s) A	See Note	1112823
NOTE: A - No longer available as a service part in the field.		

### Candy Red Pearl

Right	Left	Pair
1112805-NLA	1112806-NLA	1112824-NLA

## Shrouds and Attaching Hardware After 7/13/03

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Side Shrouds Candy Red Pearl	1	1
1		See Chart	#1 Side Shrouds Dark Blue Pearl	1	1
1		See Chart	#1 Side Shrouds Magic Black	1	1
2		1025451	Screw, Slotted Hex Washer Head (#10-32 x 5/8")	1	8
3		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	12
4		1028755	Screw, Phillips Pan Head (#10-32 x 5/8")	1	4
5		1106366-NLA	NLA - Plate, Shroud -Right	1	1
5		1106367-NLA	NLA - Plate, Shroud -Left	1	1
6	B	1106343-NLA	NLA - Hardware, Shroud, Side (1 Pair)	2	1
7		1109032-NLA	NLA - Decal, Shroud - Pronto M71	1	1
8		1104797	Decal, Drive/Push	1	1
9	C	1108464-NLA	NLA - Kit, Shroud Clamp Hardware - Right	1	1
9	C	1108465-NLA	NLA - Kit, Shroud Clamp Hardware - Left	1	1
10		1109036-NLA	NLA - Bracket, Clamp Shroud Mounting, Black - Right	1	1
10		1109037-NLA	NLA - Bracket, Clamp Shroud Mounting, Black - Left	1	1
11		1026907	Washer (1/4 x 1/2 x 1/16")	1	1
12		1025736	Locknut (1/4-20)	1	1
13		1109035-NLA	NLA - Clamp, Shroud	1	1
14		1025713-NLA	NLA - Screw, Hex Head (1/4-20 x 1")	1	1
15		1095489	Decal, Invacare 3-D Logo (2-1/2" L)	1	2
16		1117240	Screw, Thumb (#10-32 x 1/2")	1	4
17		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	4
18	D	1123240-NLA	NLA - Battery Shroud w/ Decal, Black - Rear	1	1
19		1118381-NLA	NLA - Label, Warning Battery Terminal	1	1
20	A	1123237-NLA	NLA - Battery Shroud with Decals, Black - Front	1	1

NOTE: A - This shroud encompasses the controller as well as the batteries.  
 B - Includes items 2-5  
 C - Includes items 10-14  
 D - This shroud encompasses the battery charger as well as the batteries.