

Footrest & Legrest Accessories

Heel Loops (1346)

Footplate Type	Heel Loop Single #	Heel Loop Pair #	Heel Loop with Ankle Strap - Right Single #	Heel Loop with Ankle Strap - Left Single #	Heel Loop with Ankle Strap Pair #	Attaching Hardware Pair #
Composite (1651) 14" W (16" Wide Tapered)	1028412-NLA	1035912-NLA	1035908-NLA	1035909-NLA	1035906	1035931
Composite (1651) 16" W (18" Wide Tapered)	1017621	1035917-NLA	1035919	1035920	1035918	1035931
Composite (1651) 18" W (20" Wide Tapered)	1017621	1035917-NLA	1035932-NLA	1035933-NLA	1035930	1035931
Standard Aluminum (1650) 14" Wide	1025234-NLA	1035905	1035908-NLA	1035909-NLA	1035906	1035907
Standard Aluminum (1650) 16-20" Wide	1017621	1035917-NLA	1035919	1035920	1035918	1035907
Deluxe Aluminum (1450) 16-20" Wide	1017621	1035917-NLA	1035919	1035920	1035918	1086112-NLA
X-Large Aluminum (1350) 14" Wide	1025234-NLA	1035905	1035908-NLA	1035909-NLA	1035906	1055241
X-Large Aluminum (1350) 16-20" Wide	1017621	1035917-NLA	1035919	1035920	1035918	1055241

Footrest & Legrest Accessories

Heel Loop Hardware

Part Number	Consists of:
1035931	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/4") (2) Spacer (1/4 x 11/32 x 1-3/4") (2) Hex Nut w/ Nylon Cap (1/4-20) (2) Flat Washer (1/4 x 5/8 x 1/16")
1055241	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/4"), (2) Spacer (1/4 x 11/32 x 1-3/4"), and (2) Locknut (1/4-20").
1035907	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/2"), (2) Spacer (1/4 x 11/32 x 1-3/4"), and (2) Locknut (1/4-20").
1086112-NLA	(2) Screw, Phillips Pan Head (1/4-20 x 2-1/2"), (2) Spacer (1/4 x 11/32 x 1-3/4"), and (2) Locknut (1/4-20").

Toe Loops (1345)

Footplate Type	Toe Loop w/o Hardware Single #	Toe Loop w/o Hardware Pair #	Attaching Hardware Part Number	Attaching Hardware Consists of
1651, 1650, 1450, 1350, Footplates	04572X025	1037087	1037089	(4) Screw, Phillips Pan Head (#10-32 x 7/8"), (4) Washer (1/4 x 3/4 x 3/32"), and (4) Locknut (#10-32).

Footrest & Legrest Accessories

Miscellaneous Accessories

Part Number	Description
1026987	Impact Guard (8") (IMPACTG)
See Note(s) D 1143050	Non-Removable Device (1820)
1032003	"H" Heel Strap (9342) 14" Wide
1032002	"H" Heel Strap (9342) 16-20" Wide
1035902	Calf Strap (1337)
8881006725-NLA See Note(s) A, B	Footrest Support Panel (1348) 14-16" Wide
8881006726 See Note(s) B, A	Footrest Support Panel (1348) 18" Wide
8881006727-NLA See Note(s) B, A	Footrest Support Panel (1348) 20" Wide
8881035816-NLA See Note(s) C, A	Legrest Support Panel (1348) 14-16" Wide
8881035817-NLA See Note(s) A, C	Legrest Support Panel (1348) 18" Wide
8881035818-NLA See Note(s) C, A	Legrest Support Panel (1348) 20" Wide
NOTE: A - Specify upholstery color when ordering. See page titled "Frame Finishes, Touch-Up Paints and Upholstery Colors." B - For PW93 footrest C - For P904A, PAS4A and PAL4A Legrest D - For 70 deg. Tapered Footrest.	
NOTE: These items are not shown	