


Casters, Forks, and Hardware



Casters, Forks, and Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1001843	Cap, Head Tube, Black	1	1
2		20225X000	Locknut (7/16-20)	1	2
3		1065530	Bearing, Flange (7/16")	1	1
4		1074930	Sleeve w/ Bearing, Black, Adult/Hemi	1	1
5		1065589	Fork and Stem Assembly, Adult/Hemi	1	1
6		1016032	Washer, Steel (15/32 x 5/8 x 3/64")	1	2
7		57667P007	Screw, Hex Head (7/16-20 x 3-1/4")	1	1
8		1070811	Wheel Assembly, Semi-Pneumatic, Urethane, Light Grey (8" x 1-3/4")	1	1
9		1020592	Spacer, Bearing (7/16 x 5/8 x 1-3/4")	1	1
10		1006049	Bearing (7/16")	1	2