




# Disk Style Wheel Locks (Before 4/15/07)

**NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.**



## Disk Style Wheel Locks (Before 4/15/07)

**NOTE: If your chair was built with a serial # starting with "03H" or "03I", please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 Ext. 3722.**

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C,F,F1	1110225-NLA	NLA - Wheel Lock Assembly, Disk Push to Lock, Black - Right	1	1
1	C,F,F1	1110226-NLA	NLA - Wheel Lock Assembly, Disk Push to Lock, Black - Left	1	1
1	D,F,F1	1110227-NLA	NLA - Wheel Lock Assembly, Disk Pull to Lock, Black - Right	1	1
1	D,F,F1	1110228-NLA	NLA - Wheel Lock Assembly, Disk Pull to Lock, Black - Left	1	1
2	E	1110160-NLA	NLA - Tip, Plastic Wheel Lock Handle, Black	1	1
3	A	1017677	Kit, Hardware, Wheel Lock (Pair)	1	1
4			Nut, Hex Flang (1/4-20)	1	2
5			Bolt, Carriage (1/4-20 x 1-1/4")	1	2
	B,E,G	1058900	Wheel Lock Extension Handle Assembly, Black (8")	1	1
	B,C,G	1101452	Wheel Lock Assembly, Push to Lock, Black - Right	1	1
	B,C,G	1101453	Wheel Lock Assembly, Push to Lock, Black - Left	1	1
	B,D,G	1101454	Wheel Lock Assembly, Pull to Lock, Black - Right	1	1
	B,D,G	1101455	Wheel Lock Assembly, Pull to Lock, Black - Left	1	1

NOTE: A - Includes items 4 and 5 each with a quantity of 2.  
 B - Not shown.  
 C - Standard on Adult Seat to Floor models.  
 D - Standard on Hemi Seat to Floor models.  
 E - Used on Disk-Style Wheel Locks.  
 F - Disk Style Wheel Locks are no longer available as a Service Part. Must convert to Link-Pivot Style Wheel Locks listed. Order Right/Left side to match. Can use existing attaching hardware.  
 F1 - Additional Wheel Lock information can be found under category titled "Wheel Locks - Link-Pivot Style (After 4/14/07)".  
 G - Link-Pivot Style Wheel Locks