

Heel Loops

Heel Loops and Hardware

Description	Front Rigging Type	Part Number
Heel Loops (1 Pair)	Footrests and Legrests	1035929
Heel Loop Hardware (1 Pair) See Note(s) B	Footrests	1035926 See Note(s) B
Heel Loop Hardware (1 Pair) See Note(s) A	Legrests	1035907 See Note(s) A

NOTE: A - Includes 2 ea. Spacers, Screws, and Locknuts
B - Includes 2 ea. Spacers, Screws, Washers, and Locknuts.