


90 Degree Fixed Supports



90 Degree Fixed Supports

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881029426-NL A	NLA - Tube, Support, Color	1	1
2		1001762	Screw, Socket Button Head (1/4-20 x 1-3/4")	1	1
3		1017148	Pivot, Footrest	1	1
4		1027097	Plug Button, Black (3/4")	1	1
5		1025195	Locknut (1/4-20)	1	1
6		1001762	Screw, Socket Button Head (1/4-20 x 1-3/4")	1	1
7		1025736	Locknut (1/4-20)	1	1
8	A	8881029462-NL A	NLA - Footrest Support Assembly w/ Hardware, 90 Degree, Color - Right/Left	2	1

NOTE: A - Specify frame finish when ordering.