


# Crossbraces and Hardware



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1055309	Screw, Socket Button Head (1/4-28 x 1-1/4")	1	10
2		1025765	Locknut (1/4-28)	1	10
3		See Chart	#3 Crossbraces (14-18" Wide)	1	5
3		See Chart	#3 Crossbraces (19-22" Wide)	1	5

## #3 Crossbraces (14-18" Wide)

# Crossbraces and Hardware

## #3 Crossbraces (14-18" Wide)

14" Wide	15" Wide	16" Wide	17" Wide	18" Wide
1121541	1121542	1121543	1121544	1121545

# Crossbraces and Hardware

## #3 Crossbraces (19-22" Wide)

19" Wide	20" Wide	21" Wide	22" Wide
1121546	1121547	1121548	1121549