


# Drive Wheels and Casters



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1099771-NLA	NLA - Wheel Assembly, Rear (10 x 3.50) - Black Rim	1	1
2		1088720-NLA	NLA - Tire w/ Foam Fill Insert, Rear, Grey (8" x 2")	1	1
3		1099739-NLA	NLA - Hardware, Rear Wheel Mounting Includes: 1 ea. Key 1 ea. Washer (M8 x 45 x 3) 1 ea. Washer, Locking Tab 2 ea. Washer, Spring (M6 x 10) 2 ea. Screw, Phillips Pan Head (M6 x 1.0 x 16) 1 ea. Screw, Hex Head Cap (M8 x 1.25 x 20)	1	1
4	A	1115709-NLA	NLA - Caster with Nut, Rear (5")	1	1
4	A	1115777-NLA	NLA - Caster with Nut, Rear (5")	2	1
5	A	1115710-NLA	NLA - Caster with Nut, Front (6")	1	1
5	A	1115778-NLA	NLA - Caster with Nut, Front (6")	2	1

**NOTE:** A - Casters are covered under warranty for a period of 1 year.