

Wheel Lock Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1055306-NLA	NLA - Screw, Button Head Cap (1/4-28 x 7/8")	1	3
2		1122956-NLA	NLA - Plate, Wheel Lock Clamp - Right	1	1
2		1122957-NLA	NLA - Plate, Wheel Lock Clamp - Left	1	1
3		1025765	Locknut (1/4-28)	1	5
4		1026103	Washer, Nylon (1/4 x 1/2 x 1/32")	1	6
5		1122953-NLA	NLA - Shoe, Wheel Lock	1	1
6		1020363-NLA	NLA - Spacer (1/4 x 1/2 x 1/2")	1	1
7		1018733-NLA	NLA - Link, Stop, Black	1	1
8		1055304	Screw, Socket Button Head (1/4-28 x 5/8")	1	2
9		1122966-NLA	NLA - Handle, Wheel Lock, Black - Right	1	1
9		1122967-NLA	NLA - Handle, Wheel Lock, Black - Left	1	1
10		1028411-NLA	NLA - Spacer (1/4 x 1/2 x 5/32")	1	1
11		1122819-NLA	NLA - Plate, Wheel Lock Mounting	1	1
12		1122954	Bracket, Wheel Lock Mounting, Black - Right	1	1
12		1122955-NLA	NLA - Bracket, Wheel Lock Mounting, Black - Left	1	1
13		1030950-NLA	NLA - Bolt, Carriage (1/4-20 x 1-1/2")	1	1
14		1123392-NLA	NLA - Screw, Round Head Square Neck (1/4-20 x 3/4")	1	1
15		1006490	Tip, Plastic, Black	1	1
16		1025736	Locknut (1/4-20)	1	2
17	A	1121974-NLA	NLA - Wheel Lock Assembly w/ Attaching Hardware -Right	1	1
17	A	1121975-NLA	NLA - Wheel Lock Assembly w/ Attaching Hardware -Left	1	1

NOTE: A - Includes items 1-16