


# Oxygen Holder



# Oxygen Holder

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1116948	Holder, Oxygen	1	1
2		1025736	Locknut (1/4-20)	1	4
3		1121216-NLA	NLA - Bracket, Oxygen Holder, Black	1	1
4		1025708-NLA	NLA - Bolt, Hex Head Cap (1/4-20 x 3/4")	1	2
5		1018739	Washer (1/4 x 5/8 x 3/32")	1	1
6		1024796-NLA	NLA - Screw, Hex Head (1/4-20 x 2")	1	2
7		1021866-NLA	NLA - Bag, Poly (26" L x 15" W)	1	2
8		1122135-NLA	NLA - Instruction Sheet, Oxygen Holder	1	1
9		1063960	Carton	1	1
10	A	1122318-NLA	NLA - Oxygen Holder Assembly	1	1

NOTE: A - Includes items 1-9