



# Battery Harnesses

After 10/7/03

Battery Covers are no longer used on chairs After 10/7/03



Battery Cable Harness

# Battery Harnesses

After 10/7/03

Battery Covers are no longer used on chairs After 10/7/03

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1122247	Harness, Battery w/ 75 AMP Fuse and Red Connector	1	1
2		1122246-NLA	NLA - Harness, Battery - Front	1	1
3		1122245	Harness, Battery w/ 75 AMP Fuse and Black Connector	1	1
4		1122244-NLA	NLA - Harness, Battery - Rear	1	1
5		16721X000	Locknut (1/4-20)	1	4
6		1062246	Washer, Locking Tab	1	2
7		1102337-NLA	NLA - Screw, Hex Flange Whiz Lock(1/4-20 x 5/8")	1	4
	A	98007X005	Tie Wrap (11-1/2" L)	1	6
	A	1118378-NLA	NLA - Instruction Sheet, Battery Harness	1	1
8	B	1122249-NLA	NLA - Harness, Battery w/ Terminal Hardware	1	1
9		1116381	Cable, Battery Extension, Black Connector	1	1

NOTE: A - Not Shown  
B - Includes items 1-7 and Cables Ties, Instruction Sheet.