

Side Frames with Hardware

Model K2000

Side Frame w/ Hardware	w/ Removable Arm	w/ Permanent Arm
Right	1123171-NLA See Note(s) A	1123173-NLA See Note(s) A
Left	1123172-NLA See Note(s) A	1123174-NLA See Note(s) A
NOTE: A - Complete frame assembly includes handgrip, caster journal bearings, all tube plugs/caps and the lower frame insert tube with attaching hardware.		

Model K3000

Side Frame w/ Hardware	Part Number
Right	1123191-NLA See Note(s) A
Left	1123192-NLA See Note(s) A
NOTE: A - Complete side frame assembly includes handgrip, caster journal bearings, plastic h-blocks, rivets & washers and all tube plugs/caps.	

Model K4000

Side Frame w/ Hardware	Part Number
Right	1123204-NLA See Note(s) A
Left	1123205-NLA See Note(s) A
NOTE: A - Complete side frame assembly includes caster journal bearings, plastic h-blocks, rivets & washers, all tube plugs/caps and back cane assembly with handgrip, attaching screw and knob.	

Model K7000

Side Frame w/ Hardware	Part Number
Right	1123210-NLA See Note(s) A
Left	1123211-NLA See Note(s) A
NOTE: A - Complete side frame assembly includes handgrip, caster journal bearings, all tube plugs/caps and the lower frame insert tube with attaching hardware.	