

Casters, Forks and Mounting Hardware

Model K2000, K3000, K4000

Part Number	Fork and Stem Assembly includes:
1123177-NLA	Fork & Stem Assembly (1) 8" x 1" Caster Assembly (1) Fork Stem Locknut (1) Fork Stem Washer (1)
See Note(s) A	
NOTE: A - 8x1" caster includes tire, bearings, hub spacer, axle bolt and locknut.	

Model K7000

Part Number	Fork and Stem Assembly includes:
1123214-NLA	Fork & Stem Assembly (1) 8 x 2" Caster Assembly (1) Fork Stem Locknut (1) Fork Stem Washer (1)
See Note(s) A	
NOTE: A - 8x2" caster includes tire, bearings, hub spacer, axle bolt and locknut.	