



Battery Harnesses & Hardware



chairs with MKIV Electronics



chairs with MK5 Electronics

Battery Harnesses & Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1114498-NLA	NLA - Battery Terminal Hardware (GP22, MKIV)	1	1
2	B	1115350	Harness, Wire, Battery (75 Amp) - Rear (MKIV)	1	1
3	B	1115349	Harness, Wire, Battery (75 Amp) - Front (MKIV)	1	1
4	C,D	1122814-NLA	NLA - Harness, Battery, 100 Amp w/ Terminal Hardware - MK5 (Red Connector)	1	1
NOTE: A - Battery Terminal Hardware is NLA. See item 4 for replacement. B - Individual components for the front and rear battery harnesses are not available C - Includes All battery harnesses, battery install instruction sheet, 1 ea. battery lift strap, 2 ea. washer locking tabs, 2 ea. 1/4-20 flange locknuts, 2 ea. 1/4-20 x 7/8" hex flange lock screws, 2 ea. 1/4-20 locknuts, and 6 ea. cable ties. D - Battery Extension Cable (p/n 1122199) may be needed as well. See Controller, Cables, Programmer, and Hdwr. (MK5) category.					