


A.S.B.A. Back Upholstery, Nylon



Back Upholstery w/ Dual Liners

Back Width	16"-17" High	18"-19" High	20" High
18" Wide	1045182	1045214	1068954
19" Wide	1045183	1045215	1068955
20" Wide	1045184	1045216	1068956
21" Wide	1045185	1045217	1068957
22" Wide	1045186	1045218	1068958
23" Wide	1045187-NLA	1045219-NLA	1068959-NLA
24" Wide	1045188-NLA	1045220-NLA	1068960-NLA

A.S.B.A. Back Upholstery, Nylon

Back Upholstery Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	1026298	Package, Screw with Washer (#10 x 7/8") (Pair)	2	1
	A	98007X005	Tie Wrap (11-1/2" L)	1	2
	A	See Chart	Back Upholstery w/ Dual Liners	1	1

NOTE: A - Not Shown