




Front, Top, and Rear Shrouds and Attaching Hardware



Front Shroud



Top Shroud Cover



Rear Shroud

Front, Top, and Rear Shrouds and Attaching Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1095489	Decal, Invacare 3-D Logo (2-1/2" L)	1	1
2		1070574-NLA	NLA - Stud, Wing Head (1" Long)	1	1
3		1093312-NLA	NLA - Spring, Torsion Ejector	1	1
4		1070581-NLA	NLA - Washer, Nylon Cupped (9/32 x 19/32 x 3/32")	1	1
5		1070580-NLA	NLA - Washer, Nylon Retainer (1/4 x 1/2 x 1/32")	1	1
6		1112571	Shroud with Hardware - Front	1	1
7		1070509-NLA	NLA - Fastener, Loop with Adhesive Back (1" W x 2" L)	1	4
8	D	1079203	Decal, Warning Pinch Point	1	2
9		1118397-NLA	NLA - Label, Warning Wiring Dual, 22NF Battery	1	1
10	B	1035900	Decal, Transportation	1	1
11		1122816	Shroud, Top w/ Decals	1	1
12		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	1
13	C	1025451	Screw, Slotted Hex Washer Head (#10-32 x 5/8")	1	1
14		1123225	Kit, Shroud w/ Hardware, Rear, Candy Red Pearl	1	1
14		1123226	Kit, Shroud w/ Hardware, Rear, Magic Black	1	1
14	C	1123227	Kit, Shroud w/ Hardware, Rear, Dark Blue Pearl	1	1
NOTE: A - Includes Shroud plus items 1-5 B - Includes Shroud plus items 7-10 C - Includes Shroud plus items 1, 12, and 13 D - Part #1118397 is no longer available as a service part individually. Order item 11 as a replacement.					