




Battery Harnesses - MK5 Electronics

Chairs built After 10/7/03

Battery Covers are no longer used on chairs After 10/7/03





Battery Cable Harness




Battery Lifting Strap (UI)

chairs with MK5 NX Controller



Battery Cable Harness



Battery Lifting Strap (UI)

chairs with MK5 NX-80 or MK5 EX Controller

Battery Harnesses - MK5 Electronics

Chairs built After 10/7/03

Battery Covers are no longer used on chairs After 10/7/03

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1119579-NLA	NLA - Harness, Battery Rear w/ 75 AMP Fuse	1	1
2		1119580-NLA	NLA - Harness, Battery Front w/ 75 AMP Fuse	1	1
3		16721X000	Locknut (1/4-20)	1	4
4		1062246	Washer, Locking Tab	1	2
5		1102337-NLA	NLA - Screw, Hex Flange Whiz Lock(1/4-20 x 5/8")	1	4
	A	98007X005	Tie Wrap (11-1/2" L)	1	6
	A	1118378-NLA	NLA - Instruction Sheet, Battery Harness	1	1
6	B,E	1119585-NLA	NLA - Harness, Battery, 75 Amp. w/ Terminal Hdwr, MK5 & MK6i, Black Connector	1	1
7	B	1116381	Cable, Battery Extension, Black Connector	1	1
8	C,D	1122814-NLA	NLA - Harness, Battery, 100 Amp w/ Terminal Hardware - MK5 (Red Connector)	1	1
9	C	1122199	Cable, Battery Extension, Red Connector	1	1
10	F	1127455-NLA	NLA - Strap, Battery Lifting (U1 Batteries)	1	1

NOTE: A - Not Shown
 B - Used with MK5 NX Electronics.
 C - Used with MK5 NX-80 & MK5 EX Electronics.
 D - Includes All battery harnesses, battery install instruction sheet, 1 ea. battery lift strap, 2 ea. washer locking tabs, 2 ea. 1/4-20 flange locknuts, 2 ea. 1/4-20 x 7/8" hex flange lock screws, 2 ea. 1/4-20 locknuts, and 6 ea. cable ties.
 E - Includes items 1-5 and Cables Ties, Instruction Sheet.
 F - Battery Lifting Strap is standard when chair is ordered without "batteries onboard" option