


Split Non-Adjustable Stroller Handles -1483



Split Non-Adjustable Stroller Handles -1483

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1133698	Grip, Foam Hand (9" L)	1	1
2		1056328	Plug Button, Black (1")	1	1
3	A	1061996-NLA	NLA - Handle, Stroller, Black - Right	1	1
3	A	1061997-NLA	NLA - Handle, Stroller, Black - Left	1	1
4	B	1043722-NLA	NLA - Handle, Stroller, Black - Right	1	1
4	B	1043723-NLA	NLA - Handle, Stroller, Black - Left	1	1
5		16348X000-NLA	NLA - Rivet (3/16 x 7/16")	1	2
6		1061993-NLA	NLA - Tube, Insert	1	1
7		1045470	Pin, Quick Release with Lanyard	1	1
8		1013870-NLA	NLA - Anti-Rattle (1")	1	1
9	A,C	1062566	Stroller Handle Assembly	1	1
10	B,D	1044714-NLA	NLA - Stroller Handle Assembly	1	1
NOTE: A - Used with Fixed Height Back Canes B - Used with Adjustable Height Back Canes C - Includes items 1-3,5-7 D - Includes items 1,2,4,5,7,8					