


Crossbraces - Hardware - Seat Extensions



Crossbraces - Hardware - Seat Extensions

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1029501	Plug Button, Black (1")	1	4
2		1030949	Screw, Hex Head (1/4-20 x 1-1/2")	1	4
3		1012180	Bracket, Pivot, Black	1	2
4		7930323-NLA	NLA - Package, Locknut (1/4-20)	25	1
5		1012172-NLA	NLA - Clamp, Saddle	1	2
6	A,B	1016252-NLA	NLA - Crossbrace, Black (16" W)	1	2
6	A,B	1016253	Crossbrace, Black (18" W)	1	2
6	A,C	1061299-NLA	NLA - Crossbrace, Black (16" W)	1	2
6	A,C	1061300-NLA	NLA - Crossbrace, Black (18" W)	1	2
7	D,F	1014666-NLA	NLA - Link, Pivot, Black (14 - 16" W)	1	2
		7200188-NLA	NLA - Pivot Link w/ Hardware (14 - 16" W)	1	1
7	E	1014668	Link, Pivot, Black (17 - 18" W)	1	2
8		1024793-NLA	NLA - Screw, Hex Head (1/4-20 x 1-3/4")	1	2
9		1098625	Screw, Hex Head (5/16-24 x 3-1/2")	1	1
10		16308X000	Washer, Flat (5/16 x 11/16 x 1/16")	1	2
11		1016263	Spacer, Coved (3/4 x 1 x 1/2")	1	4
12		7930313	Package, Washer, Nylon (5/16 x 7/8 x 1/8")	25	1
13		7930330-NLA	NLA - Locknut (5/16-24)	25	1
14		1037913-NLA	NLA - Bushing (5/16 x 7/16 x 9/32")	1	1
15		1025704	Washer, Flat (5/16 x 5/8 x 1/16")	1	1
16		7930333	Locknut (1/4-20)	25	1
17		1058616-NLA	NLA - Kit, Seat Extension - 1" (Pair)	2	1
17		1058617	Kit, Seat Extension - 2" (Pair)	2	1
18		1026298	Package, Screw with Washer (#10 x 7/8") (Pair)	2	1
18		1076284	Screw w/ Washer (#10 x 7/8")	1	2

NOTE: A - Decals not included with crossbraces.
 B - Before 4/1/96
 C - After 3/31/96
 D - 14-16" wide
 E - 17-19" wide
 F - Pivot Link is no longer available as a service part individually. Order part #7200188 Pivot Link w/ Hardware (not shown) as a replacement. Includes items 2-4, 7 and 8.