


Seat Pan and Hardware - A.S.B.A.



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
5	A	1070459-NLA	NLA - Kit, Seat Pan Mounting Hardware	1	1
1			Screw, Hex Head (1/4-20 x 1-5/8")	1	2
2			Screw, Hex Head (1/4-20 x 1-1/2")	1	4
3			NLA - Spacer (1/4 x 3/8 x 3/16")	1	4
4			Locknut (1/4-20)	1	6
6		1039949-NLA	NLA - Fastener, Hook Adhesive Back (2" W x 16" L)	1	2
6		1044578-NLA	NLA - Fastener, Loop Adhesive Back (2" W x 16" L)	1	2
7		See Chart	#7 Seat Pan	1	1

NOTE: A - Includes Items 1-4

Seat Pan and Hardware - A.S.B.A.

#7 Seat Pan

Chair Width	16" Deep	17" Deep	18" Deep	19" Deep	20" Deep	21" Deep	22" Deep
18" Wide	1066015-NLA	1066026	1066037-NLA	1066048-NLA	1066059-NLA	1066070-NLA	1066081-NLA
19" Wide	1066016-NLA	1066027-NLA	1066038-NLA	1066049-NLA	1066060	1066071-NLA	1066082-NLA
20" Wide	1066017-NLA	1066028-NLA	1066039-NLA	1066050-NLA	1066061-NLA	1066072-NLA	1066083-NLA
21" Wide	1066018-NLA	1066029-NLA	1066040-NLA	1066051-NLA	1066062-NLA	1066073-NLA	1066084-NLA
22" Wide	1066019-NLA	1066030-NLA	1066041-NLA	1066052-NLA	1066063-NLA	1066074-NLA	1066085-NLA
23" Wide	1066020-NLA	1066031-NLA	1066042-NLA	1066053-NLA	1066064-NLA	1066075-NLA	1066086-NLA
24" Wide	1066021-NLA	1066032-NLA	1066043-NLA	1066054-NLA	1066065-NLA	1066076-NLA	1066087-NLA