


# Front Rigging Hanger



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1106315-NLA	NLA - Brace, Arm Socket	1	1
2		1025736	Locknut (1/4-20)	1	2
3		1025195	Locknut (1/4-20)	1	1
4		1027097	Plug Button, Black (3/4")	1	1
5		1106313	Hanger, Footrest Support - Right	1	1
5		1106314	Hanger, Footrest Support - Left	1	1
6		1003545	Screw, Socket Head (1/4-20 x 1-1/4")	1	2
7		60285X006	Screw, Socket Head (1/4-20 x 2-1/4")	1	1