




Walker (WalkLite)

Model 6300-5F, 6300-ATA and 6300-JRA



**Front and Rear
Leg Extensions
Sold in Pairs Only**



**Leg Extension
Shown for
Reference Only**

Walker (WalkLite)

Model 6300-5F, 6300-ATA and 6300-JRA

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	1107793-NLA	NLA - Backrest, Nylon Upholstery	1	1
1	D	1043538-NLA	NLA - Nut, Push	1	4
1		1130784-NLA	NLA - Kit, Push Nut	10	4
2	B,H	1124254	Kit, Extension Front Leg Assembly	2	1
3		1053179	Glide	1	2
4		1040656	Tip, Rubber, Dark Grey (1-1/8")	1	2
5	C,E	1124255-NLA	NLA - Kit, Legs with Tips and Hardware (Pair)	2	1
5	C,F	1149805-NLA	NLA - Kit, Legs with Tips and Hardware (Pair)	2	1
6		1063124	Anti-Rattle, Grey (1-1/8")	1	4
7		1094945	Button, Snap	1	4
	G	1109367-NLA	NLA - Seat Hardware	1	1
8	G	1021421-NLA	NLA - Screw, Phillips Pan Head (1/4-20 x 2-1/4")	1	2
8	G	1025195	Locknut (1/4-20)	1	2
8	G	1025728	Washer (1/4 x 1/2 x 1/32")	1	2
8	G	1108477-NLA	NLA - Washer, Nylon (1/4 x 3/4 x 5/8")	1	2
9		1107737-NLA	NLA - Seat	1	1
10		1044944	Spring, Compression (31/64 x 1 x 3/64")	1	2
11		1044943	Fork	1	2

NOTE: A - Not Shown.
 B - Includes Leg Extension Tube, Locknut, Bushing, Washers, Wheel and Screw.
 C - Kit includes (2) Leg Extensions, (2) Anti-Rattles, (2) Glides, and (2) Tips.
 D - Individual item no longer available, refer to Kit Part #1130784.
 E - For use with 6300-5F and 6300-JRA.
 F - For use with 6300-ATA.
 G - Kit no longer available. Order each of item #8 in order to receive a complete kit.
 H - For use with 6300-5F, 6300-JRA and 6300-ATA.