


Battery Box, Shrouds, Chargers and Cables



Battery Box, Shrouds, Chargers and Cables

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	E	1030426-NLA	NLA - Box, Battery with Harness	1	1
1	F	1030427-NLA	NLA - Box, Battery with Harness	1	1
1	G	1033964-NLA	NLA - Box, Battery with Harness	1	1
1	H	1035477-NLA	NLA - Box, Battery with Harness	1	1
1	I	1033963-NLA	NLA - Box, Battery with Harness	1	1
1	J	1030426-NLA	NLA - Box, Battery with Harness	1	1
1	K	1030429-NLA	NLA - Box, Battery with Harness	1	1
1	L	1030428-NLA	NLA - Box, Battery with Harness	1	1
1	M	1035478-NLA	NLA - Box, Battery with Harness	1	1
1	N	1034614-NLA	NLA - Box, Battery with Harness	1	1
1	O	1030429-NLA	NLA - Box, Battery with Harness	1	1
1	P	1035222-NLA	NLA - Box, Battery with Harness	1	1
1	Q	1034862-NLA	NLA - Box, Battery with Harness	1	1
2	R	1029685-NLA	NLA - Lid, Battery Box	1	1
2	S	1030443-NLA	NLA - Lid, Battery Box	1	1
2	T	1035221-NLA	NLA - Lid, Battery Box	1	1
2	U	50000X025-NLA	NLA - Clip, Lock Over Latch	1	2
2	U	1019694-NLA	NLA - Slide Over Clip (Black)	1	4
2	U	1029687-NLA	NLA - Slide Latch, Grey	1	4
2	U	1010578-NLA	NLA - Spacer, Battery Clamp (3/16 x 5/8 x 2-5/8")	1	2
2	U	1010701-NLA	NLA - Rivet	1	4
2	U	16348X005-NLA	NLA - Washer, Steel (5/32 x 3/8 x 1/32")	1	2
2	U	47375X277-NLA	NLA - Fastener, Loop with Adhesive Back (1" W x 4" L)	1	2
3		1029686-NLA	NLA - Bottom, Box	1	1
4	V	See Page	See Footnote	1	2
5		12642X000	Screw (#8-18 x 5/8")	1	1
6	R,S	1023702-NLA	NLA - Harness, Battery with Fuse	1	1
6	T,W	1104224-NLA	NLA - Harness, Battery with Fuse	1	1
6		50000X166-NLA	NLA - Fuse - Buss 2-1/4" L (60 AMP)	1	1
6		1042740-NLA	NLA - Fuse - Shaumut 1-1/2" (30/50 AMP)	1	1
7		50000X123	Cable, Battery Jumper	1	1
8		1010681-NLA	NLA - Plug, Clamp	1	1
9		1013399-NLA	NLA - Cable, Charger Extension	1	1
10	Z	1098892-NLA	NLA - Fuse - Buss ATC-15 (15 AMP)	1	2
10	Y	1090117-NLA	NLA - Fuse - Buss Mini ATM-15 (15 AMP)	1	2
11		1095503-NLA	NLA - Package, Tie Wrap (5-5/8" L)	10	1
12		7930335-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 1/2")	25	2
13		7930332-NLA	NLA - Locknut, High Profile (#10-32)	25	2
14	D	1014758-NLA	NLA - Charger, Battery, 10 Amp. 115vac, Off Board, Dual Mode Square Plug, Lester	1	1
14	B,D	1053163-NLA	NLA - Charger, Battery, 8 Amp, 115vac, Off Board, Dual Mode, Square Plug, Lester	1	1
14	D	1085685-NLA	NLA - Charger, Battery, 8 Amp. 115vac, Off Board, Dual Mode, Square Plug w/ AC Power Line Cord, Dynamic	1	1

Battery Box, Shrouds, Chargers and Cables

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
14	B1,B2, D	1062556-NLA	NLA - Battery Charger, 8 Amp. 115vac (Off Board) Single Mode Round Plug w/ AC Power Line Cord and Square Plug Adapter Cable, CTE Corp.	1	1
15	AA	1089627-NLA	NLA - Cord, AC Power, 115VAC, 2 Prong Off-Board (11' 6" L)	1	1
15		1121086	Cord, AC Power, 115VAC, 3 Prong On-Board (6' 6-1/2" L)	1	1
16	X	1027965-NLA	NLA - Shroud, Center (14" W)	1	1
16	X	1027966-NLA	NLA - Shroud, Center (15" W)	1	1
16	X	1027967-NLA	NLA - Shroud, Center (16" W)	1	1
16	X	1027968-NLA	NLA - Shroud, Center (17" W)	1	1
16	X	1027969-NLA	NLA - Shroud, Center (18" W)	1	1
16	X	1027970-NLA	NLA - Shroud, Center (19" W)	1	1
16	X	1027971-NLA	NLA - Shroud, Center (20" W)	1	1
16	X	1027972-NLA	NLA - Shroud, Center (21" W)	1	1
16	X	1027973-NLA	NLA - Shroud, Center (22" W)	1	1
17		1028019-NLA	NLA - Decal, Center Shroud	1	1
18		1028362-NLA	NLA - Rivet, Snap	1	1
19		1010424-NLA	NLA - Cap, Anti-Wear	1	1
20		1022597	Screw, Phillips Pan Head Tap (#10-32 x 3/8")	1	2
21		1027263-NLA	NLA - Washer, Lock Internal Tooth (#10)	1	1
22		1030340-NLA	NLA - Shroud, Bottom, Mounting Bracket, Standard Height	1	1
22		1034906-NLA	NLA - Shroud, Bottom, Mounting Bracket (2" Lower)	1	1
23		1030444-NLA	NLA - Decal, Pulley/ Clutch Shroud - Action Arrow	1	1
23		1028020-NLA	NLA - Decal, Pulley/ Clutch Shroud - Action Ranger X	1	1
23		1028021-NLA	NLA - Decal, Pulley/ Clutch Shroud - Action Ranger II	1	1
24		1022649-NLA	NLA - Guide, Cable, Battery Box	1	1
25		1029291-NLA	Holder, Fuse, 1559 Dual Mode Charger	1	1
26		1028809-NLA	NLA - Timer, 1559 Dual Mode Charger	1	1
	A,C	1027962-NLA	NLA - Screw, Phillips Flat Head (10-32 x 1/2")	1	1
	A,C	1030446-NLA	NLA - Bracket, Shroud Mounting	1	1
27		1053831-NLA	NLA - Cable, Adapter	1	1
28		See Chart	# 28 Pulley/Clutch Shroud - Grey	1	1

Battery Box, Shrouds, Chargers and Cables

NOTE: A - On 2" Lower Chairs ONLY, the mounting bracket is used as a spacer to align the shroud when attaching to the chair. When using a 2" lower recliner, the screws are used to mount the bracket in a specific location prior to mounting the shroud.

AA - No longer available as a service part. Please use 1121086 as its replacement.

B - On 7/18/01, Lester 8 amp. battery charger was reinstated for service parts.

B1 - CTE charger introduced 4/1/04. Used on Gel Cell or AGM (Absorbed Glass Mat) batteries. For wet cell batteries refer to the Lester charger.

B2 - CTE charger looks exactly the same as Dynamic's charger except for the company name that is silk screened on the narrow side panel, where the charger plug is (hard wired) attached. Serial number starts with 04A00000001.

C - Not Shown

D - See Pictorial

E - Non-Recliners, 16" to 18" Wide with Sectional Back

F - Non-Recliners, 14" and 15" Wide with Sectional Back

G - Non-Recliners, 20" to 24" Wide with Permanent Uprights and Dual Crossbraces

H - 2" Lower Manual Recliners, 14" Wide

I - 2" Lower Manual Recliners, 16" Wide

J - 2" Lower Manual Recliners, 18" Wide

K - Manual Recliners, 16" to 18" Wide

L - Manual Recliners, 14" and 15" Wide

M - 2" Lower Manual Recliners, 14" Wide

N - 2" Lower Manual Recliners, 16" Wide

O - 2" Lower Manual Recliners, 18" Wide

P - Power Back Recliners, 16" to 18" Wide

Q - Labac, 16" to 18" Wide

R - Sectional Backs

S - Manual Recliners

T - Power Back Recliners

U - Part numbers 50000X025, 1019694, 1029687, 1010578, 1010701, 16348X005, and 47375X277 are available for individual sale or as part of "Miscellaneous Lid Hardware".

V - See page titled "Battery Box Hangers".

W - Labac

X - Introduced After 1/17/92

Y - After 3/10/03

Z - Before 3/11/03

Item # 28 Pulley/Clutch Shroud - Grey

Side	Non-Recline Standard Height	Non-Recline 2" Lower Height	Recliner Standard Height	Recliner Tall	Recliner 2" Lower Height
Right	1037648-NLA	1031881-NLA	1028683-NLA	1043253-NLA	1035132-NLA
Left	1037649-NLA	1031882	1028684-NLA	1043254-NLA	1035133-NLA

NOTE: Chairs manufactured after 1/17/92 went from black pulley/clutch shrouds to a grey. Shrouds should be ordered in pairs along with items 20, 21, 22 for chairs made before 1/18/92.